RESIN MODIFIED GLASS IONOMER
LUTING CEMENT

INSTRUCTIONS FOR USE

Riva Cem Automix is a self-curing, radiopaque, fluoride releasing paste/paste resin modified glass-ionomer luting material. It is indicated for permanent cementation of metal-based and strengthen-core ceramic restorations and orthodontic appliances.

Riva Cem Automix comes with the additional light tack cure option for early removal of excess cement and fluoride release for your cementation needs.

Riva Cem Automix is available in a universal light-yellow shade in double barrel syringes. The automix tip attachment offers convenience over traditional hand mixed cement systems.

INDICATIONS FOR USE:

Permanent cementation of:

Indicated material	Type of restorations
PFM	Crowns & bridges
Prefabricated/Cast	Posts
Metals	Crowns, bridges, inlays & onlays, orthodontic appliances & posts
Ceramics (high strength)	Crowns & bridges, inlays & onlays
Ceramics (low strength) ¹	Inlays
Keramik (höhere Festigkeit) ¹	Inlays
Keramik (niedrigere Festigkeit) ²	Inlays

¹High strength ceramics, e.g. Zirconia, Lithium Disilicate
²Low strength ceramics, e.g. Feldspathic Porcelains, Glass Ceramics

CONTRA-INDICATIONS:
1. Pulp capping.
2. Product may cause skin irritations to some people. In such cases, discontinue use and seek medical attention. Any persons having known resin allergies should immediately discontinue use.

INSTRUCTIONS FOR USE:
1. PREPARATION:
a) Prepare tooth using standard clinical guidelines.
b) For optimum adhesion, condition the tooth with a mild solution of polyacrylic acid (**Riva Conditioner**, SDI Ltd.). Rinse thoroughly with water.
c) In cases where pulp protection is necessary, use a calcium hydroxide.

2. MIXING PROCEDURE:
Paste / Paste ratio (ml/ml) 20/1.0
Working time (min, sec)² (from start of mixing) 1'30"

Excess cement:
- Tack light cure
- Self-curing gel-phase²
Setting time (min, sec)² (after seating the restoration) 5'00"
Temperatur Condition: (26 ± 1°C, 78 ± 1°F)
Radiopacity 1.5 mm Al

¹ Test conditions: temperature (23 ± 1°C, 73 ± 1°F), relative humidity (50 ± 10%), ISO 9917-2
² These times are approximate only, and may vary based on storage conditions, humidity, age of product and temperature.

a) Remove the syringe from the foil pack.
b) Remove the cap from the syringe.
c) Before each use, dispense a small amount to ensure even dispensing or to remove any dried out paste. Discard this material.

Automixing with SDI Riva Cem Automix tip
d) Attach a new mixing tip and secure it by turning it to the right.
e) Dispense and discard a small amount until an evenly mixed paste in a homogenous colour flows out of the tip.
Note: The mixing tip is for single use only. Discard the used mixing tip and place the original syringe cap onto the syringe. Syringe must be stored with the cap.

Handmixing with a spatula
d) Extrude required amount of pastes onto a mixing pad.
e) Place the mixing pad onto the middle of the spatula in the middle of the plunger to obtain the correct dispensing ratio of the two pastes.
f) Using a spatula, mix the two pastes together thoroughly for 15 seconds. Be careful not to incorporate air bubbles.

Note: Do not expose **Riva Cem Automix** paste directly to operatory light to avoid early setting of the material.

3. CEMENTATION PROCEDURE:
a) Ensure the internal surface of the restoration to be seated is cleaned, dried and prepared according to manufacturer's instructions.
b) Apply a thin layer of cement to the bonding surface.

c) Seat the restoration within 30 seconds of completion of mixing before any light curing is carried out.

Note: Avoid overfilling crowns.
2. Room temperatures above 23°C/73°F will reduce the working time.

Light tack cure any excess material for 5 seconds per surface using a high power LED curing light (460-480nm emission Radii-Cal, Radii Plus, Radii Xpert, SDI Ltd.).

g) Remove any excess cement using a scalpel or explorer.
Note: Excess cement can also be removed using self-curing gel stage (when cement feels rubbery), after 1 min - 30 seconds.

f) Maintain isolation for at least 5 min after seating of restoration.

STORAGE AND HANDLING:
- Avoid excess heat and direct sunlight.
- Store between 4°C and 25°C (39° - 77°F).
- Use at room temperature between 20°C and 25°C (68° - 77°F).
- If used at cooler temperature allow product to reach room temperature prior to use.
- Do not remove Riva Cem Automix syringes from foil packets until ready to use.

- For best results, use a new mixing tip. To store and avoid cross-contamination, discard the mixing tip and replace with original sealing cap securely attached to the syringe. Do not allow the pastes to dry out.

PRECAUTIONS:
- For professional use only.
- Do not take internally.
- Do not use after expiry date.

- Caution: Federal Law restricts this device to sale by or on the order of a dentist.
- SDB available at www.sdi.com.au or contact your regional representative.

FIRST AID:
- Eye (kontakt): Wash thoroughly with water and seek medical attention.
- Skin (kontakt): Remove using a cloth or sponge soaked in alcohol.
- Ingestion: Drink plenty of water/milk. Seek medical attention if symptoms persist.

Inhalation: No symptoms expected.

SHADE:
- Light Yellow

Made in Australia by SDI Limited

3-15 Brundson Street,
Bayswater, Victoria 3153

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