

# ADVANCED SHADING SYSTEM IN LUNA 2, A NEW UNIVERSAL COMPOSITE

## OBJECTIVES


To quantify and analyse the CIELAB colour parameters of new universal composite restorative Luna 2 (SDI Limited, Bayswater, Australia), and to demonstrate the relationship between shades in the VITA Classical range for ease of clinical use.

## EXPERIMENTAL METHODS

Universal composite Luna 2 was analysed in VITA Classical shades A1-A4, B1-B3, C1-C3, D2-D3. For colour analysis, discs of each shade of Luna 2 composite were prepared according to ISO-4049:2019 using a Radium Xpert curing light (SDI Limited, 440-480nm). Discs were analysed within the CIELAB colour space using an X-Rite SP-64 spectrophotometer. Additionally, Luna 2 sample kits were distributed to clinicians worldwide for in-vivo clinical evaluation and asked to respond to a survey to assess clinical performance of the restorative, including ease and accuracy of shade matching (n=28).





Figure 1: VITA Classical shade range and corresponding Luna 2 shades



## RESULTS AND DISCUSSION

Luna 2 composite discs were observed to have a close visual shade match to the VITA Classical shade guide (Figure 1).



The 1976 CIELAB system is one of the most widely used colour spaces<sup>1</sup> and represents colour within a three-dimensional space (Figure 2). In this space, data points are described by coordinates in the L\*, a\* and b\* axes. The L\* axis describes the lightness, or value, of the point, while the a\* axis represents a scale of green (-) to red (+), and the b\* axis represents a scale of blue (-) to yellow (+). The colour saturation, or chroma, is calculated according to the equation<sup>2</sup>:

$$C^*_{ab} = (a^{*2} + b^{*2})^{1/2}$$




Figure 2: CIELAB colour space, showing value and chroma of an L\*a\*b\* data point

Trends were observed within Luna 2 VITA Classical groups A, B, C, and D. Within a single shade group, value showed a consistent decrease across shades (Figure 3). Meanwhile, chroma increased across the group. This system demonstrates a simple and predictable relationship between shades within the VITA Classical range.



Figure 5: Clinical case restored using Luna 2 (images courtesy of Dr. Giovane Neri - Brazil)

As part of an in-vivo clinical evaluation of the product, 100% of the 28 surveyed clinicians responded with "Agree" or "Strongly Agree" to the statement "Selecting the shade was easy and shade matching was accurate" (Figure 4). The predictability of the Luna 2 advanced shading system allows for seamless shade matching and highly aesthetic clinical outcomes, as shown in Figure 5.

## CONCLUSION

The CIELAB colour parameters of shades in the Luna 2 VITA Classical shade range demonstrate a linear progression across shade groups, with decreasing lightness and increasing colour saturation (chroma). The predictability of this system simplifies the shade selection whilst ensuring a highly aesthetic result, supported in clinical use by 100% approval from surveyed dentists.

## REFERENCES

- <sup>1</sup>Durmus D. CIELAB color space boundaries under theoretical spectra and 99 test color samples. *Color Res Appl.* 2020; 45(5):796-802.
- <sup>2</sup>Hajira N., Mehta D., Ashwini P., Meena N., Usha HL. Influence of different enamel shades and thickness on chroma and value of dentin VITA shade: an in vitro comparative assessment study. *J Contemp Dent Pract.* 2015; 16(4):304-309.