

# ORDER DETAILS



## GS-80 CAPSULES


<b>Boxes of 50 capsules - One Spill 400mg alloy</b>	
Regular Set	4401403
Fast Set	4401402
<b>Boxes of 50 capsules - Two Spill 600mg alloy</b>	
Slow Set	4402404
Regular Set	4402303
Fast Set	4402202
<b>Boxes of 50 capsules - Three Spill 800mg alloy</b>	
Regular Set	4403303
Fast Set	4403202
<b>Boxes of 50 capsules - Five Spill 1200mg alloy</b>	
Regular Set	4405000
<b>Jars of 500 capsules - One Spill 400mg alloy</b>	
Slow Set	4421404
Regular Set	4421403
Fast Set	4421402
<b>Jars of 500 capsules - Two Spill 600mg alloy</b>	
Regular Set	4422303
Fast Set	4422202
<b>Jars of 500 capsules - Three Spill 800mg alloy</b>	
Regular Set	4423303
Fast Set	4423202
<b>Jars of 500 capsules - Five Spill 1200mg alloy</b>	
Regular Set	4425000

## GS-80 ALLOY TABLETS

Regular set 5oz	4450303
Regular set 1oz	4450505
Regular set 150g	4440303
Regular set 30g	4440304

## GS-80 ALLOY POWDER

Regular set 5oz	4460303
Regular set 30g	4470303
Regular set 50g	4480303
Regular set 250g	4960303



## EQUIPMENT

<b>Ultramat 2</b> (USA - includes power cord for USA)	5546045
<b>Ultramat 2</b> (All countries - includes power cord for the region the product is shipped to)	5546052
<b>Ultramat S</b> (All countries - includes power cord for the region the product is shipped to)	5546058

## MIXING TIMES

<b>Ultramat S</b>	
Frequency	Fixed (4600 rpm)
1, 2 & 3 Spill	8 ± 1 sec
5 Spill	6 ± 1 sec
<b>Ultramat 2</b>	
Frequency	Fixed (4550 rpm)
1, 2 & 3 Spill	8 ± 1 sec
5 Spill	6 ± 1 sec

GS-80 capsules conform to ISO 24234, ANSI/ADA Specification No.1.

\*Amalcap Plus, Contour, Stratosphere and Premium are not the registered trademarks of SDI Limited.

[1] Source-Published and SDI test data.



MADE IN AUSTRALIA  
by SDI Limited  
3-15 Brunson Street,  
Bayswater, Victoria 3153  
[www.sdi.com.au](http://www.sdi.com.au)

**AUSTRALIA** 1800 337 003  
**AUSTRIA** 00800 022 55 734  
**BRAZIL** 0800 770 1735  
**FRANCE** 00800 022 55 734  
**GERMANY** 0800 100 5759

**ITALY** 00800 0225 5734  
**NEW ZEALAND** 0800 808 855  
**SPAIN** 00800 022 55 734  
**UNITED KINGDOM** 00800 022 55 734  
**USA & CANADA** 1 800 228 5166

# SDI | GS-80

ECONOMICAL, HIGH STRENGTH,  
ADMIX AMALGAM



# HIGH STRENGTH NON-GAMMA 2 ADMIX ALLOY

The technically advanced, non-gamma 2 admix amalgam at an affordable price.

GS-80's high strength and positive dimensional change ensures longevity and patient satisfaction.

Over 200 million successful GS-80 restorations have been placed since its worldwide release in 1986.

**CHOICE OF CAPSULE SIZES AND SETTING TIMES**

**HIGH COMPRESSIVE STRENGTH**

## OPTIMAL HANDLING CHARACTERISTICS

GS-80's consistent mix, ease of loading, placement in preparations, condensability, burnishability, quality of interproximal contact, carvability, and setting characteristics ensure the final restorations meet the needs of the operator.




**POSITIVE DIMENSIONAL CHANGE**

# HIGH STRENGTH AND POSITIVE DIMENSIONAL CHANGE

## HIGH COMPRESSIVE STRENGTH


GS-80's high compressive strength minimizes the risk of fracturing the restoration<sup>[1]</sup>. High compressive strength ensures a tough and durable long lasting restoration.

### COMPRESSIVE STRENGTH @ 1 HOUR<sup>[1]</sup>




## POSITIVE DIMENSIONAL CHANGE

GS-80's positive dimensional change ensures excellent sealing qualities at the margins<sup>[1]</sup>. Negative values can contribute to post operative sensitivity.



## GREATER CHOICE OF SETTING TIMES AND LARGER CHOICE OF CAPSULE SIZES

The working and setting times allow large restorations to be condensed and carved without difficulty.




## PHYSICAL PROPERTIES

Composition:	Ag 40%, Sn 31.3%, Cu 28.7%, Hg 47.9%
Alloy Particle:	Spherical and Lathe Cut
Compressive strength @ 1 hour:	225MPa [32,625psi]
Compressive strength @ 24 hours:	510MPa [73,950psi]
Diametral tensile strength @ 1 hour:	22MPa [3,190psi]
Diametral tensile strength @ 24 hours:	54MPa [7,830psi]
Static creep @ 7 days:	0.1%
Dimensional change @ 24 hours:	+2µm/cm

## CAPSULE SYSTEM

The SDI capsule system is designed to give consistent and accurate mixes every mix.



## CONVENTIONAL CAPSULES

The plunger denotes the spill size:



The capsule body denotes the setting speed:



**NOTE:** 5 Spill only available in regular set.

## ULTRAMAT 2

The dual voltage Ultramat 2 with a high precision microprocessor controlled timer provides consistent and accurate trituration, ensuring the optimal handling and clinical characteristics of mixed materials.

