


# INSTRUCTIONS FOR CAPSULES^

- 1** Isolate and prepare tooth core.
- 2** Ensure the internal surface of the restoration to be seated is cleaned and dried. Do not dehydrate the tooth. Do not remove the smear layer. In cases where pulp protection is necessary, use a calcium hydroxide liner.
- 3** Activate the capsule and immediately mix in a triturator. Do not click before you mix.  
Press capsule down to click.  
Ten seconds in a triturator.
- 4** Immediately place into capsule applicator and click trigger until paste is seen through the clear nozzle.
- 5** Gently extrude approximately 1mm of cement onto the bonding surface.
- 6** Seat restoration within 30 seconds.
- 7** Remove excess cement at the first formation of the gel stage (approx. 2 min from start of mixing).
- 8** Maintain isolation until set is verified.


^ For powder and liquid instructions, please refer to Instructions for Use

## ORDER DETAILS



### CAPSULES, POWDER AND LIQUID

Riva Luting Capsules 50 x Riva Luting Capsules	8650008
Riva Luting Powder/Liquid Kit 35g Riva Luting Powder 25g Riva Luting Liquid Accessories	8650508
Riva Luting Powder/Liquid Triple Pack 3 x 35g Riva Luting Powder 3 x 25g Riva Luting Liquid Accessories	8650510



### APPLICATORS

Riva Applicator	5545009
Riva Applicator 2	5545013



\* Not the registered trademarks of SDI Limited.

\*\* Published and SDI Test Data.

[1] McCabe JF, Al-Naimi OT. Fluoride release into water for the Riva GIC products compared with competitor products. Experiment 2. University of Newcastle (UK); September 2004.

[2] Waldo B, Zhang P, Bennett J, Michalek SM, Katz J, Broome JC. Antimicrobial Activity of Composite-Resin and Glass-Ionomer Cements. University of Alabama at Birmingham, School of Dentistry, USA; March 2005.



06-2023

M870061 B



MADE IN AUSTRALIA  
by SDI Limited  
Bayswater, Victoria 3153  
Australia 1800 337 003  
[www.sdi.com.au](http://www.sdi.com.au)

**AUSTRIA** 00800 0225 5734  
**BRAZIL** 0800 770 1735  
**FRANCE** 00800 0225 5734  
**GERMANY** 0800 100 5759  
**ITALY** 00800 0225 5734

**NEW ZEALAND** 0508 802 557  
**SPAIN** 00800 0225 5734  
**UNITED KINGDOM** 00800 0225 5734  
**USA & CANADA** 1 800 228 5166

# **SDI** | **RIVA LUTING**

GLASS IONOMER LUTING CEMENT

CHEMICALLY  
ADHERES  
TO TOOTH  
STRUCTURE  
AND METAL  
SUBSTRATES



# SELF-CURING GLASS IONOMER LUTING CEMENT

Riva Luting is a conventional, self curing, glass ionomer luting cement designed for final cementation of metal based restorations. It chemically bonds to metal substrates and the tooth.

Riva Luting utilizes SDI's proprietary *ionglass™* filler developed by our glass technologists. *ionglass™* is a radiopaque, high ion releasing, reactive glass used in SDI's range of dental cements. Riva Luting releases substantially higher fluoride to assist with remineralization of the natural dentition.



## 17μM

LOW FILM THICKNESS

## 1


SIMPLE  
COLOR



## LOW SOLUBILITY

Riva Luting has low solubility in the oral environment. This increases the material's ability to resist disintegration and wear caused by oral acidity.

### WATER SOLUBILITY (7 MINUTES AFTER MIXING)\*\*




## SETTING TIMES

There will be minimal loss of working time when using Riva Luting. With adequate working time and a quick setting property, Riva Luting will shorten chair time and ensures early resistance.



## HIGH FLUORIDE RELEASE

Riva Luting's free movement of fluoride provides benefits to the tooth. Fluoride plays several significant roles in any caries - prevention program. These include the formation of fluorapatite, which is more acid resistant than hydroxyapatite. Fluorapatite is very important for tooth remineralization.



## ADVANTAGES

Chemically adheres to tooth structure and metal substrates

BPA & HEMA free

Very low film thickness - flows easily between restoration and tooth

Easy to clean up

Low water solubility  
- increased longevity and aesthetics

High bond strengths

High fluoride releasing and recharging ability

Bioactive proprietary *ionglass™* technology - high ion exchange

Great for PFM crowns, metal posts, orthodontic bands, implants

Light yellow shade

## INDICATIONS

Cementation of metal or porcelain fused to metal inlays, onlays and crowns

Cementation of stainless steel crowns

Cementation of posts and screws

Attachment of metal orthodontic bands

Lining

## PROVEN ANTIMICROBIAL EFFECT

Compared to other materials such as fluoride releasing resin cements, Riva Luting has proven antimicrobial activity against three cariogenic bacteria: *Streptococcus mutans*, *Streptococcus sobrinus* and *Lactobacillus* <sup>[2]</sup>.



Photos courtesy of Dr Shuichi Tsubura D.D.S, Tochigi, Japan