SDI LUNA

ENG	UNIVERSAL COMPOSITE
BUL	УНИВЕРСАЛЕН КОМПОЗИТ
CES	UNIVERZÁLNÍ KOMPOZIT
DAN	UNIVERSAL KOMPOSIT
DEU	UNIVERSAL KOMPOSIT
ELL	ΚΑΘΟΛΙΚΟ ΣΥΝΘΕΤΟ
FIN	UNIVERSAALI YHDISTELMÄMUOVI
FRA	COMPOSITE UNIVERSEL
HUN	UNIVERZÁLIS KOMPOZIT
ITA	COMPOSITO UNIVERSALE
LIT	UNIVERSALUS KOMPOZITAS
NOR	UNIVERSAL KOMPOSIT
POL	UNIVERSALNY KOMPOZYT
POR	RESINA COMPOSTA UNIVERSAL
RON	MATERIAL KOMPOZIT UNIVERSAL
SLK	UNIVERZÁLNÝ KOMPOZIT
SPA	COMPOSITE UNIVERSAL
SWE	UNIVERSAL KOMPOSIT

Made in Australia by SDI Limited
3-15 Brunson Street, Bayswater, Victoria 3153, Australia

Australia 1 800 337 003
Austria 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Italy 00800 0225 5734
Spain 00800 0225 5734
UK 00800 0225 5734
USA & Canada 1 800 228 5166

Last Revised: 08-2025

ENG INSTRUCTIONS FOR USE

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With non-stick handling, the radiopaque light cured, condensable and low-shrinkage composite, LUNA offers polishability and strength for anterior and posterior restorations alike.

COMPOSITION:

24% wt. (42% vol.) multifunctional methacrylic esters
76% wt. (58% vol.) inorganic fillers (40nm - 1.5µm)

INDICATIONS:

Direct anterior/posterior restorations and veneers. Core build ups. Indirect inlays, onlays and veneers. Splinting. Composite and porcelain repair. Sandwich technique with glass ionomers.

CONTRAINDICATIONS:

Pulp capping. Product may cause skin irritations to some people. In such cases, discontinue use and seek medical attention. Do not use in conjunction with any eugenol-containing materials.

INSTRUCTIONS FOR USE:

- Clean and isolate tooth. Rubber dam is the preferred mode of isolation.
- Prepare a conservative cavity, employing standard techniques.
- Use a calcium hydroxide liner and appropriate glass ionomer cement liner to protect deep cavity preparations in close proximity to the pulp.
- Select composite shade by matching with a clean and moist tooth using the SDI composite shade guide.
- Apply bonding system according to manufacturer's instructions or apply ZipBond Universal adhesive (SDI, Limited) or Stae single component dentin/enamel adhesive (SDI Limited).
- Place composite in increments of 2mm or less. NOTE: If using the syringe system, after each use, turn the screw plunger anti-clockwise half a turn to release any built up pressure.
- Light cure each increment for at least 10 seconds for light shades and at least 20 seconds for dark and/or opaque shades using a high powered LED light such as the Radii Plus (SDI Limited).
- Shape, polish and finish employing standard techniques.

NOTE: Instruct patient not to eat or drink anything for at least one hour after the procedure.

STORAGE AND HANDLING:

- Store at temperatures between 10° and 25°C (50°-77°F).
- Use at room temperature between 20° and 25°C (68°- 77°F).
- Syringe system: when not in use replace the cap tightly.

PRECAUTIONS:

- Avoid prolonged contact of etchant, adhesive and composite with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization. For professional use only.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiry date.
- Do not expose material to direct light.
- Caution: Federal Law restricts this device to sale by or on the order of a dentist.
- SDS available at www.sdi.com.au or contact your regional representative.

FIRST AID:

- Eye (contact):** Wash thoroughly with water and seek medical attention.
- Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- Ingestion:** Drink plenty of water/milk. Seek medical attention if symptoms persist.
- Inhalation:** No symptoms expected.

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