# SDBI polaoffice

## ADVANCED TOOTH WHITENING SYSTEM

### INSTRUCTIONS FOR USE

ENGLISH

Pola Office is a hydrogen peroxide based in-office tooth whitening system requiring minimal chair time. It is a neutral pH gel and contains desensitizers to maximize patient comfort.

#### INDICATIONS:

- Whitening of discoloured vital teeth.
- Whitening of discoloured non-vital teeth.

#### COMPOSITION:

Pola Office Liquid 35% Hydrogen peroxide  
65% Water

Pola Office Powder 73.26% Thickeners  
26.2% Catalysts  
0.04% Dye  
0.5% Desensitizing agents

Gingival Barrier 83.95% Methacrylic ester  
16% Silica  
0.04% Pigment  
0.01% Butylated hydroxy toluene

#### NET WEIGHT:

Depending on kit purchased. Please refer to outside of box:  
Pola Office: 1 syringe and 1 pot: 0.3g Bleach powder

2.25g Hydrogen peroxide liquid  
1g Gingival Barrier

Pola Office Bulk Kit: 3g Bleach powder  
22.5g Hydrogen peroxide liquid  
3g Gingival Barrier

#### PRECAUTIONS:

- For professional use only.
- Patient must wear protective eyewear.
- Operators must wear gloves, mask and protective eyewear.
- Do not use on pregnant or lactating women.
- Do not use on children under 14 years of age.
- Do not use the Gingival Barrier on any persons having known resin allergies.
- Do not use Pola Office on any persons having known peroxides allergies.
- Anyone with a history of chemical allergies is advised to undergo allergy checks before treatment.
- Do NOT anesthetize patients.
- Pola Office Liquid syringe may be under pressure – remove cap carefully.
- Do not allow mixed gel or Pola Office Liquid to come in contact with skin, eyes and soft tissues.
- Pola Office will not lighten any restorative materials.
- Do not use on patients with extremely sensitive teeth.
- SDS available at [www.sdbi.com.au](http://www.sdbi.com.au) or contact your regional representative.

#### FIRST AID:

- SKIN/TISSUES (contact):** Immediately wash thoroughly with water. Apply a slurry of sodium bicarbonate and water to the affected area. Seek medical attention if symptom persists.
- EYES (contact):** Open the eye wide and thoroughly wash for 15 minutes with running water.
- INGESTION:** Rinse mouth with water. Gargle with salt water and drink lots of milk. Seek medical attention if symptoms persist.
- INALHALATION:** no symptoms expected.

#### STORAGE AND HANDLING:

- Store in cool place (2°-25°C / 35°-77°F).
- Do not use after expiration date.

#### TIPS:

- Select appropriate candidate. Note: Not everyone is suitable for in-office whitening e.g. patient with translucent teeth or in the Vita\* C group.
- Patient must have good overall health.
- Seal all exposed roots and restorations.
- Inform patient that results are not guaranteed as the removal of stains varies between patients.
- Acid etching before treatment is NOT recommended.
- Large diameter suction tips are not recommended as they can spread the gel to other areas.
- Mild discomfort in teeth during treatment can be a common side effect.
- If patient experiences sensitivity during light curing, move the tip slightly away from the tooth. If this doesn't work, remove gel and rinse the affected area. Assess other possible causes.
- Delay any esthetic restorations for 2 to 4 weeks after treatment.
- Advise patient not to smoke, consume highly coloured, hot or cold foods or beverages for at least 2 days.
- Treat any post-operative sensitivity with a desensitizing gel or toothpaste.
- Prescribe Pola Day or Pola Night to maintain brightness if desired.
- A patient consent form is available for download at [www.sdbi.com.au](http://www.sdbi.com.au)

#### INSTRUCTIONS FOR VITAL TOOTH WHITENING:

- Determine and record pre-operative shade. Please follow this Vita\* shade guide arrangement according to degree of brightness: **B1•A1•B2•D2•A2•C1•C2•D4•A3•D3•B3•A3.5•B4•C3•A4•C4**
- Clean teeth with ONLY a flour based pumice.
- Place cheek retractors and then cover exposed lip surface with petroleum gel
- Dry teeth and apply Gingival Barrier to both arches, slightly overlapping enamel and interproximal spaces.
- Light cure in a fanning motion for 10-20 seconds until Gingival Barrier is cured.
- Open powder pot. Take one Pola Office syringe, **firmly** attach a tip, and **carefully pull back** plunger to release pressure.
- Carefully** extrude contents of syringe into the pot.
- Immediately mix using a brush applicator until gel is homogeneous.
- Apply a **thick** layer of gel to all teeth undergoing treatment.
- Leave gel on for 8 minutes (Optional: a curing light may be used at this stage, please refer to light manufacturer's recommendations)
- Suction off any surplus aspirator tip
- Complete Steps 9-11 three times (optional 4 times). Note: One pot should only be used for 2 applications or within 20 minutes of mixing.
- After the last application, suction all the gel off, then wash and apply suction.
- Remove Gingival Barrier by lifting it from an end.

#### TIPS FOR NON-VITAL TOOTH WHITENING:

- Isolate tooth effectively.
- Protect oral mucosa to prevent chemical burns.
- Verify adequate endodontic obturation.
- Use protective barrier to prevent peroxide from reaching the periodontal ligament.
- Avoid acid etching as periodontal ligament irritation may result.
- Advise the patient that a whitened non-vital tooth will not look the same as a vital tooth.
- Recall patients periodically to examine the whitened teeth (clinically & radiographically).

#### INSTRUCTIONS FOR NON-VITAL TOOTH WHITENING:

- Determine and record pre-operative shade. Please follow this Vita\* shade guide arrangement according to degree of brightness: **B1•A1•B2•D2•A2•C1•C2•D4•A3•D3•B3•A3.5•B4•C3•A4•C4**
- Isolate tooth/teeth undergoing whitening using a rubber dam.
- Isolate the dam fits tightly at the cervical margin. Note: If a rubber dam is not used, cover approximately half of the adjacent vital teeth with 3-4mm thick layer of Gingival Barrier to insulate from heat discomfort of the curing light. This is done when one tooth is being whitened non-vitaly.
- Open lingual access to pulp chamber and remove all composite, base material and gutta cerca 2-3 mm apical to the CE junction. Place a glass ionomer/zinc phosphate base, at least 1 mm thick to ensure a barrier between the "sealed" root canal and whitening gel.
- Place and light cure a thick layer of Gingival Barrier at the gingiva region of the isolated tooth.
- Mix the Pola Office gel as instructed in steps 6-8 of the vital tooth bleaching instructions.
- Place the Pola Office gel into the opened pulp chamber and onto the labial surface.
- Leave gel on/in for 8 min (Optional: curing light may be used at this stage, please refer to light manufacturer's recommendations)
- NOTE: If using a curing light and treating inner and outer surfaces, 2 curing lights can be used (one on each side of tooth) to activate Pola Office.
- Suction off the gel.
- Repeat steps 6-8 until desired shade is achieved. Note: One pot should only be used for 2 applications or within 20 minutes of mixing.
- After the last application, suction and wash the gel off.
- Remove Gingival Barrier and re-seal tooth with appropriate restorative.
- Record the final shade.

\*Vita shade guide is a trademark of VITA Zahnfabrik, H. Rauter GmbH & Co. KG, D-79713 Säckingen, Germany

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