# SDI riva protect

GLASS-IONOMER FISSURE AND TOOTH PROTECTOR INSTRUCTIONS FOR USE

ENGLISH

Riva Protect is a high-fluoride-release glass-ionomer cement, designed for surface protection and sealing of newly-erupted or caries-challenged dentition. Its viscosity and handling is designed for ease of use and placement, and its direct adhesion ensures long-term benefit. Riva Protect is available in white and pink shades. The pink shade is darker (L\* =65) for easy identification.

**INDICATIONS FOR USE:**

- Pit and fissure sealing.
- Root /tooth surface protection and hypersensitivity prevention.
- Temporary restorations in non-load-bearing applications.
- Temporary endodontic sealing.

**CONTRAINDICATIONS:**

- Pulp capping.
- Product may cause skin irritations to some people. In such cases, discontinue use and seek medical attention.

**INSTRUCTIONS FOR USE:**

**PREPARATION:** Isolate and clean tooth. Lightly air dry without dehydrating tooth surface.

<b>ACTIVATION OF CAPSULE:</b>		
Setting Plunger colour	REGULAR (White Shade: White) (Pink Shade: Pink)	FAST (White Shade: Blue) (Pink Shade: Orange)
Powder/Liquid ratio (g/g)	0.34/0.19	0.34/0.19
Mixing time (min.sec.)	10"	10"
Working time (min.sec.)	1'40"	1'05"
Initial setting time (min.sec.)	4'10"	2'00"
Finishing time @ 37°C from start of mixing (min.sec.)	5'00"	3'00"
Minimum deliverable volume per capsule: 0.13mL		
Test conditions: temperature (23 ± 1°C), relative humidity (50 ± 10%). ISO 9917:2003 (Dental water based cements) Note: Higher temperatures will shorten setting time.		

- Push the plunger firmly by hand or on a hard surface, until it is flush with the body of the capsule.

- Immediately place the capsule into the Ultramat 2 Amalgamator, or any other high-speed mixer (4000-4800rpm), and triturate for 10 seconds. Once capsule trituration is complete, immediately remove the capsule and place into the Riva applicator. NOTE: Prior to placement, click the trigger of the Riva applicator until glass ionomer paste is seen through the clear nozzle.

<b>or</b>	
<b>POWDER-LIQUID DISPENSING AND MIXING:</b>	
Powder-liquid dosage:	2 flat scoop to 3 drops
Powder/Liquid ratio (g/g)	2:1
Recommended mixing time (sec.)	30"
Working time (min.sec.)	2'00" (approximate)
Finishing time @ 37°C from start of mixing (min.sec.)	5'00" (approximate)
Mixing pad type: non-absorbent paper or glass slab	
Note: Use a cool glass slab for extended working time	
Test conditions: temperature (23 ± 1°C), relative humidity (50 ± 10%). ISO 9917:2003 (Dental water based cements)	

- Gently tap powder jar on the hand. Using the plastic scoop provided, dispense two level measures of powder into the mixing pad.
- Carefully dispense three drops of liquid next to the powder on the mixing pad. Replace cap tightly.
- Use a plastic spatula to divide the powder into two equal parts.
- Mix the liquid with one part of the powder for 10 seconds and then include the second part and continue mixing for another 15 to 20 seconds.

**PIT AND FISSURE SEALING AND ROOT/TOOTH SURFACE PROTECTION**

- Apply a polyacrylic acid conditioner such as Riva Conditioner (SDI Ltd.) to the tooth surface for 10 seconds.
- Wash and air-dry. Leave surface moist.
- Place Riva Protect directly onto the surface of the tooth.
- Spread the material onto the tooth surface or into the pits and fissures using a small brush eg POINTS (SDI Ltd.).
- When the material has lost its surface gloss, apply a thin film of light cured, unfilled resin such as Riva Coat (SDI Ltd.).
- Light cure for 20 seconds.
- Instruct the patient not to eat for at least one hour.

**TEMPORARY RESTORATION**

- Apply a polyacrylic acid conditioner such as Riva Conditioner (SDI Ltd.) to the cavity for 10 seconds.
- Wash and air-dry. Leave surface moist.
- Place Riva Protect into the cavity, taking care to avoid air entrapment under the restoration.
- Form contour.
- When the material has lost its surface gloss, apply a thin film of light cured, unfilled resin such as Riva Coat (SDI Ltd.).
- Light cure for 20 seconds.
- Instruct the patient not to eat for at least one hour.

**TEMPORARY ENDODONTIC SEALING**

- Prepare and treat pulp chamber accordingly.
- Clean and gently dry the pulp chamber.
- Fill the pulp chamber with a cotton bud.
- Place Riva Protect over the cotton bud.
- When the material has lost its surface gloss, apply a thin film of light cured, unfilled resin such as Riva Coat (SDI Ltd.).
- Light cure for 20 seconds.

**STORAGE AND HANDLING**

- Do not remove Riva Protect capsules from their foil packets until ready to use.
- Store at temperatures between 4°C - 20°C (39°F - 68°F) for capsules and 4°C - 30°C (39°F - 86°F) for powder/liquid kits.
- Use at room temperature between 20°C and 25°C (68°F - 77°F).

**PRECAUTIONS:**

- For professional use only.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiry date.
- Do not mix powder or liquid with any other brand of glass ionomer products.
- Caution: Federal Law restricts this device to sale by or on the order of a dentist.
- MSDS available at www.sdi.com.au or contact your regional representative.

**FIRST AID**

- Eye (contact):** Wash thoroughly with water and seek medical attention.
- Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- ingestion:** Drink plenty of water/milk. Seek medical attention if symptoms persist.
- Inhalation:** (powder) Remove to fresh air. If breathing is difficult, seek medical attention.

**PACKAGING:**

- 50 Capsule Boxes (Pink or White) (Regular or Fast set)
- Powder/Liquid Kits (Pink or White): 10g (9.1mL) Liquid Bottle, 15g Powder Jar, 10mL Riva Conditioner bottle, accessories
- 10g (9.1mL) Liquid Bottle Refills
- 15g Powder Jar Refills (Pink or White)

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