# STELA CAPSULE STELA AUTOMIX STELA PRIMER

ENG

## Bulk Fill Flowable Self Cure Restorative

### INSTRUCTIONS FOR USE

Stela is an easy to use bulk-fill radiopaque self-cure restorative system. The combination of Stela Primer with Stela Capsule or Stela Automix will bond to both enamel and dentine. No scrubbing nor light-curing are required for any parts of the Stela restorative system.

#### COMPOSITION OF STELA CAPSULE:

76.8 wt% (55.4 vol%) inorganic filler

#### COMPOSITION STELA AUTOMIX:

61.2 wt% (36.4 vol%) inorganic filler

#### INDICATIONS FOR USE (Stela Primer)

- Dentin and enamel bonding

#### INDICATIONS FOR USE

- Direct restoration of Class I, II, III and V cavities (ideally Class I and II)
- Base or liner
- Core build-ups

#### CONTRA-INDICATIONS:

- **Do not** use Stela Primer with other composite products as it will not cure
- **Do not** use if a dry field cannot be established
- **Do not** use for pulp capping
- **Do not** use in conjunction with any eugenol containing materials
- **Do not** use in anyone who has an allergy to acrylics or who has an allergic reactions to the material of this product.

**NOTE:** This product may cause skin irritation to some people. In such cases, discontinue use and seek medical attention.

#### PROCEDURE

##### A. CAVITY PREPARATION

1. Prepare a conservative cavity with round internal angles, employing standard techniques and instruments. Undercuts for retention are not necessary.
2. Clean and isolate tooth. Rubber dam is the preferred mode of isolation. Pumice and water or chlorhexidine solutions are the only substances recommended to clean the cavity. Residues of other substances can impair the bond strength and setting reaction.
3. In cases where pulp protection is necessary, use a calcium hydroxide liner (pulp capping) and/or appropriate glass ionomer cement (liner/base) to protect deep cavity preparations, before applying Stela Primer.
4. For Class II cavities or whenever the clinical situation requires, apply a thin, cambered, well-adapted matrix band as needed. It is recommended to fixate the matrix approximately with interdental wedges. For Classes III cavity, a plastic matrix is preferred.

##### B. STELA PRIMER APPLICATION

1. Dispense one or two drops (depending on restoration size) of Stela Primer into a mixing well. Immediately replace the cap.
2. Using a disposable brush applicator with solvent-resistant fibres (Points, SDI Limited), apply Stela Primer onto prepared cavity surfaces and leave on cavity for 5 seconds. Scrubbing is **NOT** necessary.
3. Gently blow with air (free of water/oil) until you see no movement of the primer (2-3 seconds).

**NOTE:** As Stela Primer is a thin liquid, a hard blow is NOT required. A vacuum aspirator can be used to prevent primer liquid from scattering.

**NOTE:** Do **NOT** light cure Stela Primer. It is designed to completely cure when in contact with Stela Capsule or Stela Automix paste.

##### C. MIXING AND PLACING STELA CAPSULE OR STELA AUTOMIX

The Stela restorative is available in Capsules (single dose) or Automix syringes (multiple doses). The following instruction cover both Stela presentations.

	CAPSULE	AUTOMIX
Mixing time*	10 sec	Not applicable (mixing tip)
Working time*	1 min 30 sec	1 min 30 sec
Finishing time at 37°C/99°F from start of mixing**	4 min	4 min

Type 1, Class I, polymer based restorative material according to ISO 4049

#### NOTE:

\*Test conditions: temperature (23 ± 1°C, 73 ± 2°F), relative humidity (50 ± 10%). Temperatures higher than 23 ± 1°C and 73 ± 2° F will shorten the working time.

\*\*Setting time at 37°C/99°F (ISO 4049): 1 min 30 sec. Refer to Section D for information on Finishing.

#### I. CAPSULE

1. Activate the capsule by pressing the plunger into a hard surface until it completely locks inside the body of the capsule. Important: Do not click with an applicator before you mix.
2. Immediately place activated capsule into the Ultramat amalgamator or any high-speed capsule mixer (4000 – 4800 rpm).
3. Mix for 10 seconds.
4. Immediately remove the capsule and place into a capsule applicator.
5. Click the trigger of the capsule applicator until paste is seen through the clear nozzle.

**NOTE:** Stela capsule may discharge some residual powder during the initial applicator clicks. This may produce airborne particles. Perform this action while wearing personal protective equipment (PPE) and away from patients.

6. Extrude paste into the cavity, being careful not to trap air under the restoration. Fill the entire cavity in a single step. Slightly overfill to ensure good contact with the primer at the margin.

**NOTE:** To prevent the formation of an oxygen inhibition layer, a protective gel (e.g., glycerine) can be used after placing Stela paste.

#### II. AUTOMIX

1. Remove cap. For first use only, dispense 2-3mm of material onto a pad and discard (to ensure correct ratio of extrusion). Attach a new mixing tip to the Automix syringe.
2. Dispense a small amount until an evenly mixed paste flows out of the tip. Discard this material.
3. Extrude paste into the cavity, being careful not to trap air under the restoration. Fill the entire cavity in a single step. Slightly overfill to ensure good contact with the primer at the margin.

**NOTE:** To prevent the formation of an oxygen inhibition layer, a protective gel (e.g., glycerine) can be used after placing Stela paste.

**NOTE:** The mixing tip is for single use only. Discard the used mixing tip and place the original syringe cap onto the syringe. Syringe must be stored with the cap.

#### D. FINISHING

1. Finish the restoration using standard techniques, after at least 4 minutes from the start of the mixing time
2. Instruct patient not to eat or drink anything for at least one hour after the procedure.

#### STORAGE AND HANDLING

- Avoid excess heat
- Store at temperatures between 2° and 25°C (35° - 77°F)
- For optimum freshness, refrigeration is optional. Product temperature must be between 20° and 25°C (68° and 77°F) when used.

#### PRECAUTIONS

- Do not use Stela Primer with any other products, as it will not cure.
- Avoid prolonged contact of Stela products with the skin or oral tissue, as they may cause inflammation of the oral tissues or skin sensitization
- Stela Primer: contact with eyes may cause severe eye damage. Wear eye protection to prevent injury. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Stela Primer is a flammable liquid. Avoid inhaling the vapour.
- Do not use in anyone who has an allergy to acrylics or who has an allergic reactions to the material of this product.
- For professional use only
- Keep out of reach of children and pets
- Do not take internally
- Do not use after expiry date

**CAUTION:** Federal Law restricts this device to sale by or on the order of a dentist. SDS available at [www.sdi.com.au](http://www.sdi.com.au) or contact your regional representative.

#### FIRST AID

- **Eye (contact):** Wash thoroughly with water and seek medical attention.
- **Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- **Ingestion:** Rinse mouth with water and discard. Drink plenty of water/milk. Do NOT induce vomiting. Seek medical attention.
- **Inhalation:** No symptoms expected under normal conditions of use. If affected, move to well ventilated area.

Made in Australia by SDI Limited  
3-15 Brunson Street,  
Bayswater, Victoria 3153  
Australia

Australia 1 800 337 003  
Austria 00800 0225 5734  
Brazil 0800 770 1735  
France 00800 0225 5734  
Germany 0800 1005759  
Italy 00800 0225 5734  
Spain 00800 0225 5734  
United Kingdom 00800 0225 5734  
USA & Canada 1 800 228 5166  
[www.sdi.com.au](http://www.sdi.com.au)

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