SDI SUPER ETCH SUPER ETCH LV

E
N
G
L
I
S
H

INSTRUCTIONS FOR USE

Super Etch and Super Etch LV (low viscosity) 37% phosphoric acid etchant for etching enamel and dentin.

The Super Etch syringe facilitates controlled direct extrusion for faster procedures and easier access into difficult to reach areas. SDI's inorganic filler particles wash away without leaving a residue.

Super Etch LV is specialized for sealants. Super Etch LV's greater flowability is formulated for pit and fissure etching. Super Etch LV attains greater surface contact for deeper etchant penetration into the pits and fissures. The Super Etch LV tip's diameter opening of only 0.41mm, enhances penetration into the smallest pits and fissures; crevasses that will not admit the wider dispensing tips of other etchants.

Indications

1. Etching of enamel
2. Etching of dentine

Contraindications

1. Soft tissues.
2. If the treatment site cannot be etched.
3. Patients who are allergic to phosphoric acids and blue dyes.

INSTRUCTIONS:

1. Isolate tooth. Prepare cavity.
2. Clean the surface to be etched using a non-fluoridated, oil-free prophylaxis paste, or a pumice and water slurry.
3. Rinse thoroughly with water.
4. Thoroughly dry the surface to be etched with dry, oil-free air.
5. **CAUTION:** DO NOT dispense Super Etch LV without a dispensing tip on. Only use the pre-bent tips supplied. Before dispensing, slightly pull the plunger backward (~1mm) and then push forward.

Note: on first usage, or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the etchant's viscosity and rate of extrusion.

6. Leave on:
 - a. enamel for at least 30 seconds.
 - b. dentin and glass ionomer liner for 20 seconds
7. Rinse well with water for 5-10 seconds.
8. Remove excess water.
9. Apply bonding agent, such as **Stae**, or pit and fissure sealant, such as **Conseal F**, according to manufacturer's instructions.

WARNINGS & PRECAUTIONS:

- **CORROSIVE** - Contact can cause pain, redness, burns, and blistering. Avoid contact with oral tissues, eyes and skin.
- For professional use only.
- Keep out of reach of children.
- Do not take internally.
- When not in use replace the bottle cap tightly.
- Use at room temperature.
- Do not use after expiry date.
- When not in use keep the cap tightly sealed.
- Federal Law restricts this device to sale by or on the order of a dentist.
- SDS available at www.sdi.com.au or contact your regional representative.

FIRST AID:

- **Eye contact:** Hold eyelids apart and flush the eye continuously with running water. Continue flushing until advised to stop by a Poisons Information Centre (e.g. Australia 13 11 26; New Zealand 0800 764 766) or a doctor, or for at least 15 minutes.
- **Skin contact:** Remove contaminated clothing and flush skin and hair with running water.
- **Ingestion:** Do NOT induce vomiting. Rinse mouth. Drink plenty of water. Seek medical attention.
- **Inhalation:** No symptoms expected unless heated or misted.

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