SDI

aura

ULTRA UNIVERSAL
RESTORATIVE MATERIAL
INSTRUCTIONS FOR USE

PAK41333 A

SDI

aura

ULTRA UNIVERSAL
RESTORATIVE MATERIAL
INSTRUCTIONS FOR USE

INTRODUCTION

aura ultra universal restorative system is a the easy solution to creating as close to nature anterior and posterior restorations. Based on the natural make-up and colours of the dentine and enamel, **aura** uses a logical yet innovative approach to shade-matching and layering in one single composite system.

INDICATIONS

	Enamel Shade	Dentine Shade
Anterior restorations	•	•
Posterior restorations	*	•
Indirect inlays, onlays & veneers	*	•
Core build ups		•
Splinting	•	•
Sandwich technique with glass ionomer		•

*For use as thin overlay over Dentine Chroma shade.

CONTRA-INDICATIONS

- Pulp capping.
- Do not use in conjunction with any eugenol containing materials.
- Anyone with allergy to acrylics

SHADE RANGE

Enamel shades:

- **aura** enamel shades are microfilled composites based on different colours of natural enamels.
- E1 has translucency, high opalescence and is whiter. It is designed to emulate young enamel.
- E2 has intermediate translucency, is neutral in colour and less opalescence. It is designed to replace adult enamel.
- E3 has a yellow tint with high translucency and is suitable for replacing aged enamel.

Note: **aura** enamel shades are not radiopaque


Dentine shades:

- **aura** dentine shades are nanohybrid composites based on single opacity, single hue (base colour similar to the "A" series of the Vita Classical Shade Guide), fluorescent with a broad range of chroma. "DC" stands for "dentine chroma" and is equally spaced out in ascending chroma value. A bleached dentine, Db, is also available for bleached teeth.

SHADE SELECTION

Note: Shade selection must be carried out before isolation and ideally should be done under day light using the **aura** shade guide provided.

1. Select the dentine shade (Db / DC1 / DC2 / DC3 / DC4 / DC5 / DC6 / DC7) by matching the shade tab to the cervical area where the enamel has minimal thickness.
2. Select the enamel shade (E1 / E2 / E3) by matching the shade tab to the incisal of the neighbouring teeth or tooth to be restored.

PLACEMENT PROCEDURE:

1. Isolate tooth.
2. Prepare cavity using standard technique for composite resin restorations.
Note: In cases where pulp protection is necessary, use a calcium hydroxide liner.
3. Apply a dentin/enamel bonding agent such as SDI GO! as per manufacturer's instructions.
4. Place **aura** Dentine / Enamel shade in increments of no more than 2mm.
5. Light cure each increment using a high powered LED curing light (460 - 480nm wavelength) for:
- 10 seconds for E shades (Enamel Shades)
- 20 seconds for Db, DC (Dentine Shades)
6. Finish restoration using standard technique.

STORAGE AND HANDLING

- Store at temperatures between 10° and 25°C (50°-77°F).
- Use at room temperature between 20° and 25°C (68°- 77°F).
- Syringe system: when not in use replace the cap tightly.

PRECAUTIONS

- Avoid prolonged contact of **aura** with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization.
- For professional use only.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiry date.
- Do not expose material to direct light.
- Caution: Federal Law restricts this device to sale by or on the order of a dentist.
- SDS available at www.sdi.com.au or contact your regional representative.

FIRST AID

- **Eye (contact):** Wash thoroughly with water and seek medical attention.
- **Skin (contact):** Remove using a cloth or sponge soaked in alcohol.
- Wash thoroughly with water.
- **Ingestion:** Rinse mouth thoroughly. Drink plenty of water/milk. Seek medical attention if symptoms persist.
- **Inhalation:** No symptoms expected.

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