SDI | STAE

ENG	DENTIN / ENAMEL SINGLE COMPONENT ADHESIVE TOTAL ETCH ADHESIVE
DEU	EINKOMPONENTEN TOTAL ETCH ADHÄSIV-SYSTEM FÜR DENTIN/SCHMELZ
FRA	ADHÉSIF MONOCOMPOSANT À MORDANÇAGE DENTINE/EMAIL
ITA	ADESIVO SMALTO/DENTINA MONOCOMPONENTE TOTAL ETCH
KOR	상아질/법랑질용 1성분 토탈 에칭 접착제 시스템
POL	JEDNOSKŁADNIKOWY SYSTEM WIĄZĄCY ZE SZKLIWEM I ZEBINĄ
POR	ADESIVO SIMPLIFICADO PARA ESMALTE/ DENTINA TOTAL ETCH
SPA	ADHESIVO GRABADOR TOTAL DE UN SOLO COMPONENTE PARA LA DENTINA Y EL ESMALTE
TUR	DENTIN/ENAMEL İÇİN TEK BİLEŞENLİ TOPLAM ETCH ADEZİV SİSTEMİ

	Made in Australia by SDI Limited 3-15 Brunson Street, Bayswater, Victoria 3153, Australia		
	Australia 1 800 337 003	SDI Germany GmbH	
	Austria 00800 0225 5734	Harsestrasse 95	
	Brazil 0800 770 1735	51149 Cologne, Germany	
	France 00800 0225 5734	T: +49 2203 9255 0	
	Germany 0800 1005759		
	Italy 00800 0225 5734		
	Spain 00800 0225 5734		
	UK 00800 0225 5734		
	USA & Canada 1 800 228 5166	Last Revised: 08-2025	
	www.sdi.com.au		

PAK0372 L

ENG INSTRUCTIONS FOR USE

Stae is a fluoride releasing single component dentin / enamel total - etch adhesive system, designed for direct bonding. Stae completely and homogeneously infiltrates the hybrid layer to ensure superior bonding to the tooth.

COMPOSITION:

Acrylic monomer
Acetone
Fluoride
Stabilizer

INDICATIONS:

- I) All direct restorative composites
- II) Compomers
- III) Composite / Ceramic+ / Metal+ / Porcelain+ repairs.
 - * Bonding surface requires pre-treatment with a silane primer.

CONTRAINDICATIONS:

- I) Not to be used in combination with chemically-cured composite resin.
- II) Do not use for pulp capping.
- III) Any persons having known resin allergies should refrain from using Stae.
- IV) Not to be used in combination with eugenol materials or hemostatics containing ferric compounds.
- V) Not to be used in combination with desensitising agents.

INSTRUCTIONS FOR USE:

- Clean and isolate tooth. Rubber dam is the preferred mode of isolation.
- Prepare a conservative cavity, employing standard techniques.
- Use a calcium hydroxide liner and appropriate glass ionomer cement liner to protect deep cavity preparations in close proximity to the pulp.
- Thoroughly dry the bonding surfaces of the tooth with dry, oil-free air.
- Etch with Super Etch 37% phosphoric acid (SDI Limited) for at least 20 seconds. NOTE: Enamel subjected to fluoridation should be etched for 90 to 120 seconds. IMPORTANT: If using an etch syringe, ensure that the dispensing tip hub is **firmly** attached to the syringe by twisting the hub securely onto the syringe. On first usage or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the etchant’s viscosity and rate of extrusion.
- Wash thoroughly with water.
- Remove excess water. Keep moist. Avoid contamination e.g. saliva.
- Apply Stae to saturate all internal surfaces.
- Blow gently with dry, oil-free air to evaporate solvent. Leave surface glossy.
- Light cure for 10 seconds.

NOTE: do not polymerize Stae by itself in a mixing well; it will be a non – polymerized looking material. Stae contains a unique hydrophobic and hydrophilic monomer that does not cure when exposed to oxygen. The product will polymerize once covered by any lightcure composite restorative material.
- Place composite in increments of 2mm or less.

STORAGE AND HANDLING:

Store at temperatures between 10° and 25°C (50°-77°F). Use at room temperature between 20° and 25°C (68°- 77°F). Shelf life: 3 years

PRECAUTIONS:

Avoid prolonged contact of etchant, adhesive and composite with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization. For professional use only. Keep out of reach of children. Do not take internally. Do not use after expiry date. Do not expose material to direct light. Caution: Federal Law restricts this device to sale by or on the order of a dentist. SDS available at www.sdi.com.au or contact your regional representative.

FIRST AID:

Eye (contact): Wash thoroughly with water and seek medical attention.
Skin (contact): Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
Ingestion: Drink plenty of water/milk. Seek medical attention if symptoms persist.
Inhalation: No symptoms expected.

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