| | |
|---|--|
| SDI | RIVA SILVER |
| ENG | SILVER-ALLOY REINFORCED GLASS IONOMER RESTORATIVE MATERIAL |
| DAN | SOLVFORSTÆRKT GLASIONOMER FYLINGSMATERIALE |
| DEU | SILBERVERSTÄRKTER GLAS-IONOMER-FÜLLUNGSZEMENT |
| FRA | CEMENT VERRE IONOMÈRE D'ORTURATION RENFORCE A L'ARGENT MATERIALE VETRO-IONOMERICO PER RESTAURO RINFORZATO CON ARGENTO |
| ITA | 중합된 글라스 이온오머 복합 재료 |
| KOR | SDABRU SUSTIPRITAS STIKLO IONOMERIS UZPILDAS |
| LIT | GLASJONOMER D VYKELNIN WZMONONNY OPIKAM SREBRA |
| POL | MATERIAŁ DO RESTAUROACJI DE LIGA DE PRATA COM IONOMERO DE VIDRO |
| POR | IONOMERO DE VIDRO RESTAURADOR REFORZADO CON ALEACION DE PLATA |
| SPA | GUMIS LE GUCLENDIRILIMIS CAM IONOMER KADEM MADDESİ |
| TUR | |
| PAK411110N | |
|  | Made in Australia by SDI Limited 3-15 Brunston Street, Bayswater, Victoria 3153, Australia |
|  | Australia 1 800 337 003 Austria 0800 0225 5734 Brazil 0800 770 1735 France 0800 0225 5734 Germany 0800 1005759 Italy 0800 0225 5734 Spain 0800 0225 5734 UK 0800 0225 5734 USA & Canada 1 800 228 5166 |
|  | SDI Germany GmbH Hansesstraße 85 51149 Cologne, Germany T: +49 2203 9255 0 |
|  | CE |
|  | MD |
|  | PDR |
|  | 0297 |
|  | Rx |
|  | only |
|  | capsule |
|  | powder/liquid |
|  | 41°C 39°F |
|  | 20°C 68°F |
|  | 41°C 39°F |
|  | 30°C 86°F |

ENG INSTRUCTIONS FOR USE

Riva Silver is the ideal radiopaque, silver-alloy-reinforced glass ionomer restorative material. With excellent tooth adhesion, sustained fluoride release and high compressive strength, Riva Silver is useful for a wide range of restorative applications. Shade: light grey.

INDICATIONS

- Small class I and II restorations.
- Deciduous teeth restorations.
- Geriatric restorations.
- Intermediate restoratives and base material for class I and II cavities using the sandwich technique.
- Class V restorations in non-aesthetic situations.
- Core buildups.
- Temporary fillings.
- Restorations where radiopacity is desirable.
- Restorative in the field using the ART technique.
- Temporary endodontic fillings.

CONTRAINDICATIONS

- Pulp capping
- Product may cause skin irritations to some people. In such cases, discontinue use and seek medical attention.

PROCEDURE

A. CAVITY PREPARATION:

- Clean and isolate tooth.
 - Note: A conservative cavity should be prepared, employing standard techniques and instruments. Undercuts are not necessary.
- In cases where pulp protection is necessary, use a calcium hydroxide liner.

B. CONDITIONING / ETCHING:

- Apply Riva Conditioner to the prepared surfaces and leave in cavity for 10 seconds.

OR

- Apply Super Etch 37% Phosphoric Acid to the prepared surfaces and leave in cavity for 5 seconds.
 - Rinse thoroughly with water.
 - Remove excess water. Keep moist. Avoid contamination.
- IMPORTANT:** Be careful not to completely dry the surface. Surface should be glistening.

C. MIXING AND PLACING

| | |
|--|------------|
| POWDER / LIQUID SYSTEM | |
| Powder/Liquid ratio (g/g) | 0.72/0.13 |
| Mixing time (min. sec.) | 10 seconds |
| Working time (min. sec.) | 1'30" |
| Initial setting time (min. sec) | 5'30" |
| Finishing time @ 37°C from start of mixing (min. sec) | 6'00" |
| Minimum delivery/capsule | 0.15mL |
| Test conditions: temperature (23 ± 1°C), relative humidity (50 ± 10%), ISO 9917-1 (Non translucent radiopaque glass polyalkenoate cement). Minimum net contents per capsule: 0.72g powder and 0.14g (0.12mL) liquid. | |

- Activate the capsule by pushing the plunger until it is flush with the body.
 - Important: Do not click with applicator before you mix.**
- Immediately place the capsule into the Ultramat 2 amalgamator, or any other suitable mixer, and triturate for 10 seconds.
- Immediately remove the capsule and place into the Riva applicator.
- Click the trigger of the Riva applicator until glass ionomer paste is seen through the clear nozzle.
- Fill Riva Silver into the cavity, being careful not to trap air under the restoration. **Note:** At 23°C/74°F working time will be about 90 seconds (1'30") from the start of mixing. In situations above this temperature, working times will shorten. Adhesion strength will decrease if material is manipulated after this time.
- Immediately apply Riva Coat to all exposed surfaces of restoration.
- Light cure for 10 seconds.
- At least 6 minutes after the start of trituration, finish the restoration using standard techniques.
- Apply another layer of Riva Coat.
- Light cure for 10 seconds.
- Instruct patient not to eat for at least one hour after procedure.

| | |
|-------------------------------------|------------------------|
| POWDER / LIQUID SYSTEM | |
| MIXING OF LIQUID AND POWDER: | |
| Powder-liquid dosage: | 1 flat scoop to 1 drop |
| Powder-liquid ratio: | 0.300g/0.043g (7:1) |
| Recommended mixing time: | 30 secs. |

- Gently tap Riva Silver powder bottle on the hand.
- Dispense one level measure of powder onto the mixing pad/ glass slab. Replace cap tightly.
- Dispense one drop of Riva Silver onto the mixing pad/glass slab next to the powder. Replace cap tightly.
- Divide the powder into two equal parts.
- Mix the liquid with one part of the powder for 10 seconds using a plastic spatula, then include the second part and mix for another 15 to 20 seconds.

STORAGE AND HANDLING

- Do not remove Riva Silver capsules from their foil packets until ready to use.
- Store at temperatures between 4°C - 20°C (39°F - 68°F) for capsules and 4°C - 30°C (39°F - 86°F) for powder/liquid kits.
- Use at room temperature between 20°C and 25°C (68°F - 77°F).

PRECAUTIONS

- For professional use only.
- Keep out of reach of children.
- Do not use after expiry date.
- Do not mix powder or liquid with any other brand of glass ionomer products.
- Caution: Federal Law restricts this device to sale on or by order of a dentist.
- SDS available at www.sdi.com.au or contact your regional representative.

FIRST AID

Eye (contact): Wash thoroughly with water and seek medical attention.

Skin (contact): Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.

Ingestion: Drink plenty of water/milk. Seek medical attention if symptoms persist.

Inhalation: No Symptoms expected (capsules). Remove to fresh air (powder).

PRODUCT RANGE

Riva Silver Capsules
50 x Riva Silver Capsules
Re-order: 8670008

Riva Silver Powder/Liquid Kits
10g (8.9mL) Riva Silver Liquid Bottle
32g Riva Silver Powder jar
accessories
Re-order: 8670508

Riva Silver Powder refills
32g Riva Silver Powder jar
Re-order: 8670108

Riva Silver Liquid refills
10g (8.9mL) Riva Silver Liquid Bottle
Re-order: 8670900

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C. BLANDING OG PLACERING

| | |
|--|------------|
| CAPSULE SYSTEM | |
| Powder/Liquid ratio (g/g) | 0.72/0.13 |
| Mixing time (min. sec.) | 10 seconds |
| Working time (min. sec.) | 1'30" |
| Initial setting time (min. sec) | 5'30" |
| Finishing time @ 37°C from start of mixing (min. sec) | 6'00" |
| Minimum delivery/capsule | 0.15mL |
| Test conditions: temperature (23 ± 1°C), relative humidity (50 ± 10%), ISO 9917-1 (Non translucent radiopaque glass polyalkenoate cement). Minimum net contents per capsule: 0.72g powder and 0.14g (0.12mL) liquid. | |

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