SDI CONSEAL F CONSEAL

E
N
G
L
I
S
H

PIT AND FISSURE SEALANT INSTRUCTIONS FOR USE

Conseal f and Conseal are the lowest viscosity sealants with the greatest flowability. Their optimal flow characteristics achieve superior homogeneous penetration. The lower shrinkage of SDI's unique resin system results in a tighter seal that optimizes the mechanical retention.

Conseal f is fluoride releasing and combines clinically proven SDI sealant technology with an intensive initial fluoride boost in conjunction with long term fluoride treatment. The unique Conseal F filler component, with its blend of particles and high surface-to-area volume ratios maximizes the sealant's fluoride release.

INDICATIONS:

- Sealing of teeth with sound but vulnerable pit and fissure sites.
- Sealing of teeth with minimal pit and fissure caries.

CONTRA-INDICATIONS:

- Do not apply over carious areas.
- Do not apply to unetched areas.
- Product may cause skin irritations to some people. In such cases, discontinue use and seek medical attention.

PRINCIPAL COMPONENT:

Multifunctional methacrylic ester.

APPEARANCE:

A low viscosity resin which changes to a solid mass after light curing.

INSTRUCTIONS FOR USE:

- Clean the occlusal surface of teeth with sound but vulnerable pit and fissure sites by using a non-fluoridated, oil-free prophylaxis paste. Remove minimal pit and fissure caries with small round burs in a slow speed handpiece. Do not extend into the sound portion of the fissure system. Rinse thoroughly with water. Isolate tooth.
- Dry the surface to be etched with clean, dry, oil-free air.
- Etch occlusal surface with 37% phosphoric acid such as Super Etch or Super Etch Low Viscosity (LV), (SDI Limited) for at least 30 seconds but no more than 60 seconds. IMPORTANT: If using an etch syringe, ensure that the dispensing tip hub is **firmly** attached to the syringe by twisting the hub securely onto the syringe. On first usage or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the material viscosity and rate of extrusion.
- Wash thoroughly with water.
- Completely dry tooth with clean, dry oil-free air for 15 seconds. IMPORTANT: A dry etched enamel surface should have a matt white appearance. If this is not so or the area has become contaminated with saliva or blood, repeat steps 2-5.
- Apply sealant. IMPORTANT: If using the syringe system, ensure that the dispensing tip hub is **firmly** attached to the syringe by twisting the hub securely onto the syringe. On first usage or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the material viscosity and rate of extrusion. Place the tip against the surface to be coated and slowly extrude the sealant. **Note: After each use, slightly pull back the syringe plunger to prevent dripping.**
- Light cure sealant for 20 seconds using SDI's Radii Plus LED Curing Light or any other visible light curing device (440 - 480 nm wavelength). Further increments may be added provided proper isolation and dryness are maintained. Wipe off the oxygen inhibited layer and ensure no voids.

STORAGE AND HANDLING:

- Store at temperatures between 10° and 25°C (50°-77°F).
- Use at room temperature between 20° and 25°C (68°- 77°F).
- Exposure to temperatures outside the recommended range in this instruction sheet may shorten shelf life.
- Exposure to ambient and operating lights may cause premature curing of the sealants.

PRECAUTIONS:

- For professional use only.
- Keep out of reach of children.
- Shake sealant bottles before use.
- When not in use replace the cap tightly of syringes and bottles.
- Avoid prolonged contact of etchant and sealant with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization.
- Do not take internally.
- Do not use after expiry date.
- Do not expose material to direct light.
- Caution: Federal Law restricts this device to sale by or on the order of a dentist.
- SDS available at www.sdi.com.au or contact your regional representative.
- This product conforms to ISO 6874 (type 2 light cured pit and fissure sealant).

FIRST AID:

- Eye (contact):** Wash thoroughly with water and seek medical attention.
- Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- Ingestion:** Drink plenty of water/milk. Seek medical attention if symptoms persist.
- Inhalation:** No symptoms expected.

PACKAGING:

- Conseal F (white): syringe kit, bulk syringe kit, bottle kit, bottle refill, complet refill.
- Conseal Clear: syringe kit, bottle refill
- Conseal Light Grey: bottle kit, bottle refill

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