SDI TRAY FABRICATION

- ENG TRAY FABRICATION
- DEU ANFERTIGUNG DER ZAHNSCHIENE
- FRA FABRICATION DES GOUTTIÈRES
- ITA PREPARAZIONE DELLE MASCHERINE
- NLD TRAY FABRICATIE
- POR FABRICAÇÃO DE MOLDEIRAS
- SPA FABRICACIÓN DE LA CUCHARILLA

PAK41275 D

Made in Australia by SDI Limited
3-15 Brunson Street, Bayswater, Victoria 3153, Australia

Australia 1 800 337 003
Austria 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Italy 00800 0225 5734
Spain 00800 0225 5734
UK 00800 0225 5734
USA & Canada 1 800 228 5166

EU REP

SDI Germany GmbH
Hansestrasse 85
51149 Cologne, Germany
T: +49 2203 9255 0

Last Revised: 11-2025




ENG INSTRUCTIONS FOR USE

Material for fabrication of mouth tray. The mouth tray delivers topical dental treatments.

REQUIRED MATERIALS:

- Teeth model – trim base of model to the gingival margin
- SDI Tray material
- SDI Spacer Resin
- Dental Curing Light
- Vacuum former

PROCEDURE

1. Apply approximately 0.5 mm of Spacer Resin to the teeth on the model selected for treatment. Apply to about 2.5 - 3.5 mm from the gingival line. 
2. Light cure (Radii Xpert or Radii Plus (SDI Limited)) each tooth for 20 seconds. 
3. Firmly attach a piece of tray material into the holder. Note: We recommend the tray holder to be fitted with a rubber gasket. This will prevent slipping of the material during heating.
4. Position the tray holder for heating. 
5. Place the model on the centre of the base.
6. Position the pre-heated heating element over the tray. 
7. Heat until the material starts to sag. Turn the heat off. Note: The thickness of the tray will depend on a number of factors, e.g. heating time, position of heating element above tray etc.
8. Lower the tray onto the model. 
9. Apply vacuum for approximately 10 seconds. 
10. Allow the model to cool before removing from vacuum former.
11. Carefully remove the tray away from the impression. 
12. Trim the tray to the gingival margin.

CAUTION: Ensure that the dispensing tip hub of the Spacer Resin is firmly attached to the syringe by twisting the hub securely onto the syringe.

Precautions: Avoid prolonged contact of Spacer Resin with skin when dispensing, as it may cause inflammation or skin sensitization. Any persons having known resin allergies should immediately discontinue the use of Space Resin. SDS for Spacer Resin is available at www.sdi.com.au or contact your regional representative.

NOTE: On first usage, or after prolonged storage, extrude a small amount of Spacer Resin onto a mixing pad for familiarity with the viscosity and rate of extrusion. Use at room temperature. Do not expose material to direct light.

Spacer Resin 2 syringe refill
Reorder 7700012
2 x 1g Spacer Resin syringes accessories

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2. Fotocure (Radii Xpert or Radii Plus (SDI Limited)) cada diente durante 20 segundos.