INSTRUCTIONS FOR USE

ULTRA UNIVERSAL RESTORATIVE MATERIAL
ULTRA UNIVERSAL RESTAURATIONS MATERIAL
MATÉRIEL DE RESTAURATION ULTRA UNIVERSEL
MATERIALE DA RICOSTRUZIONE ULTRA UNIVERSALE
MATERIAL RESTAURADOR UNIVERSAL ULTRA
MATERIAL RESTAURADOR UNIVERSAL
ULTRA UNIWERSALNY MATERIAL ODTWÓRCZY
ULTRA UNIVERSALT FYLLNINGSMATERIAL
ULTRA UNIVERSAL FYLDNINGSMATERIALE
LOISTAVA UNIVERSAALI PAIKKAMATERIAALI
UNIVERSELT RESTAURERINGSMATERIALE
ΕΥΡΕΙΑΣ ΧΡΗΣΗΣ ΥΛΙΚΟ ΑΠΟΚΑΤΑΣΤΑΣΗΣ
GÜÇLENDİRİLMİŞ UNIVERSAL DOLGU MADDESİ
ULTRA UNIVERZÁLNÍ VÝPLŇOVÝ MATERIÁL
COMPOZIT UNIVERSAL PENTRU OBTURAJI ȘI RESTURĂRI
УЛТРА УНІВЕРСАЛІЕН РЕСТОРАТІВЕН МАТЕРІАЛ
ULTRA UNIVERZÁLIS TÖMŐANYAG
ULTRA UNIVERZALNI RESTORATIVNI MATERIAL
超级通用修复材料

Made in Australia by SDI Limited
Bayswater, Victoria 3153
Australia 1 800 337 003
Austria 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Italy 00800 0225 5734
Spain 00800 0225 5734
United Kingdom 00800 0225 5734
USA & Canada 1 800 228 5166
www.sdi.com.au

EC

REP

SDI Germany GmbH
Hansestrasse 85
51149 Cologne
Germany
T+49 2203 9255 0

PAK41333BR B

SDI

AURA

ULTRA UNIVERSAL
RESTORATIVE MATERIAL
INSTRUCTIONS FOR USE

ENGLISH

INTRODUCTION

aura ultra universal restorative system is a the easy solution to creating as close to nature anterior and posterior restorations. Based on the natural make-up and colours of the dentine and enamel, **aura** uses a logical yet innovative approach to shade-matching and layering in one single composite system.

INDICATIONS

	Enamel Shade	Dentine Shade
Anterior restorations	•	•
Posterior restorations	*	•
Indirect inlays, onlays & veneers	*	•
Core build ups		•
Splinting	•	•
Sandwich technique with glass ionomer		•

*For use as thin overlay over Dentine Chroma shade.

CONTRA-INDICATIONS

- Pulp capping.
- Do not use in conjunction with any eugenol containing materials.
- Allergien with allergy to acrylics

SHADE RANGE

Enamel shades:

- **aura** enamel shades are microfilled composites based on different colours of natural enamels.
- E1 has translucency, high opalescence and is whiter. It is designed to emulate young enamel.
- E2 has intermediate translucency, is neutral in colour and less opalescence. It is designed to replace adult enamel.
- E3 has a yellow tint with high translucency and is suitable for replacing aged enamel.
Note: **aura** enamel shades are not radiopaque

Dentine shades:

- **aura** dentine shades are nanohybrid composites based on single opacity, single hue (base colour similar to the “A” series of the Vita Classical Shade Guide), fluorescent with a broad range of chroma. “DC” stands for “dentine chroma” and is equally spaced out in ascending chroma value. A bleached dentine, Db, is also available for bleached teeth.

SHADE SELECTION

Note: Shade selection must be carried out before isolation and ideally should be done under day light using the **aura** shade guide provided.

1. Select the dentine shade (Db / DC1 / DC2 / DC3 / DC4 / DC5 / DC6 / DC7) by matching the shade tab to the cervical area where the enamel has minimal thickness.
2. Select the enamel shade (E1 / E2 / E3) by matching the shade tab to the incisal of the neighbouring teeth or tooth to be restored.

PLACEMENT PROCEDURE:

1. Isolate tooth.
2. Prepare cavity using standard technique for composite resin restorations.
Note: In cases where pulp protection is necessary, use a calcium hydroxide liner.
3. Apply a dentin/enamel bonding agent such as SDI Go! as per manufacturer's instructions.
4. Place **aura** Dentine / Enamel shade in increments of no more than 2mm.
5. Light cure each increment using a high powered LED curing light (460 - 480nm wavelength) for:
 - 10 seconds for E shades (Enamel Shades)
 - 20 seconds for Db, DC (Dentine Shades)
6. Finish restoration using standard technique.

STORAGE AND HANDLING

- Store at temperatures between 10° and 25°C (50°-77°F).
- Use at room temperature between 20° and 25°C (68°- 77°F).
- Syringe system: when not in use replace the cap tightly.

PRECAUTIONS

- Avoid prolonged contact of **aura** with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization.
- For professional use only.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiry date.
- Do not expose material to direct light.
- Caution: Federal Law restricts this device to sale by or on the order of a dentist.
- SDS available at www.sdi.com.au or contact your regional representative.

FIRST AID

- **Eye (contact):** Wash thoroughly with water and seek medical attention.
- **Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- **Ingestion:** Rinse mouth thoroughly. Drink plenty of water/milk. Seek medical attention if symptoms persist.
- **Inhalation:** No symptoms expected.

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