
## INSTRUCTIONS FOR USE

ENGLISH

Thank you for purchasing the Radii-Cal CX LED curing light. The Radii-Cal CX is a powerful, collimated cordless, lightweight curing light using an LED for the polymerization of restorative materials that contain the photoinitiator Camphorquinone. The Radii-Cal CX produces a strong blue light in the wavelength range of 440 – 480 nm, the relevant range for Camphorquinone containing products.

## PRODUCT CONTENTS LIST

- Handpiece
- Charger
- Plug pack
- 100 barrier sleeves
- 3 lens caps
- Light shields (5 Pack)

## UNPACKING AND INITIAL CHARGING

The Radii-Cal CX has been carefully packed to withstand reasonable damage during transit. Ensure all parts contained within the box match the product contents list above.

- 1) Remove all parts from the box.
- 2) Place the plug pack cable into the charger.
- 3) Plug the plug pack into an available power outlet and turn power outlet on. The charger will conduct a brief LED test: the charge light LED should momentarily show red/green and the 5 blue radiometer LEDs should turn on/off.
- 4) Place the curing light handpiece into the charger. The handpiece can be rotated 360° in the charger. Depending on the length of time from manufacture to package opening, the charge remaining in the battery can vary. The charge light will show the below three conditions.

Charge Light Color	Battery condition
Red	discharged
Flashing Green	charging (OK to use)
Green	charged

To ensure the highest level of performance, allow the Radii-Cal CX to fully charge prior to first use. On average, it will take approximately 1-3 hours to initially charge the Radii-Cal CX.

## OPERATION

## • Turning unit on and off

Simply pick up the curing light handpiece, press the on/off button once to activate and de-activate the light. Audible beep facilitate accurate time measurement of curing. Replace the Radii-Cal CX back in the charger when not in use.

## • Audible beeps

An audible beep will be heard after the first 15 seconds (5 second ramp, 10 seconds full cure) one audible beep after a further 10 seconds. The full cycle will continue as follows;

Audible beeps	Time between beeps (seconds)	Total elapsed time (seconds)
1 short beep	15 sec	15 sec
2 short beeps	10 sec	25 sec

At the end of 25 seconds one beep will be heard before the blue LED is automatically switched off.

## • Lens Cap

The Radii-Cal CX has been supplied with a lens cap fitted. The lens cap should be inspected prior to each use for scratches, cracks or foreign matter that may reduce the effectiveness of the LED light source. Generally lens caps should be replaced every 2 weeks. To remove the lens cap carefully unscrew the old lens cap counter-clockwise and fit the new lens cap by tightening clockwise.

## • Built in Radiometer

The Radii-Cal CX's built-in radiometer provides a simple test to check the functionality of the unit. The radiometer is calibrated for the Radii-Cal CX's power and collimated lens. With the Radii-Cal CX on, carefully position the lens cap on the built-in radiometer. 4 or 5 lights should illuminate indicating a fully functional unit. If only 1,2 or 3 lights illuminate, it could be due to the following reasons:

- 1) There could be some material on the lens cap such as composite: - the lens cap should be replaced.
- 2) The lens cap could be damaged or cracked: - it should be replaced.
- 3) The battery could be extremely low - the unit should be recharged.

In a severe case the light could be malfunctioning: - the unit should be serviced.

**Note:** Portable radiometers can give inaccurate and inconsistent results. These results can vary by as much as 300 mW/cm² or more.

## • Disposable barrier sleeves

Transparent, disposable barrier sleeves are supplied to protect the handpiece and prevent cross contamination. Testing showed no significant difference in spectral and energy output when operating with or without a barrier sleeve in place.

Slide the barrier sleeve over the end of the handpiece before placing the orange light shield on the end of the handpiece.

**It is essential to use disposable barrier sleeves to prevent any liquids getting on the handpiece. Evidence of liquids entering handpiece causing damage due to not using barrier sleeves will void warranty.**

## • Automatic shut off function

Should the Radii-Cal CX LED become too hot due to excessive use (after about 4 minutes of continuous irradiation), the unit will automatically shut off and be un-useable for about 1 minute until it cools down.

## • Changing the battery pack

Carefully unscrew the battery pack counter-clockwise to remove and replace. A single beep will be heard when the battery pack is successfully connected.

## • Cleaning the Radii-Cal CX:

The ideal way to protect patients from cross-infection is to utilise Radii-Cal CX barrier sleeves. Barrier sleeves are also essential to keep the Radii-Cal CX clean. When cleaning the Radii-Cal CX handpiece, ensure control section, lens cap and battery pack are connected as one unit.

- LIQUID CONTACT WITH TERMINALS OF HANDPIECE AND CHARGER WILL CAUSE DAMAGE VOIDING WARRANTY.
- LIQUID CONTACT WITH LED WILL CAUSE DAMAGE VOIDING WARRANTY.
- LIQUIDS MUST NOT ENTER ANY PART OF THE HANDPIECE OR CHARGER.
- DO NOT SPRAY LIQUIDS DIRECTLY ONTO, OR APPLY COPIOUS AMOUNTS OF LIQUIDS TO THE RADII-CAL CX.
- DO NOT AUTOCLAVE THE RADII-CAL CX.

## General Surface Cleaning

General surface cleaning can be conducted with anti-microbial surface disinfectant such as Glutaraldehyde, Chlorhexidine gluconate, and 70% Isopropyl alcohol. When applying the disinfectant to any component of the device, spray the disinfectant onto a piece of cloth and wipe over the **external** areas to be cleaned.

## • Curing Times

Curing times will differ for different formulations of restorative materials. Follow the manufacturers instructions for recommended curing times. Always bench test new materials before use in-vivo.

## TROUBLE SHOOTING

## Problem

	Causes	Solution
• Charge light flashes red	Possible battery fault	Replace battery
• Handpiece blue LED flashes then goes off	Battery charged too low	Re-charge battery
• No lights appear on the charger when the hand piece is placed on the charger	1. Check terminals on the battery are clean and free of debris 2. Check power is getting to charger, turn power outlet on. The charge light LED should momentarily show red/green and the 5 blue radiometer LEDs should turn on/off. 3. Check the power outlet is turned on, if so use a different power outlet. 4. Check the plug pack is correctly connected to the charger. 5. Check that the hand piece is placed into the charger correctly. 6. If the problem still exists after checking 1, 2, 3 and 4 above the charger may be faulty. Return the unit to the supplier for service.	

When returning faulty units please make sure that all parts are returned together in the original packaging.

## • Disposal of battery pack

As a means of protecting the environment, the Radii-Cal CX battery pack is fitted with a Lithium Ion battery pack. Dispose of battery pack in accordance with local legal regulations.

## SAFETY AND PRECAUTIONS

Use of the Radii-Cal CX is restricted to qualified medical professionals in accordance with the safety instructions below:

Radii-Cal CX LED curing light is classified as a Risk Group 2 device according to IEC 62471-1.

The safe operation of this unit depends on strictly following the operating instructions in this manual. The Manufacturer accepts no liability for any damage resulting from the use of this unit for any other purpose than the polymerization of dental materials.

- Do not expose the Radii-Cal CX battery pack to prolonged temperature above 40 degrees Celsius as this may affect its overall estimated working life. The equipment can not be exposed or immersed in water or wet locations. The equipment does not have protection against liquid penetration.
- Equipment not suitable for use in the presence of flammable anesthetic mixtures with air oxygen or nitrous oxide.

## CAUTION!

Radii-Cal CX LED curing light produces high-output curing energy. Prolonged exposure to the eyes, exposed skin, oral mucosa, or tissues near the pulp can result in irritation or injury. It is important to observe the following precautions and procedures:

**Positioning:** When curing, position the LED Light attachment tip directly over the cure material. Do not place light directly on or toward unprotected gingiva or skin.

**Eye exposure:** Do not look directly at the light emitted from this curing unit. Do not operate this device without suitable eye protection for the operator, assistant, and patient. Select eye protection that filters blue-violet and ultraviolet light. Securely attach Light Shield before each use. Avoid prolonged exposure of unprotected eyes to light scatter or reflected light. If light exposure causes blurring/staring of vision or visual after-image that persists, consult a physician.

**Skin exposure:** Avoid prolonged exposure to skin. If skin rash, sensitization or other reaction occurs, discontinue use and seek medical attention.

**Oral mucosa exposure:** Avoid prolonged exposure to oral soft tissues. Prolonged exposure may cause soft tissue irritation or burns. Most minor reactions require only thorough cleansing and palliative treatment. If irritation persists, seek medical attention.

## Do not use the Radii-Cal CX in patients or by users with:

- heart pacemaker implants who have been advised to be cautious in regards to their exposure to small electrical devices.
- a history of photo biological reactions (including individuals with urticaria solaris or erythropoietic protoporphyria) or who are currently on photo-sensitizing medication (including 8-methoxysoralen or dimethylchlorotetraacycline).
- a history of cataract surgery. These people may be particularly sensitive to the exposure to light and should be discouraged from Radii-Cal CX treatment unless adequate safety measures, such as the use of protective goggles to remove blue light, are undertaken.
- a history of retinal disease. These people should seek advice from their ophthalmologist prior to operating the unit. In operating the Radii-Cal CX unit, this group of individuals must take extreme care and comply with all safety precautions (including the use of suitable light-filtering safety goggles).

## WARRANTY PERIOD

SDI Ltd., the Manufacturer, extends a two year warranty against defects in materials or workmanship to the original purchaser of this product.

The Manufacturer agrees to correct any defects which develop within the warranty period, either by repair or replacement, at its option. This warranty is valid providing factory inspection indicates that any such defect developed during normal and proper use subject to the conditions below.

## WARRANTY CONDITIONS - PLEASE READ CAREFULLY

- Please fill out the warranty card included in kit and send promptly back to SDI with proof of purchase. Failure to do so may void your warranty.
  - Alternatively goto <http://www.sdi.com.au/warranty> to register your warranty.
- a) Claims for damage in shipment should be filed promptly with the transportation company.  
b) All shipments claimed defective can only be returned to the Manufacturer with the written consent of the Manufacturer. All returned products must be accompanied by a full description of discrepancy or malfunction.  
c) Only SDI authorized service personnel are allowed to carry out repairs to the Radii-Cal CX. Manufacturer shall be released from all obligations under this warranty if repairs or modifications are made by persons other than Manufacturer's own authorized service personnel.  
d) Use only SDI parts exclusively to replace defective components. The product warranty does not cover any damage resulting from the use of third-party replacements parts.

e) **Manufacturer shall be released from all obligations under this warranty in the event of improper installation; damage due to casualty; or obvious misuse including but not limited to damage from liquids and autoclaving.**

f) Manufacturer shall in no event be liable for any consequential damages arising from the misuse of the product.

g) Only the above stated warranty shall be applicable, and all other warranties, expressed or implied, are hereby excluded.

All correspondence concerning the product should specify the model and serial number.

## TECHNICAL DATA

Handpiece (battery, control section and standard LED attachment)

wavelength range:	440 nm – 480 nm
peak wavelength range:	460 nm
RAMP to full intensity:	5 seconds
light intensity:	1200 mW/cm²(peak)
duration of continuous use with fully charged battery:	20 seconds
total continuous run time:	3 hours, 20 minutes
dimensions:	22mm diameter x 243mm long
weight:	144g / 5.1 ounces
battery:	Li-ion 7.4V/14wh
operating temperature:	10°C - 40°C (50°F - 104°F)
operating Relative humidity:	10% - 85%

## Plug Pack

AC input:	100-240V 0.8A, 50/60 Hz
DC output:	12V 1.66A, 20.0W MAX
dimensions:	34mm high x 41mm wide x 62.5mm
long:	
weight:	104g / 3.65 ounces
operating temperature:	10°C - 40°C (50°F - 104°F)
operating Relative humidity:	10% - 85%

## Charger

time to charge battery:	60 – 180 minutes
operating temperature:	10°C - 40°C (50°F - 104°F)
operating Relative humidity:	10% - 85%
inserted:	
dimensions:	255mm
long:	70mm wide x 40mm high x 130mm
weight:	84g / 3 ounces

According to classification IEC 60601-1

## REPLACEMENT PARTS AND ACCESSORIES

Replacement parts and accessories which can be installed without special training or equipment, can also be obtained from your authorized SDI dealer.

Replacement lens caps (5 Pack)  
Barrier sleeves (1000 pack)  
Light shield (5 Pack)

Re-order Code  
5600310  
5600055  
5600056

## Transport and storage environment:

Ambient temperature: 0°C - 40°C (32°F - 104°F)

Relative humidity: 10% - 85%

Atmospheric pressure: 0.5 atm - 1.0 atm (500 hPa to 1060 hPa)

CAUTION  
(refer to accompanying instructions)

Type B Applied Part  
(Handpiece)

HIGH INTENSITY LIGHT  
Do not look directly into the light source.  
The tip of the LED attachment is the only emission aperture that relates to the above hazard warnings.

Made in Australia by SDI Limited  
3-15 Brundson Street  
Bayswater, Victoria 3153, Australia  
Australia 1 800 337 003  
Austria 00800 0225 5734  
Brazil 0800 770 1735  
France 08000 0225 5734  
Germany 0800 1005759  
Italy 0800 0225 5734  
Spain 08000 0225 5734  
United Kingdom 00800 0225 5734  
USA & Canada 1 800 228 5166  
[www.sdi.com.au](http://www.sdi.com.au)

EC | REP

SDI Germany GmbH  
Hansestrasse 85  
51149 Cologne  
Germany  
+49 2203 9255 00

Last Revised: 05-2024

Table 201 - Guidance and MANUFACTURER'S declaration - ELECTROMAGNETIC EMISSIONS - for all equipment and systems

Guidance and manufacturer's declaration - electromagnetic emissions		
The Radii Cal CX is intended for use in the electromagnetic environment specified below. The customer or the user of the Radii Cal CX should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Radii Cal CX uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Radii Cal CX is suitable for use in all establishments, including domestic establishments. The Radii Cal CX is powered by 2 x 4.2V DC Lithium ion batteries thus Harmonic emission and Voltage fluctuation/flicker emissions are not applicable.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Table 202 - Guidance and MANUFACTURER'S declaration - Electromagnetic IMMUNITY - for all equipment and systems

Guidance and manufacturer's declaration - electromagnetic immunity			
The Radii Cal CX is intended for use in the electromagnetic environment specified below. The customer or the user of the Radii Cal CX should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 204 - Guidance and MANUFACTURER'S declaration - Electromagnetic IMMUNITY - for all equipment and systems that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration - electromagnetic immunity			
The Radii Cal CX is intended for use in the electromagnetic environment specified below. The customer or the user of the Radii Cal CX should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the TENS PRO 900, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range .b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Radii Cal CX is used exceeds the applicable RF compliance level above, the Radii Cal CX should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or re-locating the Radii Cal CX.			

Table 206 - Recommended separation distances between portable and mobile RF communications equipment and the equipment and systems - for equipment and systems that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the Radii Cal CX		
The Radii Cal CX is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Radii Cal CX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Radii Cal as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0,01	0.12	0.23
0,1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

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