

wave

FLUORIDE RELEASING FLOWABLE COMPOSITE

FLUORIDE FREISETZENDES FLIESSFÄHIGES KOMPOSIT

COMPÓSITO DE ALTA FLUIDEZ COM FLÚOR

COMPOSITE FLUIDO QUE LIBERA FLUOR

COMPOSITE A FAIBLE VISCOSITÉ DÉGAGEANT DU FLUOR

COMPOSITO FLUIDO CON RILASCIO DI FLUORO

FLUORIDE RELEASING FLOWABLE COMPOSIT

FLUORAFGIVENDE FLYDENDE KOMPOSIT

FLYTENDE FLUORAVGIVENDE KOMPOSITT

FRUORIA VARAUTTAVA FLOW -YHDISTELMÄMUOVI

ΡΕΥΣΤΗ ΣΥΝΘΕΤΗ ΡΗΤΙΝΗ Η ΟΠΟΙΑ ΑΠΕΛΕΥΘΕΡΩΝΕΙ ΦΘΟΡΙΟ

MATERIAŁ KOMPOZYTOWY TYPU FLOW UWALNIAJĄCY FLUOR

FLUORID LEADÓ FOLYÉKONY KOMPOZIT

FLUORIIDI ERALDAV VOOLAV KOMPOSIIT

FLUORID-SPROŠČUJOČI TEKOČI KOMPOZIT

FLUORU IZDALOŠS PLŪSTOŠAIS KOMPOZĪTS

SKYSTAS TAKUMO FLUORĄ ATPALAI DUOJANTIS KOMPOZITAS

TEKUTÝ KOMPOZITNÍ MATERIÁL UVOLŇUJÍCÍ FLUORIDY

TEKUTÝ KOMPOZIT UVOĽŇUJÚCI FLUORID

ウエーブ フッ素徐放型フロアブルコンポジット

可释氟流动性树脂



SDI

a flowable for all your clinical situations

wave




**THE DENTAL
ADVISOR**
++++^{1/2}

The Dental Advisor.
October 1998,
Vol. 15, No. 8.

For use in just about every restoration

- | | |
|---|---|
| • Class V restorations | • Temporary crowns |
| • Tunnel / air abrasion preparations | • Cement porcelain/ceramic/ composite veneers |
| • Pit and fissure sealants | • Splint mobile teeth |
| • Minimal Class I, II, III, IV restorations | • Attach fiber bridges |
| • Radiopaque cavity liners | • Restorations in deciduous teeth |
| • Cervical abrasion/erosion lesions | • Repair porcelain restorations |
| • Abfraction lesions | • Block out undercuts in inlay/onlay/crown preparations |
| • Implant seals | • Cover stains |
| • Incisal abrasions | • Adhering tooth jewelry |
| • Repair enamel defects | • Ceramic inserts |
| • Minor core build-ups | |

A choice of 3 viscosities



Flow at 90° angle for 2 minutes

The wide number of viscosities available in the Wave range ensures the clinician has a choice for all different clinical conditions.

SDI Flowable: Compare to:

Wave	UltraSeal XT Plus* Premise Flowable*
Wave MV	Filtek Supreme Plus Flowable* Heliomolar Flow* Tetric Flow*
Wave HV	Revolution Formula 2* Flow-It! ALC* Point 4 Flowable* Filtek Flow* AeliteFlo LV* Esthet-X Flow* AeliteFlo*


Nanofilled

Wave, Wave MV and Wave HV contain specially treated nano-fillers to maximize polishability, wear resistance and strength. Ease of extrusion is improved with the addition of these specially treated nano-fillers.

Filler loading


Wave, Wave MV and Wave HV all have excellent wear profiles due to their high filler loading level:

% by weight



Fluoride release

Fluoride's cariostatic effect enhances remineralization and inhibits enamel demineralization. In an aqueous environment, the fluoride ions in Wave, Wave MV and Wave HV diffuse from the resin into the surrounding tooth. This ionic movement is caused by oral fluid passing in and out of the resin and tooth, acting as a carrier for the fluoride ions.



Natural fluorescence

Wave, Wave MV and Wave HV contain a fluorescent agent that ensures perfectly natural looking teeth under ultraviolet lights.

No one will know!



Wave on natural teeth under UV light


Radiopacity

Wave, Wave MV and Wave HV are all visible on x-rays. This allows easy identification of secondary caries at the interface.

Compressive strength

Wave's, Wave MV's and Wave HV's high compressive strength values expand the suitability for multiple applications.

MPa @ 1 hour**



UDMA


Wave's UDMA resin system minimizes shrinkage and reduces post operative sensitivity.

Specialty shades

Wave's 14 shades are ideal for all anterior and posterior restorations: A1, A2, A3, A3.5, B1, B2, B3, C2, C3, D3, OA2, incisal, gingival and bleach. Gingival shade is ideal for restoring cervical lesions in elderly patients. For teeth lighter than B1, Bleach shade is available. A semi-translucent shade, Incisal, is perfect as a final layer for composite restorations.

Tip length comparison

Wave's complet tip length is 12.2mm, 40% longer than the alternative of 8.7mm. Wave's longer tip length allows more direct placement of the material into difficult to reach areas.



instructions:


- 1 Clean and isolate tooth
- 2 Remove excess water. Preferably keep moist, but a dry surface is ok
- 3 Apply Go! to saturate all internal surfaces and leave in place for 20 seconds. Or, apply other bonding agent according to manufacturer's instructions



- 4 Blow firmly with dry, oil free air for 5 seconds. Leave surface glossy




- 5 Light cure for 10 seconds



- 6 Directly inject Wave in increments of 2mm or less in:


- 6.1 Class V restorations **or**

- 6.2 Tunnel preparations **or**



- 6.3 Pit and fissure sealants **or**

- 6.4 Minimal Class I, Class II, Class III and Class IV restorations, or other indications as required



- 7 Cure Wave for 20 seconds in increments of 2mm

- 8 Polish and finish

wave



Syringe refills

1 x 1g syringe, 5 single use disposable tips

Shade	Wave	Wave MV	Wave HV
A1	7510102	8310102	8210102
A2	7510203	8310203	8210203
A3	7510304	8310304	8210304
A3.5	7510405	8310405	8210405
B1	7511101	8311101	8211101
B2	7511211	8311211	8211211
B3	7511314		
C2	7512223	8312223	8212223
C3	7512324		
D3	7513334		
OA2	7510204	8310204	8210204
Incisal	7514142		
Gingival	7514143		
Bleach	7514147	8314147	8214147

Disposable tips for Wave syringes

20 x single use disposable tips (20 gauge)

Reorder 7500001

200 x single use disposable tips (20 gauge)

Reorder 8100062

Disposable tips for Wave MV & HV syringes

20 x single use disposable tips (18 gauge)

Reorder 8100065

Wave syringe introductory kit

4 x 1g syringes A2, B3, C3, OA2

20 single use disposable tips

Reorder 7500000

Wave MV syringe introductory kit

4 x 1g syringes A2, B2, C2, OA2

20 single use disposable tips

Reorder 8300000

Wave HV syringe introductory kit

4 x 1g syringes A2, B2, C2, OA2

20 single use disposable tips

Reorder 8200000

Wave syringe bulk kits

10 x 1g syringes

40 single use disposable tips

A1 7500011

A2 7500012

A3 7500013

Complet refills

20 x 0.25g complets

Shade	Wave	Wave MV
A1	7520001	8320001
A2	7520002	8320002
A3	7520003	8320003
A3.5	7520004	

* Not the registered trademarks of SDI Limited.

** Source-Published and SDI test data.



Made in Australia by SDI Limited

Bayswater, Victoria 3153

Australia 1 800 337 003

Austria 00800 2255 734

Brazil 0800 770 1735

France 00800 022 55 734

Germany 0800 100 5759

Ireland 01 886 9577

Italy 800 780625

Spain 00800 022 55 734

United Kingdom 00800 022 55 734

USA & Canada 1 800 228 5166

www.sdi.com.au

