DOSAGE:

- Each Step 1 & Step 2 capsules can treat up to 5 teeth per patient.
- Limit of one (1) Step 1 & Step 2 capsule per visit. If needed a second application of Riva Star may be administered after 1 week.

PRECAUTIONS:

- For professional use only.
- Do not repeat treatment on pregnant or lactating women.
- Solution from **silver** capsule will stain any surface if solution from **green** capsule is not applied effectively. (See HANDLING STAINS.)
- Apply Gingival Barrier before applying solution from **silver** capsule to avoid soft tissue burn.
- Avoid contact of solution from **silver** capsule with oral tissue, eyes and skin as it is corrosive.
- Avoid contact of solution from **green** capsule with eyes.
- Keep out of reach of children.
- Do not ingest.
- Wear gloves and eye protection when handling.
- Single use only.
- Do not use after expiry date.
- Solution from **silver** capsule contains 0.06g / mL fluoride ions.
- SDS is available at www.sdi.com.au or contact your SDI representative.

STORAGE AND HANDLING:

- Do not pierce foil of **silver** capsule or **green** capsule until ready to use.
- Store in a cool and dark place at temperatures between 2° and 8°C (35° - 45°F).
- Use at room temperature between 20° and 25°C (68° - 77°F).

FIRST AID:

- Eyes (contact):** Hold eyelids apart and flush the eye continuously with running water. Continue flushing until advised to stop by a Poisons Information Centre or a medical doctor, or for at least 15 minutes. Seek medical attention if symptom persists.
- Skin/oral soft tissue (contact with solution from silver capsule):** Immediately apply a generous amount of solution from **green** capsule until creamy precipitate turns clear. Flush with running water. Seek medical attention if symptom persists.
- Ingestion:** Not expected due to small dose/brush. If affected, do NOT induce vomiting. Give a glass of milk or water. Seek medical attention if symptom persists.
- Inhalation:** Not expected due to small dose/brush. If affected, remove to fresh air.

HANDLING STAINS:

- Skin:** Apply solution from **green** capsule with gentle rubbing until stain disappears. Rinse with water.(If nothing is done, stain will gradually disappear after a few days).
- Clothing:** Soak in Napisan* (or other fabric stain remover) overnight and wash per normal.
- Hard surfaces:** Immediately apply household bleach (sodium hypochlorite) then wipe off with water. For stubborn stains, saturate a paper napkin with household bleach and place over affected area. Leave overnight and mop up with water.


Napisan* is not the registered trademark of SDI Limited.

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