SDI RIVA STAR AQUA

ENG	TOOTH DESENSITISING AGENT
BUL	ДЕСЕНЗИТИРАЩ АГЕНТ ЗА ЗЪБИ
CES	PROSTŘEDEK K DESENBILIZACI ZUBU
DAN	MIDDEL MOD TANDFØLSOMHED
DEU	DESENSIBILISIERUNGSMITTEL FÜR ZÄHNE
ELL	ΠΑΡΑΓΟΝΤΑΣ ΑΠΕΥΑΙΣΘΗΤΟΠΟΙΗΣΗΣ ΔΟΝΤΙΩΝ
FIN	AINE HAMPAAN DESENSITISOINTIIN
FRA	AGENT DÉSENSIBILISANT DENTAIRE
HUN	FOG ÉRZÉKENYSÉG MEGSZÜNTETŐ
ITA	AGENTE DESENSIBILIZZANTE DENTALE
KOR	치아 탈감작제
LIT	DANTŲ JAUTRUMĄ MAŽINANTI MEDŽIAGA
NOR	DESENSIBILISERINGSMIDDEL FOR TENNER
POL	PREPARAT ZMNIEJSZAJĄCY NADWRAŹLIWOŚĆ ZĘBÓW
POR	AGENTE DESSENSIBILIZANTE DENTAL
RON	AGENT DESENSIBILIZANT DENTARE
SLK	ZUBNÝ DEAKTIVÁTOR CITLIVOSTI
SPA	AGENTE DESENSIBILIZANTE DEL DIENTE
SWE	DESENSIBILISERANDE MEDEL FÖR TÄNDER

Made in Australia by SDI Limited
3-15 Brunsdon Street, Bayswater, Victoria 3153, Australia

Australia 1 800 337 003
Austria 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Italy 00800 0225 5734
Spain 00800 0225 5734
UK 00800 0225 5734
USA & Canada 1 800 228 5166

EU REP
SDI Germany GmbH
Hansestrasse 85
51149 Cologne, Germany
T: +49 2203 9255 0

P/AK41451 A

Last Revised: 08-2026

2°C 35°F 77°F Allergies & Precautions 0297 MD Rx only Paper

ENG INSTRUCTIONS FOR USE

DESCRIPTION
Riva Star Aqua is an ammonia-free two-step kit containing Silver Fluoride (AgF) and potassium iodide (KI) solutions to desensitise teeth.

Step 1 is 38% (w/v) silver fluoride (0.06g / mL fluoride ions) in an aqueous solution.

Step 2 is an aqueous potassium iodide solution.

INTENDED USE
Riva Star Aqua is used by dental professionals to desensitise cervical tooth hypersensitivity.

CONTRAINDICATIONS
- Do not use on patients with known thyroid conditions.
- Do not use on patients who have a history of allergies to silver, fluoride, potassium and/or iodine.
- Do not use on pregnant or lactating women.

Note: If in doubt, consult a doctor or an allergy specialist before use.

SIDE EFFECTS
1. Contact with soft tissues may cause inflamed tissues / aphthous ulcers in the mouth and temporary staining.
2. Permanent staining of unwound tooth structure.

INTERACTIONS
Residual precipitate or moisture may reduce the bond strength of resin-based adhesives.

A. INSTRUCTIONS FOR USE

TOOTH PREPARATION
1. Using a non-fluoride prophylaxis paste in a rubber cup, clean the tooth/teeth to be treated and abutting teeth.
2. Assess condition of teeth. Only teeth with sound dentine and enamel should be treated.
3. Isolate area to be treated with cotton rolls.
4. Dry teeth to be treated and adjoining gingival tissues.
5. Protect the soft tissue by applying petroleum jelly, cocoa butter and/or rubber dam.

B. APPLICATION DOSAGE
Capsule : Each Step 1 and Step 2 capsule can treat up to five (5) teeth per patient.

Bottle : One (1) drop of Riva Star Aqua Step 1 and two (2) drops of Riva Star Step 2 can treat up to five (5) teeth per patient.

Application is via bottle (See B1) or capsule (See B2)

1. Application using Riva Star Aqua bottles:
I. Riva Star Aqua Step 1 bottle: Dispense one (1) drop of the solution (blue label) on a non-absorbent mixing pad or dappen dish.
II. Carefully apply solution to treatment site only; a medium sized micro brush, such as Points (SDI Limited), is recommended.
III. Riva Star Step 2 bottle: Immediately after, dispense two (2) drops of solution (green label) onto fresh non-absorbent mixing pad.
IV. Apply the solution to treatment site, repeating application until the creamy white precipitate turns clear (if this precipitate does not turn clear, staining may result). A medium sized micro brush, such as Points (SDI Limited) is recommended.
Note : Also apply Step 2 solution to any other tooth surfaces that may have come into contact with Step 1 solution.
V. Blot dry.
VI. Remove all protective/isolation materials used in the mouth.
VII. Clean dappen dishes and discard used mixing pads and micro-brushes in accordance with local regulations.

NOTE: Capsules and micro brushes are for single use only. Reuse will result in cross contamination.

STORAGE AND HANDLING

BOTTLE
- Do not dispense either solution from bottle until ready to use. Store bottles upright in the packaging provided.
- Store at temperatures between 2 – 25 °C (35 – 77 °F).
- Use at room temperature between 20 - 25 °C (68 - 77 °F). If stored at cooler temperature allow product to reach room temperature before use.
- Replace cap on bottle and seal tightly immediately after use.

CAPSULE
- Do not pierce foils of capsules until ready to use.
- Store at temperatures between 2 – 25 °C (35 – 77 °F).
- Use at room temperature between 20 - 25 °C (68 - 77 °F). Allow product to reach room temperature before use.

PRECAUTIONS
- Patient must be informed of potential tooth staining prior to treatment.
- Due to potential staining, it is recommended that anterior teeth not be treated.
- Avoid contact of Riva Star Aqua solution with oral soft tissue, eyes and skin.
- Use on soft dentine may cause staining.
- Solution from Riva Star Aqua Step 1 will stain surfaces. Wear gloves and eye protection when handling (see "HANDLING STAINS").
- Possible patient discomfort due to metallic taste.
- Capsules are single use only.
- For professional use only.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiry date.
- SDS available at sdi.com.au or contact your regional representative.

If a serious incident has occurred due to this product, please report this incident to the manufacturer of this product and to the competent authority in your state or country.

FIRST AID
Eye (contact): Open the eye wide and thoroughly wash for at least 15 minutes with running water. Seek medical attention.
Skin (contact): Immediately and thoroughly wash with water. Seek medical attention if symptoms persist.
Ingestion: Do NOT induce vomiting. Rinse mouth with water. Seek medical attention if symptoms persist
Inhalation: Not expected under normal conditions of use. If inhalation occurs, seek medical attention.

HANDLING STAINS
- **Skin:** The stain should disappear gradually after a few days.
- **Clothing:** Soak contaminated clothing in a diluted fabric stain remover and wash per normal.
- **Hard surfaces:** Immediately apply household bleach (sodium hypochlorite) then wipe off with water. For stubborn stains, saturate a paper napkin with household bleach and place over affected area. Leave overnight and mop up with water.

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