SDI POLA OFFICE+

ENG IN OFFICE TOOTH COLOUR RESTORATION SYSTEM
BUL СИСТЕМА ЗА ВЪЗСТАНОВЯВАНЕ НА ЦВЕТА НА ЗЪБИТЕ В ОФИСА
CES SYSTÉM PRO OBNOVU BARVY ZUBNÍ SKLOVINY V ORDINACI
DAN SYSTEM TIL GENOPRETNING AF TANDFARVE PÅ KONTORET
DEU PRAXISSYSTEM FÜR DIE RESTAURATION DER ZAHNFARBE
ELL ΠΡΟΪΟΝ ΛΕΥΚΑΝΣΗΣ ΔΟΝΤΙΩΝ ΓΙΑ ΤΟ ΟΔΟΝΤΙΑΤΡΕΙΟ
FIN VASTAANOTOLLA KÄYTETTÄVÄ HAMPAAN VÄRIN PALAUTUSJÄRJESTELMÄ
FRA SYSTÈME EN CLINIQUE DE RESTAURATION DE LA COULEUR DES DENTS
HUN RENDELŐI FOGSZÍN-HELYREÁLLÍTÓ RENDSZER
KOR 치과에서 시술하는 치아 미백 장비
ITA SISTEMA DI RIPRISTINO DEL COLORE DENTALE IN STUDIO
LIT DANTŲ SPALVOS ATSTATYMO SISTEMA KLINIKOJE
NOR SYSTEM FOR RESTAURERING AV TANNFARGE PÅ KONTORET
POL SYSTEM ODTWARZANIA KOLORU ZĘBÓW W GABINECIE
RON SISTEMA DE RESTAURAŢIA CULORII DINŢILOR ÎN CABINET
SLK SYSTÉM NA AMBULANTNÉ OBNOVENIE FARBY ZUBOV
SPA SISTEMA DE RESTAURACIÓN DEL COLOR DE LOS DIENTES EN LA CONSULTA
SWE SYSTEM FÖR ÅTERSTÄLLNING IAV TANDFÄRG I PRAKTIKEN

SDI Limited, 3-15 Brunsdon Street, Bayswater Victoria 3153, Australia

Austria 00800 0225 5734
France 00800 0225 5734
Germany 0800 1005759
Italy 0800 0225 5734
Spain 00800 0225 5734

EU REP
SDI Germany GmbH
Hansestrasse 85
51149 Cologne, Germany
T: +49 2203 9225 0

Last Revised: 05-2026

2°C 35°F | 8°C 46°F | Allergies & Precautions | CE 0297 | MD | PAK41446 B

ENG INSTRUCTIONS FOR USE

DESCRIPTION
Pola Office+ is a hydrogen peroxide based dental device for professional in-office use, supplied as a dual-component gel in a double-barrel syringe, mixed during application. The active ingredient of the mixed gel, hydrogen peroxide, penetrates the tooth structure and oxidises chromogenic compounds responsible for intrinsic discolouration, supporting restoration of natural tooth colour in clinically indicated cases.

INTENDED USE
Pola Office+ is used to restore the natural tooth colour in eternally discoloured vital and non-vital teeth, which can only be done safely in the controlled environment of a dental practice in the case of difficult access and limited patient compliance.
Pola Office+ is a hydrogen peroxide-based system for restoring tooth colour in the dental office with minimal time.
The gel is pH neutral and contains potassium nitrate for maximum patient comfort. It is mixed in a dual-chamber syringe and can be easily and precisely applied directly to the tooth surface to be treated.

COMPOSITION
The product consists of a dual barrel syringe. One barrel contains a clear gel formulation with 37.5% hydrogen peroxide. The second barrel contains a blue gel formulation without hydrogen peroxide.

When mixed, the combined gel contains: AQUA, HYDROGEN PEROXIDE (30.4 %w/w), PVP, SODIUM HYDROXIDE, ETIDRONIC ACID, POTASSIUM NITRATE, SODIUM STANNATE, ULTRAMARINES (CI77007).

INDICATIONS FOR USE
- Pola Office+ 37.5% is used for the controlled, medically indicated restoration of natural tooth colour in discoloured, vital and devitalised teeth, whose natural or normal colour cannot be restored by professional tooth cleaning.
- Medical restoration of natural tooth colour is indicated when pathological tooth discolouration affects the patient's health-related quality of life.
- Medical restoration of natural tooth colour is also an alternative to much more invasive conservative or prosthetic treatments of the affected teeth (fillings, crowns, veneers) and should therefore always be the first choice.
- The causes of such internal discoloration are, for example: discolouration caused by tetracycline/minocycline, mild to severe dental fluorosis (TF index > 2), devitalized teeth (pulp necrosis), e.g. due to trauma, genetic or age-related degenerative discolouration, porphyria or jaundice, strong deposits of food dyes with penetration into the tooth substance (internal discoloration due to external dyes).
- The product supports the restoration of natural tooth colour in single and multiple teeth, including endodontically treated de-vital teeth, including intracoronal treatment.

CONTRADICATIONS
- Do not use on pregnant or lactating women.
- Do not use on children under 18 years of age.
- Do not use on persons having known sulphur allergies
- Do not use Pola Office+ on any persons having known peroxide allergies.
- Do not use on patients with extremely sensitive teeth.

SIDE EFFECTS
- Possible transient soft tissue irritation if the Gingival Barrier is not applied correctly to prevent product from contacting the soft tissue.
- Some patients may experience soft tissue swelling, redness, tenderness. These side effects are generally temporary. Clinicians should monitor for these potential reactions and manage appropriately.
- Ensure appropriate gingival protection prior to application.

INTERACTIONS
- Pola Office+ does not alter the colour of existing restorative materials (e.g. composites, glass ionomer cements, ceramics); colour mismatch may occur following treatment.
- Residual oxygen after treatment may temporarily inhibit polymerisation of resin-based materials (e.g. adhesives, composites). Delay restoration placement for 2 to 4 weeks after treatment.
- Concomitant use with other peroxide-containing products may increase the risk of tooth sensitivity and soft tissue irritation.
- Post-treatment exposure to chromogenic substances (e.g. coffee, tea, red wine, tobacco) may affect colour stability.

INSTRUCTIONS FOR USE
VITAL TEETH
1. Seal all exposed roots and restorations.
2. Determine and record pre-operative shade (recommended).
3. Clean teeth with ONLY a flour-based pumice according to manufacturer's instructions.
4. Place cheek retractors and then cover exposed lip surface with petroleum jelly.
5. Dry teeth to be treated and adjoining gingival tissues.
6. Protect the soft tissue by applying appropriate gingival protection, such as Gingival Barrier (SDI Limited), according to manufacturer's instruction.
7. Firmly attach a mixing nozzle to the Pola Office+ syringe away from patient. Dispense a small amount of gel on to a mixing pad until a uniform gel is extruded.
8. Using the nozzle as a guide, directly apply a thin layer of gel to all teeth undergoing treatment. A thin layer will help prevent the gel from running. **NOTE:** It is best not to lift the nozzle when moving from tooth to tooth. Use a continuous motion.
9. Leave gel on teeth for 8 minutes (Optional: A curing light may be used at this stage, following manufacturer's instructions).
10. Suction off using a surgical aspirator tip according to the manufacturer's instructions, taking care to retain gingival protection. Note: It is normal for a residual layer of gel to remain on the surface of the teeth.
11. Repeat steps 8-10 up to 4 times.
12. After the last application, suction off remaining gel, then wash with water.
13. Remove gingival protection.
14. Record the final shade (recommended).

INSTRUCTIONS FOR USE NON-VITAL TEETH
1. Radiographically verify adequate endodontic obturation.
2. Determine and record pre-operative shade (recommended).
3. Isolate tooth/teeth undergoing treatment using a rubber dam according to manufacturer's instructions.
4. Open lingual access to pulp chamber and remove all restorative material 2-3 mm apical to the cementoenamel junction. Place barrier material (e.g., glass ionomer cement), into the prepared canal at least 1 mm thick to prevent product from penetrating into the root canal.
5. Place the Pola Office+ gel into the opened pulp chamber and onto the labial surface.
6. Leave gel in place for 8 minutes.
7. Remove the gel with suction, employing standard techniques and instruments.
8. Repeat steps 4-6 until adequate outcome is achieved (up to 4 applications).
9. After the last application is removed, wash the treatment area(s) employing standard techniques and instruments to ensure all residual gel is removed.
10. Seal chamber access cavity with temporary restorative material.
11. Remove rubber dam.
12. Recall patient after 14 days, remove temporary restoration and replace with permanent restoration.
13. Record the final shade (recommended).

Additional notes
Pola Office+ syringe and tips are for single use only.
Reuse will cause patient cross contamination.

STORAGE AND HANDLING
- Store product in the refrigerator 2 - 8 °C (35 - 46 °F).
- Use product at room temperature.
- Do not remove syringe from foil until ready to use – Pola Office+ is heat and light sensitive.
- Do not use if the foil appears to be bloated.

⚠ High concentration hydrogen peroxide (mixed gel concentration 30.4%). May cause tooth sensitivity, soft tissue irritation or chemical burns.

PRECAUTIONS
- For professional use only. Use only by trained dental professionals.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiration date.
- Pola Office+ syringes may be under pressure – remove cap and attach nozzle with caution
- Avoid contact with soft tissues, skin and eyes.
- Patient must wear protective eyewear.
- Operators must wear gloves and protective eyewear.
- Do not allow Pola Office+ gel to come in contact with skin, eyes and soft tissues.
- SDS available at sdi.com.au or contact your regional representative.

FIRST AID
Eye (contact): Open the eye wide and thoroughly wash for 15 minutes with running water. Seek medical attention.
Skin (contact): Immediately wash thoroughly with water. Seek medical attention if symptom persists.
Ingestion: Do not induce vomiting. Rinse mouth with water. Contact local poisons center.
Inhalation: No symptoms expected under normal conditions of use. If affected, move to well-ventilated area.

If a serious incident has occurred due to this product, please report this incident to the manufacturer of this product and to the competent authority in your state or country.

Dispose in accordance with local regulations.

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