



# AURA EASY

UNIVERSAL COMPOSITE



4

SHADES

SIMPLE,  
BRILLIANT  
RESTORATIONS,  
EVERY TIME!



# BEAUTIFUL RESTORATIONS JUST GOT EASIER

## EASY ELEGANT RESTORATIONS

Its unique blend of filler and resin work in harmony to deliver the strength and esthetics needed to provide long-lasting restorations.

## EASY HANDLING FOR EFFORTLESS PLACEMENT


Optimized, non-sticky consistency makes it simple to sculpt and pack.

## EASY SHADE MATCHING FOR MOST OF YOUR POSTERIOR AND ANTERIOR RESTORATIONS

Eliminate the guesswork from finding the right shade for your restorations by having to choose from just 4 universal shades.

With shades that are equidistantly placed on the color saturation spectrum and easily equate to a VITA® shade, Aura Easy makes it simpler for you to match shades for the majority of your daily needs.

Aura Easy simplifies your procedures, minimizes your inventory and makes it super easy for your practice.



Approximate equivalence to VITA® A1-D4® shade guide

## ACHIEVABLE ESTHETICS CASE STUDIES

**DR TERRY WONG** (AUSTRALIA)



**DR HUGH FLAX** (USA)



# IT WILL TRULY ALLOW YOU TO DELIVER MORE WHILE STOCKING LESS

## EASY TO STAND THE TEST OF TIME

With a 3mm depth of cure, high compressive and flexural strengths and the benefits of the proprietary low-shrinkage resin to offset marginal stress, you get restorations that are tough and robust.

### LOW VOLUMETRIC SHRINKAGE\*

Aura Easy delivers reduced stress to provide long-lasting restorations.

### EXCELLENT FLEXURAL STRENGTH\*

Aura Easy resists major deformations without fracturing.

### HIGH COMPRESSIVE STRENGTH\*

Aura Easy resists the high forces of compression in the mouth.

2.21%

115 MPa

375 MPa

\*SDI internal data



## ADVANTAGES

4 shades cover 16 VITA® A1-D4® shades<sup>†</sup>

Increased wear resistance and strength

Smooth high gloss finish every time

Non-sticky consistency makes it simple to sculpt

Low volumetric shrinkage to offset marginal stress

4 universal shades to minimize inventory

No slumping and holds its shape

Available in syringes and complets

<sup>†</sup>Internal test - 2021

## INDICATIONS

Anterior restorations

Posterior restorations

Indirect inlays, onlay and veneers

Core buildup

Splinting



Sandwich technique with glass-ionomer



» After a couple days with Aura Easy in my practice, it was easy to see that the four shades were all I needed. Aura Easy's handling and shade matching makes it easy to create quick and beautiful restorations. «

Dr LORI TROST DMD [USA]

# ORDER DETAILS



## SYRINGES

Syringe Refill $\text{\textcircled{a}}\text{e1}$ - 4g	8560020
Syringe Refill $\text{\textcircled{a}}\text{e2}$ - 4g	8560021
Syringe Refill $\text{\textcircled{a}}\text{e3}$ - 4g	8560022
Syringe Refill $\text{\textcircled{a}}\text{e4}$ - 4g	8560023

## COMPLETS

Complet Refill $\text{\textcircled{a}}\text{e1}$ - 0.25g	8565020
Complet Refill $\text{\textcircled{a}}\text{e2}$ - 0.25g	8565021
Complet Refill $\text{\textcircled{a}}\text{e3}$ - 0.25g	8565022
Complet Refill $\text{\textcircled{a}}\text{e4}$ - 0.25g	8565023

## AURA PRODUCT RANGE

Aura Dentin Chroma  
 Aura Enamel  
 Aura Bulk Fill  
 Aura Easy  
 Aura Easyflow



YOUR SMILE. OUR VISION.



02/2024

M100263 G



+DO36M1002631-

MADE IN AUSTRALIA  
 by SDI Limited  
 Bayswater, Victoria 3153  
 Australia  
[www.sdi.com.au](http://www.sdi.com.au)

**AUSTRALIA** 1800 337 003  
**AUSTRIA** 00800 0225 5734  
**BRAZIL** 0800 770 1735  
**FRANCE** 00800 0225 5734  
**GERMANY** 0800 100 5759

**ITALY** 00800 0225 5734  
**NEW ZEALAND** 0800 808 855  
**SPAIN** 00800 0225 5734  
**UNITED KINGDOM** 00800 0225 5734  
**USA & CANADA** 1 800 228 5166