



AURA EASY

UNIVERSAL COMPOSITE



4

SHADES

SIMPLE,
BRILLIANT
RESTORATIONS,
EVERY TIME!



BEAUTIFUL RESTORATIONS JUST GOT EASIER

EASY ELEGANT RESTORATIONS

Its unique blend of filler and resin work in harmony to deliver the strength and aesthetics needed to provide long-lasting restorations.

EASY HANDLING FOR EFFORTLESS PLACEMENT


Optimised, non-sticky consistency makes it simple to sculpt and pack.

EASY SHADE MATCHING FOR MOST OF YOUR POSTERIOR AND ANTERIOR RESTORATIONS

Eliminate the guesswork from finding the right shade for your restorations by having to choose from just 4 universal shades.

With shades that are equidistantly placed on the colour saturation spectrum and easily equate to a VITA® shade, Aura Easy makes it simpler for you to match shades for the majority of your daily needs.

Aura Easy simplifies your procedures, minimises your inventory and makes it super easy for your practice.



Approximate equivalence to VITA® A1-D4® shade guide

ACHIEVABLE AESTHETICS

CASE STUDIES

DR TERRY WONG (AUSTRALIA)



DR HUGH FLAX (USA)



IT WILL TRULY ALLOW YOU TO DELIVER MORE WHILE STOCKING LESS

EASY TO STAND THE TEST OF TIME

With a 3mm depth of cure, high compressive and flexural strengths and the benefits of the proprietary low-shrinkage resin to offset marginal stress, you get restorations that are tough and robust.

LOW VOLUMETRIC SHRINKAGE*

Aura Easy delivers reduced stress to provide long-lasting restorations.

EXCELLENT FLEXURAL STRENGTH*

Aura Easy resists major deformations without fracturing.

HIGH COMPRESSIVE STRENGTH*

Aura Easy resists the high forces of compression in the mouth.

2.21%

115MPa

375MPa

*SDI internal data



ADVANTAGES

4 shades cover 16 VITA® A1-D4® shades

Increased wear resistance and strength

Smooth high gloss finish every time

Non-sticky consistency makes it simple to sculpt

Low volumetric shrinkage to offset marginal stress

4 universal shades to minimise inventory

No slumping and holds its shape

Available in syringes and composites

INDICATIONS

Anterior restorations

Posterior restorations

Indirect inlays, onlay and veneers

Core buildup

Splinting



Sandwich technique with glass-ionomer



» After a couple days with Aura Easy in my practice, it was easy to see that the four shades were all I needed. Aura Easy's handling and shade matching makes it easy to create quick and beautiful restorations. «

Dr LORI TROST DMD

ORDER DETAILS



SYRINGES

Aura Easy Syringe Kit 4x 4g syringe 1 x Aura Easy shade guide	8565112
Syringe Refill ae1 - 4g	8560020
Syringe Refill ae2 - 4g	8560021
Syringe Refill ae3 - 4g	8560022
Syringe Refill ae4 - 4g	8560023

COMPLETS

Aura Easy Complet Kit 20x4 0.25g complets 1 x Aura Easy shade guide	8565113
Complet Refill ae1 - 0.25g	8565020
Complet Refill ae2 - 0.25g	8565021
Complet Refill ae3 - 0.25g	8565022
Complet Refill ae4 - 0.25g	8565023

AURA PRODUCT RANGE

Aura Dentine Chroma
 Aura Enamel
 Aura Bulk Fill
 Aura Easy
 Aura Easyflow



YOUR SMILE. OUR VISION.



11-2021

M100262 G



+D036M1002621Z

MADE IN AUSTRALIA
 by SDI Limited
 Bayswater, Victoria 3153
 Australia
www.sdi.com.au

AUSTRALIA 1800 337 003
AUSTRIA 00800 0225 5734
BRAZIL 0800 770 1735
FRANCE 00800 0225 5734
GERMANY 0800 100 5759

ITALY 00800 0225 5734
NEW ZEALAND 0800 808 855
SPAIN 00800 0225 5734
UNITED KINGDOM 00800 0225 5734
USA & CANADA 1 800 228 5166