

NOW **eASIER** THAN EVER

4 SHADES =

Simple, brilliant restorations, every time!



Your Smile. Our Vision.
www.sdi.com.au



aura eASY

ultra universal restorative material


Beautiful restorations

eASY handling for effortless placement

Optimised, non-sticky consistency makes it simple to sculpt and pack.

eASY elegant restorations

Its unique blend of filler and resin work in harmony to deliver the strength and aesthetics needed to provide long-lasting restorations.



eASY to stand the test of time

With a 3mm depth of cure, high compressive and flexural strengths and the benefits of the proprietary low-shrinkage resin to offset marginal stress, you get restorations that are tough and robust.

2.21%

Low volumetric shrinkage

aura eASY delivers reduced stress to provide long-lasting restorations.

115
MPa

Excellent flexural strength

aura eASY resists major deformations without fracturing.

375
MPa

High compressive strength

aura eASY resists the high forces of compression in the mouth.



“After a couple days with **aura** in my practice, it was easy to see that the four shades were all I needed. **aura's** handling and shade matching makes it easy to create quick and beautiful restorations.

Dr Lori Trost, DMD

Restorations just got easier

eASY shade matching for most of your posterior and anterior restorations

Eliminate the guesswork from finding the right shade for your restorations by having to choose from just 4 universal shades.

With shades that are equidistantly placed on the colour saturation spectrum and easily equate to a Vita shade, **aura eASY** makes it simpler for you to customise shades for the majority of your daily needs.

aura eASY simplifies your procedures, minimises your inventory and makes it super easy for your practice.



It will truly allow you to deliver more while stocking less.

aura eASY it's the easy choice for your practice

- Increased wear resistance and strength
- Smooth high gloss finish every time
- Non-sticky consistency makes it simple to sculpt
- No slumping and holds its shape
- Quick and easy placement in seconds

A1/B1 ≈ ae1
A2/B2 ≈ ae2
A3/B3 ≈ ae3
A3.5/B4 ≈ ae4

Approximate equivalence to VITA classical shade guide



Dr Terry Wong (Australia)



Dr Hugh Flax (USA)

Ordering information

PART NUMBER	QTY	SYRINGES
8565112	4x1	aura eASY Syringe Kit 4g syringe (ae1, ae2, ae3, ae4)
	1	aura eASY shade guide
		aura eASY Syringe Refills
8560020	1	ae1 – 4g
8560021	1	ae2 – 4g
8560022	1	ae3 – 4g
8560023	1	ae4 – 4g



PART NUMBER	QTY	COMPLETS
8565113	20x4	aura eASY Complet Kit 0.25g complets (ae1, ae2, ae3, ae4)
	1	aura eASY shade guide
		aura eASY Complet Refills
8565020	20	ae1 – 0.25g
8565021	20	ae2 – 0.25g
8565022	20	ae3 – 0.25g
8565023	20	ae4 – 0.25g



Your Smile. Our Vision.
www.sdi.com.au

Australia 1 800 337 003
Austria 00800 0225 5734
Brazil 0800 770 1735
France 00800 0225 5734
Germany 0800 1005759
Ireland 00800 0225 5734

Italy 00800 0225 5734
New Zealand 0800 734 034
Spain 00800 0225 5734
United Kingdom 00800 0225 5734
USA & Canada 1 800 228 5166

09-01-2017 M100262 B



+D036M1002621Z