

SDI | RIVA GUIDE



NOW WITH
EASY
ACTIVATION
CAPSULE

RIVA SELF CURE + RIVA SELF CURE HV

BULK FILL GLASS HYBRID
RESTORATIVE SYSTEM

- *ionglass*TM technology
- Strong dentin replacement
- High fluoride release
- Chemically adheres to tooth structure
- BPA and HEMA free
- Suitable for bulk filling
- No shrinkage

Special features for Riva Self Cure High Viscosity

- Very easy to shape and contour
- Non sticky

Riva Self Cure compares to:

Fuji IX GP* Fuji IX GP Extra* Ketac-Molar* Equia Forte* Equia Forte HT*

Riva Self Cure:

Boxes of 50 capsules
Powder / liquid kits

Shades: A1, A2, A3, A3.5, A4
Regular Set and Fast Set

Riva Self Cure HV:

Boxes of 50 capsules
Shades: A1, A2, A3, A3.5



light safe tip!

NOW WITH
EASY
ACTIVATION
CAPSULE

RIVA LIGHT CURE + RIVA LIGHT CURE HV

LIGHT CURED, RESIN REINFORCED GLASS
IONOMER RESTORATIVE MATERIAL

- *ionglass*TM technology
- Long term wear resistance due to high surface hardness
- Chemically adheres to tooth structure
- BPA free
- High fluoride release
- Orange light safe tip prevents premature setting
- Fabulous aesthetics
- Highest colour stability to tea and coffee

Special features for Riva Light Cure High Viscosity

- Easy to manipulate and sculpt

Riva Light Cure compares to:

Fuji II LC* Photac-Fil Quick*

Riva Light Cure:

Boxes of 50 capsules
Powder / liquid kits

Shades: A1, A2, A3, A3.5

Riva Light Cure HV:

Boxes of 50 capsules

Shades: A1, A2, A3, A3.5

RIVA LIGHT CURE HV



RIVA SILVER

MERCURY FREE SELF CURING CONVENTIONAL
GLASS IONOMER RESTORATIVE MATERIAL

- *ionglass*TM technology
- Silver reinforced
- Very high radiopacity
- Chemically adheres to tooth structure
- Easy handling
- No shrinkage
- Great for core build ups
- High early strength - perfect for children and challenging patients
- High fluoride releasing and recharging ability
- High wear resistance

Riva Silver compares to:

Miracle Mix* Ketac Silver*

Riva Silver:

Boxes of 50 capsules
Powder / liquid kits
Powder refills

Shade: Light grey

*Not a registered trademark of SDI Limited.



SDI | YOUR SMILE. OUR VISION.



+0036M861001/1

MADE IN AUSTRALIA
by SDI Limited
Bayswater, Victoria 3153, Australia
www.sdi.com.au
AUSTRALIA 1800 337 003

AUSTRIA 00800 0225 5734
BRAZIL 0800 770 1735
FRANCE 00800 0225 5734
GERMANY 0800 100 5759
IRELAND 00800 0225 5734

ITALY 00 353 (83) 806 0228
NEW ZEALAND 0800 808 855
SPAIN 00800 0225 5734
UNITED KINGDOM 00800 0225 5734
USA & CANADA 1 800 228 5166



RIVA GUIDE



RIVA PROTECT

SELF CURING CONVENTIONAL GLASS IONOMER USED AS A SEALANT, LINER AND SURFACE PROTECTOR

- *ionglass™* technology
- BPA and HEMA free
- Works in a moist environment
- Perfect for caries-challenged patients
- Great protector for partially erupted teeth
- ACP nanotechnology - supercharges remineralisation
- Super high fluoride releasing and recharging ability
- Chemically adheres to tooth structure
- Assists in arresting to tooth structure

Riva Protect compares to:
Fuji VII* Fuji Triage*

Riva Protect:
Boxes of 50 capsules
Shades: White and Pink
Regular Set and Fast Set



RIVA LUTING

SELF CURING CONVENTIONAL GLASS IONOMER LUTING CEMENT

- *ionglass™* technology
- Chemically adheres to tooth structure and metal substrates, such as stainless steel crowns, PFM crowns and orthodontic bands
- BPA and HEMA free
- Very low film thickness - flows easily between restoration and tooth
- Easy to clean up
- Low water solubility - increased longevity and aesthetics

Riva Luting compares to:
Fuji I* Ketac Cem*

Riva Luting:
Boxes of 50 capsules
Powder / liquid kits
Powder refills
Liquid refills



RIVA LUTING PLUS

SELF CURING RESIN MODIFIED GLASS IONOMER LUTING CEMENT

- *ionglass™* technology
- Chemically adheres to tooth and substrates
- BPA free
- High fluoride releasing and recharging ability
- High flexural strength - reduced cracking and failure in flexion
- Clinically insoluble
- High acid resistance
- High bond strengths
- Bioactive proprietary *ionglass™* technology - high ion exchange

Riva Luting Plus compares to:
Fuji Plus* Relyx Luting Cement*

Riva Luting Plus:
Boxes of 50 capsules
Powder / liquid kits
Powder refills
Liquid refills

*Not a registered trademark of SDI Limited.



RIVA APPLICATOR PRO

EASY. STRONG. CLEAN.

- Compatible with SDI and other market leading glass ionomer capsules¹
- Contoured design for strength and comfort
- Engineered stainless steel mechanical parts and durable mechanical casing
- Autoclavable and ultrasonic bath tolerant²
- No lubrication required

RIVA COAT

- Prevents saliva and water contamination to exposed surfaces of glass ionomer cements during the first 24 hours immediately after placement.
- Useful for sealing composites & compomers
- 5mL bottle

RIVA CONDITIONER

- Partially removes the smear layer
- Optimises adhesion of the GIC
- 22% polyacrylic acid
- 10mL bottle

^{1,2} Internal SDI test data