

# SDI | RIVA GUIDE



NOW WITH  
EASY  
ACTIVATION  
CAPSULE

## RIVA SELF CURE + RIVA SELF CURE HV

BULK FILL GLASS HYBRID  
RESTORATIVE SYSTEM

- *ionglass™* technology
- Strong dentin replacement
- High fluoride release
- Chemically adheres to tooth structure
- BPA and HEMA free
- Suitable for bulk filling
- No shrinkage

### Special features for Riva Self Cure High Viscosity

- Very easy to shape and contour
- Non sticky

**Riva Self Cure compares to:**  
Fuji IX GP\* Fuji IX GP Extra\* Ketac-  
Molar\* Equia Forte\* Equia Forte HT\*

### Riva Self Cure:

Boxes of 50 capsules  
Powder / liquid kits

Shades: A1, A2, A3, A3.5, A4  
Regular Set and Fast Set

### Riva Self Cure HV:

Boxes of 50 capsules  
Shades: A1, A2, A3, A3.5



light safe tip!

NOW WITH  
EASY  
ACTIVATION  
CAPSULE

## RIVA LIGHT CURE + RIVA LIGHT CURE HV

LIGHT CURED, RESIN REINFORCED GLASS  
IONOMER RESTORATIVE MATERIAL

- *ionglass™* technology
- Long term wear resistance due to high surface hardness
- Chemically adheres to tooth structure
- BPA free
- High fluoride release
- Orange light safe tip prevents premature setting
- Fabulous aesthetics
- Highest colour stability to tea and coffee

### Special features for Riva Light Cure High Viscosity

- Easy to manipulate and sculpt

**Riva Light Cure compares to:**  
Fuji II LC\* Photac-Fit Quick\*

### Riva Light Cure:

Boxes of 50 capsules  
Powder / liquid kits

Shades: A1, A2, A3, A3.5

### Riva Light Cure HV:

Boxes of 50 capsules

Shades: A1, A2, A3, A3.5

RIVA LIGHT CURE HV



## RIVA SILVER

MERCURY FREE SELF CURING CONVENTIONAL  
GLASS IONOMER RESTORATIVE MATERIAL

- *ionglass™* technology
- Silver reinforced
- Very high radiopacity
- Chemically adheres to tooth structure
- Easy handling
- No shrinkage
- Great for core build ups
- High early strength - perfect for children and challenging patients
- High fluoride releasing and recharging ability
- High wear resistance

### Riva Silver compares to:

Miracle Mix\* Ketac Silver\*

### Riva Silver:

Boxes of 50 capsules  
Powder / liquid kits  
Powder refills

Shade: Light grey

\*Not a registered trademark of SDI Limited.

SDI | YOUR SMILE. OUR VISION.

01/2024

M861002 H



\*D038M8610021\*

MADE IN AUSTRALIA  
by SDI Limited  
Bayswater, Victoria 3153  
Australia  
[www.sdi.com.au](http://www.sdi.com.au)

**AUSTRALIA** 1800 337 003  
**AUSTRIA** 00800 0225 5734  
**BRAZIL** 0800 770 1735  
**FRANCE** 00800 0225 5734  
**GERMANY** 0800 100 5759

**ITALY** 00800 0225 5734  
**NEW ZEALAND** 0800 808 855  
**SPAIN** 00800 0225 5734  
**UNITED KINGDOM** 00800 0225 5734  
**USA & CANADA** 1 800 228 5166

# SDI | RIVA GUIDE



## RIVA PROTECT

SELF CURING CONVENTIONAL GLASS IONOMER USED AS A SEALANT, LINER AND SURFACE PROTECTOR

- *ionglass™* technology
- BPA and HEMA free
- Works in a moist environment
- Perfect for caries-challenged patients
- Great protector for partially erupted teeth
- ACP nanotechnology - supercharges remineralisation
- Super high fluoride releasing and recharging ability
- Chemically adheres to tooth structure
- Assists in arresting to tooth structure

**Riva Protect compares to:**  
Fuji VII\* Fuji Triage\*

### Riva Protect:

Boxes of 50 capsules

Shades: White and Pink

Regular Set and Fast Set



## RIVA LUTING

SELF CURING CONVENTIONAL GLASS IONOMER LUTING CEMENT

- *ionglass™* technology
- Chemically adheres to tooth structure and metal substrates, such as stainless steel crowns, PFM crowns and orthodontic bands
- BPA and HEMA free
- Very low film thickness - flows easily between restoration and tooth
- Easy to clean up
- Low water solubility - increased longevity and aesthetics

**Riva Luting compares to:**  
Fuji I\* Ketac Cem\*

### Riva Luting:

Boxes of 50 capsules  
Powder / liquid kits

2021  
**REALITY**  
MATTERS & REVIEWS  
Four Star Award  
★★★★



## RIVA LUTING PLUS

SELF CURING RESIN MODIFIED GLASS IONOMER LUTING CEMENT

- *ionglass™* technology
- Chemically adheres to tooth and substrates
- BPA free
- High fluoride releasing and recharging ability
- High flexural strength - reduced cracking and failure in flexion
- Clinically insoluble
- High acid resistance
- High bond strengths
- Bioactive proprietary *ionglass™* technology - high ion exchange

**Riva Luting Plus compares to:**  
Fuji Plus\* Relyx Luting Cement\*

### Riva Luting Plus:

Boxes of 50 capsules  
Powder / liquid kits



## RIVA APPLICATOR 2

SLEEK AND SOPHISTICATED

- Fits comfortably in smaller hands
- Autoclavable - easy to clean
- No sharp edges - hand contoured design for maximum comfort
- Capsule lock provides controlled extrusion
- Robust & long lasting
- Easy oral access
- Superbly engineered stainless steel mechanical parts



## RIVA COAT

- Prevents saliva and water contamination to exposed surfaces of glass ionomer cements during the first 24 hours immediately after placement.
- Useful for sealing composites & compomers
- 5mL bottle

## RIVA CONDITIONER

- Partially removes the smear layer
- Optimizes adhesion of the GIC
- 22% polyacrylic acid
- 10mL bottle