



SDI | RADII PLUS ATTACHMENTS

YOUR CHOICE OF ACCESSORIES

FULL ARCH BLEACH LED ATTACHMENT

FOR IN-OFFICE BLEACHING OF FULL ARCH

DIAGNOSTIC LED ATTACHMENT

SEE YOUR PATIENTS IN A DIFFERENT LIGHT!

ORTHODONTIC CURING LED ATTACHMENT

FAST AND EFFICIENT BONDING OF BRACKETS

CONTRA ANGLE CURING LED ATTACHMENT

A DECISIVE STEP FORWARD TO YOUR CURING NEEDS

RIGHT ANGLE CURING LED ATTACHMENT

GREAT FOR OCCLUSAL CURING

Please note: Your Radium Plus "+" LED curing light will work only with the new Radium Plus "+" attachments.

Look for the small "+" sign on the Control section which will indicate the improved Radium Plus.

FULL ARCH BLEACH LED ATTACHMENT

FOR IN-OFFICE BLEACHING

- Beeps every minute
- Heats to 37°C / 98.76°F



order: 5600092

POLA WHITENING STAND

SMART DESIGN & LIGHTWEIGHT

- Extremely stable and easy to manoeuvre



order: 7700989

CONTRA ANGLE LED ATTACHMENT

A DECISIVE STEP FORWARD TO YOUR CURING NEEDS

- Just like your handpiece
- Permits access to areas difficult to cure



order: 5600276

RIGHT ANGLE LED ATTACHMENT

CURING TIP

- Ease of access
- Great for occlusal curing



order: 5600260

ORTHODONTIC LED ATTACHMENT

FAST AND EFFICIENT BONDING OF BRACKETS

- Beeps every 10 seconds
- Runs continuously for 5 minutes



order: 5600259

DIAGNOSTIC LED ATTACHMENT

EXAMINATION USING TRANSILLUMINATION

- Lights up those hard to see areas to assist in identification of cracks & decay



order: 5600268

LED RADIOMETER

INSTANT DIGITAL INTENSITY READING

- Specially calibrated for any LED curing light



order: 5600028

CURING TACKING TIP

GREAT FOR SEATING OF VENEERS

- Place over the lens cap when curing
- 1mm and 2mm openings
- Autoclavable



order: 5600079

DIAGNOSTIC ADAPTER TIP

PRODUCES A FINE BEAM OF TRANSILLUMINATING LIGHT

- Replace lens cap with this adapter tip when using the diagnostic attachment
- Autoclavable



order: 5600078