

**BULK FILL GLASS HYBRID RESTORATIVE SYSTEM** 



# THE IONGLASS LONG TERM BREAKTHROUGH



Drawing on 50 years of dental research experience, SDI's *ionglass™* filler was developed specially for Riva and is manufactured in Australia by our glass experts.

The *ionglass*™ technology mimics the natural tooth structure, guaranteeing restorations will withstand long term mastication forces and patients will be satisfied with the treatment.

### **BIOACTIVE GLASS HYBRID**

Riva Self Cure contains *ionglass*™ filler, made of a unique blend of different sizes of ultrafine highly reactive glass particles.

This bioactive proprietary hybrid glass enables long term restorations, being outstanding in handling, strength and aesthetics.

The *ionglass*™ filler contains fluoride and strontium ions which synergistically enhance tooth bioremineralisation.


(Source: Thuy et al. Arch Oral Biol. 2008 Nov; 53(11):1017-22)

Superior biocompatibility ensures healthier teeth, the restoration will last longer and the patient will be more satisfied.

### SUPERIOR FLUORIDE RELEASE

The extra strong hybrid matrix provides great fluoride release<sup>†</sup>, increasing defense against caries, while having very low wear to enable the perfect long term bulk fill restoration every time.

### CUMULATIVE FLUORIDE RELEASE FOR RIVA SELF CURE



<sup>†</sup> Source: Internal test data Dec 2020

### RIVA SELF CURE IS BIOMIMETIC


Riva Self Cure replaces dentine, imitating some of its natural features. The *ionglass™* technology for Riva Self Cure mimics dentine's hydrophilic characteristic and coefficient of thermal expansion. To make the restoration time-efficient, no etching or adhesive is required. This dramatically reduces chair time for the dentist and the risk of sensitivity for the patient. A patient with a nice smile and no sensitivity has a better quality of life.

### LONG TERM ACID RESISTANCE

The *ionglass*™ technology stimulates the biomineralisation of dentine by over saturating the minerals in its organic matrix. The result is the formation of fluorapatite crystals, which are more resistant to acid attacks than the natural tooth, benefiting patients during bacterial and acidic challenges.

The ISO Standard 9917-1 dictates that a material cannot have acid erosion of more than 0.17mm. Riva Self Cure will resist acidic oral environments.

### DEPTH OF ACID EROSION AFTER 24 HOURS



Source: Internal test data Nov 2020

\*Not a registered trademark of SDI

# BEAUTIFUL HANDLING

Riva Self Cure extrudes easily into the cavity with the innovative *ionglass™* particles enabling creamier, non-stick and smoother handling. Polishability is made simple and long-term surface integrity is maintained. For an exceptional packable, high viscosity material, use Riva Self Cure HV.

Riva Self Cure has been independently evaluated as having superior handling. The unique formulation is less stringy or sticky, making placement quicker and easier.








Equia Forte HT



Riva SC



Riva SC HV (SDI)

Source: Internal test data Nov 2020 \*Not a registered trademark of SDI

 $\longrightarrow$  Riva Self Cure has the edge on glass strength and fluoride release.  $<\!<$ 

Professor Edward Lynch, University of Warwick, UK Listed in the World's Top 100 Doctors in Dentistry - January 2020. Winner of The Global Outreach Healthcare Award - June 2020.

# EASY ACTIVATION

The new Easy Activation capsules are now 70% easier to activate. This saves you time and effort when preparing a restoration.




### EASIER EXTRUSION

Working Time (s)

Riva Self Cure is 70% easier to extrude<sup>†</sup>. This minimises hand tremors and muscle stress, giving you more stability and control when placing restorations.


# CHOICE OF SETTING TIME FOR ALL TREATMENTS

Different clinical scenarios demand different setting times. Riva Self Cure Fast Set offers the fastest setting time for when a guick procedure is required. If more time is needed, Riva Self Cure offers the ideal time-frame for working with comfort and confidence.



Source: Internal test data Nov 2020 Not a registered trademark of SDI

# AVERAGE EXTRUSION FORCE



<sup>†</sup>Source: Internal test data Nov 2020 \*Not a registered trademark of SDI

### VISCOSITY OPTIONS


Riva Self Cure is available in a range of viscosity and transparency options to suit your placement preferences.



# BEST IN CLASS AESTHETICS

### OPTIMAL TRANSLUCENCY

Riva Self Cure is formulated to give the best aesthetic result<sup>†</sup>. The T-shade range provides higher translucency and better optical properties. This chameleon effect effortlessly blends with surrounding tooth structure and dentine.




 $^\dagger$ Source: Internal test data, Nov 2020

\*Not a registered trademark of SDI.

### HIGHEST RADIOPACITY

Riva Self Cure *ionglass*<sup>TM</sup> offers significantly higher radiopacity than other restorative options of its category<sup>†</sup>. On radiographs, the contrast between Riva Self Cure and dentine is evident on the tooth-restoration interface. This helps to diagnose and treat patients faster and more easily.

### RADIOPACITY



>>> Versatility and easy handling are key benefits of Riva Self Cure for use in Paediatric Dentistry. <<

Dr Daniela P. Raggio,
Associate Professor - University of Sao Paulo
Editor in Chief International Journal of Paediatric Dentistry (UK)





Photos courtesy of Dr. Marcelo Alves and Dr. Alvaro Junqueira - Brazil





Photos courtesy of Dr. Marcelo Alves and Dr. Alvaro Junqueira - Brazil

# A HYBRID RESTORATIVE WITH MORE INDICATIONS







Class I

Class I

- Class II (see technique guide on SDI website)
- Class V
- Core build up
- Deciduous teeth restorations
- Geriatric restorations
- Temporary fillings
- ✓ ART (Atraumatic Restorative Treatment)
- ✓ SMART (Silver Modified Atraumatic Restorative Treatment) e.g. Riva Star, SDI Limited.
- Dentine replacement
- Sandwich technique restorations
- Molar Incisal hypomineralisation [MIH] (intermediate approach)
- Minimally invasive dentistry

# THE UNIQUE BIOMIMETIC FOR LONG TERM RESTORATIONS

Formulated using patented *ionglass*™ technology, Riva Self Cure has restored the smiles of over 50 million people across 100 countries.

The next generation of glass ionomer, this high performance system combines a glass hybrid restorative with a strong resin coat to maximise the longevity and the strength of your patient's restorations. Riva Self Cure has earned its reputation through superior fluoride release, strength and great handling.

Perfect for a wide range of clinical situations, Riva Self Cure is the bulk fill solution you need to save you time and to mimic some of the best characteristics of the tooth.



LONG TERM STRENGTH FOR BULK FILLS Clinicians Report Jan 2020, Vol 13, Issue 1	370 <sub>MPa</sub>
SUPERIOR WEAR RESISTANCE Internal test data, November 2020	$99_{\mu m}$
SUPERIOR LONG TERM FLUORIDE RELEASE Internal test data, November 2020	3.01 <sub>µg/cm²</sub>
LONG TERM ACID RESISTANCE Internal test data, November 2020	$0.07_{\text{mm}}$
HIGH RADIOPACITY	290%


RADIOPACITY Internal test data. November 2020

# STRONGER FOR LONG TERM

Like dentine, Riva Self Cure has a very high compressive strength and a superior wear resistance, ensuring it will withstand long term mastication forces.


### HIGH STRENGTH

### COMPRESSIVE STRENGTH



Source: An independent, non-profit, dental education and product testing foundation, Clinicians Report, Jan 2020. Full report available from SDI upon request Not a registered trademark of SDI

# SUPERIOR WEAR RESISTANCE




Source: Internal test data Nov 2020 \*Not a registered trademark of SDI. †Tested with ACTA (Academic Centre for Dentistry Amsterdam) wear machine proven to be a realistic simulation of clinical wear <sup>1</sup>de Gee, A. and Pallay, P., 1994. Occlusal wear simulation with the ACTA wear

machine. Journal of Dentistry, 22, pp.S21-S27.

# INSTRUCTIONS<sup>^</sup>







Activate the capsule and mix it in a triturator for 10 seconds, before placing it into the applicator.







 For powder and liquid instructions, please refer to Instructions for Use

 ${\boldsymbol \smallfrown}$  For powder and liquid instructions, please refer to Instructions for Use



RIVA SELF CURE CAPSULES	
Riva Self Cure Regular Set 50 capsules	
A1	8600001
A2	8600002
A3	8600003
A3.5	8600004
A4	8600005
B2	8600007
В3	8600008
Riva Self Cure Fast Set 50 capsules	
A1	8605001

A1	8605001
A2	8605002
A3	8605003
A3.5	8605004
A4	8605005
B2	8605007
B3	8605008
T-A2	860 0011
T-A3	8600012
T-A3.5	8600013

### Riva Self Cure Assorted Kits 10 capsules each: A1, A2, A3, A3.5 and B2 $\,$

Regular Set	8610000
Fast Set	8620000

### Riva Self Cure High Viscosity

50 capsules

A1	8630001
A2	8630002
A3	8630003
A3.5	8630004



#### RIVA SELF CURE POWDER/LIQUID KITS

6.9mL (8g) Riva Self Cure Liquid
15g Riva Self Cure Powder
Accessories

A1	8610501
A2	8610502
A3	8610503
A3.5	8610504

#### 9.5mL (11g) Riva Self Cure Liquid 15g Riva Self Cure Powder Accessories

T-A2	8610521
T-A3	8610522
T-A3.5	8610523

# 5.2mL (6g) Riva Self Cure Liquid bottle 10g Riva Self Cure Powder bottle Accessories

A2	8610612
A3	8610613



### POWDER

Riva Self Cure Powder 15g Riva Self Cure Powder

8610103 АЗ



# RIVA SELF CURE LIQUID

6.9mL (8g) Bottle for Regular Set shades	8610900
9.5ml (11g) Bottle for T-shades	8610905



AUULUUUNILU	
Riva Coat	8610001
5mL bottle refill	
Riva Conditioner	8620001
10mL bottle refill	



#### **APPLICATORS**

Riva Applicator	5545009
Riva Applicator 2	5545013



### **EQUIPMENT**

Illtramat 2	
Utti ai i i at L	

High speed multi-use triturator

### Ultramat S

High speed multi-use triturator



**(€** 0297



MADE IN AUSTRALIA by SDI Limited Bayswater, Victoria 3153 Australia www.sdi.com.au

**AUSTRALIA** 1800 337 003 **AUSTRIA** 00800 0225 5734 **BRAZIL** 0800 770 1735 FRANCE 00800 0225 5734 **GERMANY** 0800 100 5759

**ITALY** 00800 0225 5734 **NEW ZEALAND** +64 21 400 323 **SPAIN** 00800 0225 5734 **UNITED KINGDOM** 00800 0225 5734 **USA & CANADA** 1 800 228 5166