

STELA

THE FUTURE OF COMPOSITES

- ✓ Unlimited depth of cure (self cure)
- ✓ Gap-free interface
- ✓ Only 2 steps: 15 sec prep
- ✓ High strength
- ✓ Choose: capsule or automix




25%* OFF

STELA AUTOMIX

AUTOMIX REFILL KIT - **NOW \$150.61**
AUTOMIX INTRO KIT - **NOW \$182.43**



 Images courtesy Dr Alex Olivaldo (Brazil)



The application of Stela is very easy and after 4 minutes the restoration had the same shade as the tooth. Great polishability!

PROF DR ROCIO LAZO - [PERU]



STELA CAPSULE REFILL 50 CAPS	MRP \$291.98
STELA CAPSULES INTRO KIT	MRP \$145.54

TAKE-HOME WHITENING



POLA DAY

HYDROGEN PEROXIDE BASED ADVANCED TRAY TOOTH WHITENING SYSTEM

- Available in **6% 7.5% and 9.5% hydrogen peroxide** applications
- Fast acting : from 30 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximise patient comfort



POLA NIGHT

CARBAMIDE PEROXIDE BASED ADVANCED TRAY TOOTH WHITENING SYSTEM

- Available in **10% 16% 18% and 22% carbamide peroxide** applications
- From 45 minutes to overnight wear
- High viscosity, neutral pH gel contains Potassium Nitrate to maximise patient comfort



POLA DAY CP

SUPER FAST TRAY TOOTH WHITENING SYSTEM!

- Available in **35% carbamide peroxide** applications
- Super fast : Just 15 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximise patient comfort



POLA FOR ALIGNERS

STRAIGHTEN AND BRIGHTEN TOGETHER

- Available in **10% carbamide peroxide**
- The greatest patient comfort in a whitening kit
- Contains Potassium Nitrate and fluoride to reduce sensitivity for long term use



POLA LIGHT

COMPLETE YOUR WHITENING OFFERING WITH POLA LIGHT

- Available in **6% and 9.5% hydrogen peroxide**
- No custom trays required
- Reduced patient consult time
- Entry level whitening option
- Great front desk conversation starter



DENTIST PRICE **MRP \$114.61**
RECOMMENDED PATIENT PRICE \$189.00

POLA LUMINATE

BRUSH ON & GO!

- Available in **6% hydrogen peroxide**, low viscosity fast drying gel
- No trays required - quick and easy professional brush on tooth whitener
- Apply twice a day for just 30 minutes



DENTIST PRICE **MRP \$28.80**
RECOMMENDED PATIENT PRICE \$75.00



BUY 4
GET 1
FREE^x

POLA 4 X 1.3g MINI KIT	MRP \$41.47
POLA 10 X 1.3g SYRINGE KIT	MRP \$101.57
POLA 50 X 3g DISPENSER PACK	MRP \$564.30
POLA FOR ALIGNERS	MRP \$87.10

IN-OFFICE WHITENING



POLA RAPID

Pola Rapid is an in office whitening treatment that safely whitens teeth in a 24 minute application time.

SUPERIOR HANDLING AND FAST REMOVAL

Pola Rapid's higher viscosity gel enables faster application and removal due to its specially formulated non-stick attributes.

Pola Rapid can be extruded quickly and easily and remains on teeth when applied.

ENHANCED BLUE GEL FOR FASTER APPLICATION

Pola Rapid uses an advanced blue gel for ease of visibility and application. The easy to apply non-stick gel enables faster application and removal.

MINIMAL CHAIR TIME WITH IMMEDIATE RESULTS

A unique **38% hydrogen peroxide** formulation that releases active ingredients immediately. In only one visit, Pola Rapid initiates a powerful whitening treatment designed to whiten teeth in 24 minutes. No pre-desensitiser step is required.



SUPER. FAST.
WHITENING.

POLA RAPID 3 PATIENT KIT	MRP \$473.78
POLA RAPID 3 PATIENT KIT WITH OPTRAGATE	MRP \$504.63
POLA RAPID 1 PATIENT KIT	MRP \$167.93
POLA RAPID 1 PATIENT KIT WITH OPTRAGATE	MRP \$179.67



CASE STUDIES




Dr Miles Cone
USA



BEFORE



AFTER



Dr Chanel Beydoun
Australia



BEFORE



AFTER

COMPOSITE AND GIC PROMOTIONS

only available if ordered direct from SDI

It's never been more important to

Buy Australian

than right now

australianmade.com.au



COMPOSITES

LUNA 2 | AURA EASY | AURA BULK FILL

\$70[×]

PER
SYRINGE

\$80[×]

PER COMPLET
REFILLS

MUST BE PURCHASED
IN MULTIPLES OF 10




LUNA 2

UNIVERSAL COMPOSITE

- Improved non-stick, non-slump handling and reliable mechanical properties
- BPA free
- Very high radiopacity of 250% Al

LUNA 2 4g SYR	MRP \$117.70
LUNA 2 20 X 0.25g Complets	MRP \$152.83
LUNA 2 INTRO KIT	MRP \$572.54



 Images courtesy of Dr Fady Toussoun (Australia)




AURA EASY

EASIER AND FASTER AESTHETICS FOR ANTERIOR AND POSTERIOR

- 4 shades cover 16 VITA classical A1-D4® shades
- Increased wear resistance and strength
- Smooth high gloss finish every time
- Non-sticky consistency makes it simple to sculpt
- Low volumetric shrinkage to offset marginal stress

AURA EASY 4g Syr	MRP \$123.16
AURA EASY 20 X 0.25g Complets	MRP \$138.55



 Images courtesy of Dr Terry Wong (Australia)



AURA BULK FILL

SPEEDY AND SIMPLE POSTERIOR COMPOSITE RESTORATIONS

- Natural-looking finish with just one universal shade
- High depth of cure without the need for a capping layer
- Easy placement in a matter of seconds

AURA BULK FILL 4g SYR REFILLS	MRP \$123.87
AURA BULK FILL 20 X 0.25g Complets	MRP \$139.35



 Images courtesy of Dr Bob Lowe (USA)

GLASS IONOMERS

BUY 2 BOXES OF 50 CAPSULES GET A RIVA APPLICATOR PRO

ALL RIVA AND STELA - ANY MIX OF 50 CAPSULE BOXES

FREE[×]



SAVE
\$370.66*



RIVA SELF CURE/HV

BULK FILL GLASS HYBRID RESTORATIVE SYSTEM

- Contains SDI's proprietary **ionglass™**, a high ion releasing, radiopaque, bioactive glass that remineralises and prevents secondary caries formation
- Very high fluoride releasing and recharging ability

RIVA SELF CURE	MRP \$270.15
----------------	--------------



RIVA LIGHT CURE/HV

LIGHT CURED RESIN REINFORCED GLASS IONOMER RESTORATIVE MATERIAL

- Contains SDI's proprietary **ionglass™**, a high ion releasing, radiopaque, bioactive glass that remineralises and fights against secondary caries formation

RIVA LIGHT CURE	MRP \$301.70
-----------------	--------------



RIVA LUTING

- Chemically adheres to tooth structure and metal substrates
- Very low film thickness - flows easily between restoration and tooth
- Exceptional bond strengths
- BPA & HEMA free

RIVA LUTING	MRP \$289.39
-------------	--------------



RIVA LUTING PLUS

- Chemically adheres to tooth and substrates
- Exceptional bond strengths
- High flexural strength - reduced cracking and failure in flexion

RIVA LUTING PLUS	MRP \$308.78
------------------	--------------



RIVA PROTECT

- ACP nanotechnology - supercharges remineralisation
- Very high fluoride releasing and recharging ability to protect teeth from acid attacks

RIVA PROTECT	MRP \$217.45
--------------	--------------



RIVA SILVER

- Silver reinforced
- Ideal for abfraction lesion restorations
- Ideal for bulk filling
- No shrinkage

RIVA SILVER	MRP \$478.70
-------------	--------------



A NEXT GENERATION WATER BASED SOLUTION TO ARREST CARIES



CASE STUDY:

The presented case involves a 19 year old male who attended the practice for treatment. He had recently finished high school and began a part time job. Finances were an issue. However, he was motivated to improve his oral health and would like to maintain his teeth for the long term.

On examination there were several teeth with primary and secondary caries. 16 had an O restoration which was failing. Leaking margins and secondary caries were present.

Radiographs, pulp sensibility testing with cold, percussion and periodontal assessment were carried out. Tooth 16 was asymptomatic to all tests. Periodontal health was deemed to be good. Some gingival inflammation was noted.

A significant amount of time was spent discussing general prevention strategies through daily F- exposure, brushing and flossing techniques. Further reinforced with education on caries and periodontal disease processes.

16 was treatment planned and patient was offered a number of restorative options. Patient chose to manage caries and have restoration.

Treatment was carried out on 16 under local anaesthetic and rubber dam. Upon removal of previously placed restoration extent of secondary caries was greater than radiologically expected. Explained to patient the impact of this on the tooth and prognosis for restoration. Patient was provided with option of removal of as much caries as possible and application of 2 step Riva Star

Aqua (Silver Fluoride and Potassium Iodide) as part of caries management process.

Patient was provided risks and benefits. Main benefit was the potential for arresting caries into the long term and preserving maximum tooth structure. 16 is a critical tooth for masticatory function, occlusal scheme management and aesthetics. Main risks discussed were staining of tooth and possibility caries may not be completely arrested. Also explained initially there would be a temporary restoration overlying treated caries.

Patient was provided a consent form, discussed in detail any concerns. Patient noted the ability to arrest caries and preserve tooth was the main driver for going through with treatment. In particular due to his young age and needing to utilise this tooth for a significant amount of time into the future.

16 O caries was removed using slow speed round burs, Riva Star Aqua Step 1 (Silver Fluoride) applied (Image 2), rubbed in with microbrush until colour change noted on all aspects of dentine (Image 3), Riva Star Step 2 (Potassium Iodide) followed and application time was based on the white exudate (Image 4) being removed completely (approximately 2 minutes), colour change was noted after this step, RIVA SC was placed overlying the cavity. (Image 5)

The patient returned after 3 weeks for a review. Pulp sensibility testing with cold, percussion and periodontal assessment were carried out. Tooth 16 was asymptomatic to all tests. Periodontal health was deemed to be good. Gingival inflammation had also resolved. There was very minimal staining

noted. The patient provided he was not aware of any major aesthetic changes to the tooth. Patient was quite happy with the treatment and pleased the tooth did not have any caries remaining in tooth.

Patient was again reviewed a week later and nil symptoms were noted. Patient remained pleased with treatment.

CONCLUSION:

Riva Star Aqua is an excellent addition to a clinician's tool kit. With appropriate training and education, it is a valuable option for preserving teeth in dentition for longer. Patients appreciate the opportunity of saving their teeth and having more time to work towards improved oral health.



About the Author

Dr. Vaibhav Garg has worked in Oral Health and Public Health since 2012. His career working in low-resource, rural, and remote settings in Australia and overseas has provided insight into the value of innovations for communities most in need of dental care.




Image 1 Tooth prior application of Riva Star Aqua Step 1



Image 2 Application of Riva Star Aqua Step 1 on page 2




Image 3 Post application of Riva Star Aqua Step 1 - Colour change is obvious with caries present.




Image 4 Riva Star Step 2 applied. White precipitate has turned clear.




Image 5 Application of Riva GIC



Image 5 Application of Riva GIC



Image 6 Completed restoration



Image 6 Completed restoration

THE NEXT GENERATION SDF SYSTEM



SDF INNOVATION: RIVA STAR

Drawing on 50 years of dental research experience and proudly developed in SDI's state of the art manufacturing facility in Australia, Riva Star is a silver diamine fluoride solution system used to desensitise tooth pain immediately.

This next generation SDF treatment delivers Minimal Intervention Dentistry that adopts the philosophy of integrating prevention, remineralisation and minimal intervention.

RIVA STAR BOTTLE KIT	MRP \$144.36
RIVA STAR CAPSULE KIT	MRP \$131.23

NEXT GENERATION SDF: RIVA STAR AQUA

SDI continues to innovate SDF solutions by releasing Riva Star Aqua, a next generation water based silver fluoride solution that is also effective in arresting carious lesions as an SDF solution¹, with additional patient benefits.

Similar to Riva Star (SDF), Riva Star Aqua [AgF] is also a non-invasive 2 step patented system used to desensitise tooth pain in a few drops. Riva Star Aqua is an aqueous silver fluoride solution without the Ammonia base.

RIVA STAR AQUA BOTTLE KIT	MRP \$158.79
RIVA STAR AQUA CAPSULE KIT	MRP \$144.36

\$115^{*}

BOTTLE KITS

\$99^{*}

CAPSULE KITS

ADVANTAGES	RIVA STAR	RIVA STAR AQUA
Two years and beyond desensitising effect	✓	✓
Staining reduction / elimination	✓	✓
Deep penetration and prevention of biofilm	✓	✓
Antimicrobial activity	✓	✓
Increasing bond strengths of GICs	✓	✓
Storage requirements	Bottles: Room temperature Capsules: Refrigerated	✓ Room temperature
Soft tissue irritation	Mild discomfort	✓ No discomfort
Odour / taste	Slight ammonia taste and odour	✓ No ammonia taste or odour
Non-invasive / MID	✓	✓

¹ Turton B, Horn R, Durward C. Caries arrest and lesion appearance using two different silver fluoride therapies on primary teeth with and without potassium iodide : 12 month results. Clin Exp Dent res 2020;1-11

EQUIPMENT

RADII XPERT \$1500*

SMARTER DESIGN.
OPTIMAL CURE EVERY TIME.

Radii Xpert's unique target assist technology enables you to correctly position the light and be confident of the angulation before curing. The collimated beam ensures light intensity is maintained on the restoration being cured. Achieve a predictable cure, every time with Radii Xpert.

360°
ROTATABLE HEAD

185G
LIGHT WEIGHT

2.5 HOUR
RUNTIME



PURCHASE A
RADII XPERT AND
RECEIVE A **FREE**
STELA INTRO KIT



STELA AUTOMIX INTRO KIT MRP \$243.24

SAVE
\$1104.58*

RADII XPERT
MRP \$2361.34



ACCESSORIES

POINTS

DISPOSABLE BRUSH APPLICATORS

- Points disposable brush applicators are used to quickly transfer small amounts of material into areas that are difficult to reach. Points are resistant to common solvents in dental materials.
- Solvent resistant, non-absorbent, non-linting fibers
- Boxes of 400 in medium, fine or superfine

BOX OF 400 MRP \$51.91
POINTS DISPENSER MRP \$43.07



SDI

YOUR
SMILE.
OUR
VISION.

Call 1800 337 003
info@sdi.com.au
www.sdi.com.au
www.polawhite.com.au

New South Wales

Elli Mouhtouris 0400 849 017
Aurelie Lequimener 0418 210 341
Susanna Kitchen 0400 305 512

Western Australia

Toni Doherty 0448 086 925

South Australia/Northern Territory

Jade Knowles 0417 080 771

Victoria

Jo Coyle 0419 313 765
Shelby MacDonald 0408 174 626

Queensland

Natalie Johnson 0409 807 290
Dora Zanin 0428 910 805

*Offers are valid 1st January - 28th February 2026. Offers cannot be combined. All prices include GST. Free goods will be of equal or lesser value. Purchases must be on one order to qualify for free goods. The above-mentioned offers are not applicable on corporate, government accounts or those customers who have received discounted pricing. Offers are valid in Australia only.