



HIGHLIGHTS

September | October 2025



NEW



RIVA APPLICATOR PRO

EASY. STRONG. CLEAN.

The Riva Applicator Pro represents SDI's third generation of applicators, setting a new standard in compatibility, reliability, and hygiene. Designed for seamless use and high performance, it makes extruding capsules simple and easy.

The Riva Applicator Pro is a durable, universal capsule applicator designed for dentists. Compatible with SDI and other glass ionomer capsules, its advanced lock system provides reliable extrusion. Featuring an ergonomic, maintenance-free design, the applicator withstands drops up to 1.5 meters and requires no lubrication. Fully autoclavable and easy to clean, maintaining a professional, hygienic appearance. An essential tool for modern dental practices.

Compatible with SDI and other glass ionomer capsules

Contoured design for strength and comfort

Engineered stainless steel mechanical parts and durable mechanical casing

Autoclavable and ultrasonic bath tolerant

RIVA APPLICATOR PRO

MRP \$353.01

Take-home whitening

4+1^x OR 10+3^x

BUY 4 GET 1 FREE

BUY 10 GET 3 FREE

POLA DAY

6%, 7.5%, and 9.5% HYDROGEN PEROXIDE

- Fast acting : 30 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximise patient comfort



POLA NIGHT

10%, 16% 18% and 22% CARBAMIDE PEROXIDE

- From 45 minutes to over night wear
- High viscosity, neutral pH gel contains Potassium Nitrate to maximise patient comfort



POLA DAY CP

35% CARBAMIDE PEROXIDE

- Super fast : Just 15 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximise patient comfort



POLA FOR ALIGNERS

STRAIGHTEN AND BRIGHTEN TOGETHER

- Available in **10% CP**
- Contains Potassium Nitrate and fluoride to reduce sensitivity for long term use



POLA 4 X 1.3g MINI KIT	MRP \$41.47
POLA 10 X 1.3g SYRINGE KIT	MRP \$101.57
POLA 50 X 3g DISPENSER PACK	MRP \$564.30
POLA FOR ALIGNERS	MRP \$87.10



P
O
L
A

POLA LIGHT

6% and 9.5% HYDROGEN PEROXIDE

- No custom trays required
- Reduced patient consult time
- Entry level whitening option
- Great front desk conversation starter



DENTIST PRICE MRP \$109.15
RECOMMENDED PATIENT PRICE \$189.00

POLA LUMINATE

6% HYDROGEN PEROXIDE

- No trays required - quick and easy professional brush on tooth whitener
- Apply twice a day for just 30 minutes



DENTIST PRICE MRP \$30.03
RECOMMENDED PATIENT PRICE \$75.00

In-office whitening

POLA RAPID

38% HYDROGEN PEROXIDE

- Fast and simple to use:
- From 3 x 8 minute applications
- Non-stick gel for faster application and removal
- Enhanced blue gel for superior application visibility



POLA RAPID 3 PATIENT KIT	MRP \$464.48
POLA RAPID 3 PATIENT KIT W/O	MRP \$494.73
POLA RAPID 1 PATIENT KIT	MRP \$164.63
POLA RAPID 1 PATIENT KIT W/O	MRP \$176.14

POLA RAPID
1 AND 3 PATIENT KITS

2+1^x
BUY 2 GET 1 FREE



SDI | STELA



THE FUTURE OF COMPOSITES

- ✓ Unlimited depth of cure (self cure)
- ✓ Gap-free interface
- ✓ Only 2 steps: 15 sec prep
- ✓ High strength
- ✓ Choose: capsule or automix



Cavity prepared and ready to be isolated with rubber dam



Final aspect after finishing and polishing

Images courtesy of Prof Dr Rocio Lazo - Peru

25%^{*}
DISCOUNT



The application of Stela is very easy and after 4 minutes the restoration had the same shade as the tooth. Great polishability!

PROF DR ROCIO LAZO - [PERU]



STELA AUTOMIX 1x 8g SYR REFILL	MRP \$191.26
STELA AUTOMIX INTRO KIT	MRP \$231.66
STELA CAPSULE REFILL 50 CAPS	MRP \$278.08
STELA CAPSULES INTRO KIT	MRP \$138.61

GIC & Composite promotion

only available if ordered direct from SDI


COMPOSITES

LUNA 2 | AURA EASY | AURA BULK FILL

MUST BE PURCHASED IN MULTIPLES OF 10

\$70^{*}
PER SYRINGE

\$80^{*}
PER COMPLETE REFILLS



LUNA 2

UNIVERSAL COMPOSITE

- Improved non-stick, non-slump handling and reliable mechanical properties
- BPA free
- Very high radiopacity of 250% Al

LUNA 2 4g SYR	MRP \$115.56
LUNA 2 20 X 0.25g Complets	MRP \$144.86
LUNA 2 INTRO KIT	MRP \$542.70



BEFORE



AFTER



Images courtesy of Dr Fady Toussoun (Australia)



AURA EASY

EASIER AND FASTER AESTHETICS FOR ANTERIOR AND POSTERIOR

- 4 shades cover 16 VITA classical A1-D4® shades
- Increased wear resistance and strength
- Smooth high gloss finish every time
- Non-sticky consistency makes it simple to sculpt
- Low volumetric shrinkage to offset marginal stress

AURA EASY 4g Syr	MRP \$116.75
AURA EASY 20 X 0.25g Complets	MRP \$131.33




BEFORE



AFTER



Images courtesy of Dr Terry Wong (Australia)



AURA BULK FILL

SPEEDY AND SIMPLE POSTERIOR COMPOSITE RESTORATIONS

- Natural-looking finish with just one universal shade
- High depth of cure without the need for a capping layer
- Easy placement in a matter of seconds

AURA BULK FILL	
4g SYR REFILLS	MRP \$117.42
20 X 0.25g Complets	MRP \$132.09



BEFORE



AFTER



Images courtesy of Dr Bob Lowe (USA)



GLASS IONOMERS



RIVA SELF CURE/HV

BULK FILL GLASS HYBRID RESTORATIVE SYSTEM

- Contains SDI's proprietary **ionglass™**, a high ion releasing, radiopaque, bioactive glass that remineralises and prevents secondary caries formation
- Very high fluoride releasing and recharging ability

RIVA SELF CURE MRP \$270.15



— RIVA LIGHT CURE HV —



RIVA LIGHT CURE/HV

LIGHT CURED RESIN REINFORCED GLASS IONOMER RESTORATIVE MATERIAL

- Contains SDI's proprietary **ionglass™**, a high ion releasing, radiopaque, bioactive glass that remineralises and fights against secondary caries formation

RIVA LIGHT CURE MRP \$301.70

BUY 2 BOXES OF 50 CAPSULES [Riva/Stela]
GET A RIVA APPLICATOR PRO

FREE*



SAVE
\$353.01*

RIVA LUTING

- Chemically adheres to tooth structure and metal substrates
- Very low film thickness - flows easily between restoration and tooth
- Exceptional bond strengths
- BPA & HEMA free

RIVA LUTING MRP \$289.39



RIVA PROTECT

- ACP nanotechnology - supercharges remineralisation
- Very high fluoride releasing and recharging ability to protect teeth from acid attacks

RIVA PROTECT MRP \$217.45



RIVA LUTING PLUS

- Chemically adheres to tooth and substrates
- Exceptional bond strengths
- High flexural strength - reduced cracking and failure in flexion

RIVA LUTING PLUS MRP \$308.78



RIVA SILVER

- Silver reinforced
- Ideal for abfraction lesion restorations
- Ideal for bulk filling
- No shrinkage

RIVA SILVER MRP \$478.70



THE NEXT GENERATION SDF SYSTEM



SDF INNOVATION: RIVA STAR

Drawing on 50 years of dental research experience and proudly developed in SDI's state of the art manufacturing facility in Australia, Riva Star is a silver diamine fluoride solution system used to desensitise tooth pain immediately.

This next generation SDF treatment delivers Minimal Intervention Dentistry that adopts the philosophy of integrating prevention, remineralisation and minimal intervention.

RIVA STAR BOTTLE KIT	MRP \$144.36
RIVA STAR CAPSULE KIT	MRP \$131.23

NEXT GENERATION SDF: RIVA STAR AQUA

SDI continues to innovate SDF solutions by releasing Riva Star Aqua, a next generation water based silver fluoride solution that is also effective in arresting carious lesions as an SDF solution¹, with additional patient benefits.

Similar to Riva Star [SDF], Riva Star Aqua [AgF] is also a non-invasive 2 step patented system used to desensitise tooth pain in a few drops. Riva Star Aqua is an aqueous silver fluoride solution without the Ammonia base.

RIVA STAR AQUA BOTTLE KIT	MRP \$158.79
RIVA STAR AQUA CAPSULE KIT	MRP \$144.36



\$99^x

BOTTLE KITS

\$115^x

CAPSULE KITS



¹ Turton B, Horn R, Durward C. Caries arrest and lesion appearance using two different silver fluoride therapies on primary teeth with and without potassium iodide : 12 month results. Clin Exp Dent res 2020;1-11

ADHESIVES

THE ALL-IN-ONE ADHESIVE SOLUTION FOR EVERY SITUATION

UNIVERSAL ADHESIVE

- Consistent bond strengths to dentine and enamel
- Wet and dry dentine bonding
- Proprietary combination of adhesive monomers including MDP
- Predictable bond strengths to indirect substrates
- Effective hybridisation

ZIPBOND REFILL 5mL BOTTLE	MRP \$155.19
---------------------------	--------------



DISCOUNTED PRICE
\$99.00^x

SAVE
\$56.19^{*}



CEMENTS

RIVA CEM AUTOMIX

RESIN MODIFIED GLASS IONOMER LUTING CEMENT

SUPERIOR BOND, HIGH STRENGTH AND LESS WASTE

Riva Cem Automix is a self-curing, radiopaque, fluoride releasing paste/paste resin modified glass-ionomer luting cement. It is indicated for the permanent cementation of metal and ceramic restorations. The automix tip attachment offers convenience over traditional hand mixed cement systems.

SPECIAL PRICE OF
\$99.00*

Limit 2 per clinic




Less
Cement
Waste

SAVE
\$21.07*



» The dispensing tips are great and allow me to use less cement, which reduces waste. «

 Dr Byron Davis, DDS, MS Okatie, SC



ACCESSORIES

POINTS

DISPOSABLE BRUSH APPLICATORS

- Points disposable brush applicators are used to quickly transfer small amounts of material into areas that are difficult to reach. Points are resistant to common solvents in dental materials.
- Solvent resistant, non-absorbent, non-linting fibers
- Boxes of 400 in medium, fine or superfine

BOX OF 400	MRP \$47.19
POINTS DISPENSER	MRP \$42.23





Dr. Fady Tossoun

SMILE TRANSFORMATIONS WITH SIMPLE DENTISTRY

At this seminar you will learn:

- Complex treatment planning
- Importance of communication with patients. Handling objections and complications
- Clinical tips and tricks on how to use composite resin to enhance the smiles of your every day patients
- Etch and bonding techniques
- Importance of understanding tooth form and anatomy
- How occlusion effects our results
- The importance of material choice



To register
Scan QR code

Melbourne

Wednesday 3rd September 2025

Cost : \$99



Dr. Vaibhav Garg

THE USE OF MINIMALLY INVASIVE TECHNIQUES IN RESTORATIVE DENTISTRY

For dentists, hygienists & therapists:

- A bespoke facilitated and interactive workshop for general dentists on minimally invasive dentistry and clinical procedures with Q&A.
- Restorative procedures carried out in remote communities using composites and glass ionomer materials.
- The benefits of using Silver Fluoride (AGF) and Silver Diamine Fluoride (SDF) plus the differences in these two minimally invasive restorative tools.



To register
Scan QR code



To register
Scan QR code

Mecure Hotel Perth

Wed 10th Sept 2025

Cost : \$99

Rydges Hobart

Tuesday 25th Nov 2025

Cost : \$99



Dr. James Bartalotta

MASTERING ANTERIOR AESTHETIC COMPOSITES: TECHNIQUES AND APPLICATIONS

At this seminar you will learn:

- Refined Composite Placement
- Material Mastery
- Hands-On Practice
- Complex Case Management



To register
Scan QR code

Amora Hotel Richmond

Wed 3rd December 2025

Cost : \$99



Dr Isabel Holmes

WHITEN AND RESTORE LIKE A PRO

For dentists, hygienists & therapists:

- Chairside and take-home whitening options overview
- Useful tips on how to customise whitening for your patients
- A step-by-step guide on in-chair whitening
- Composite Restorations
- Glass Ionomer Cement Restorations
- Cementing Crowns (simplified)
- Treatment Planning TIPS



To register
Scan QR code

Rydges Esplanade Resort Cairns

Thursday 6th December

Cost : \$99

SDI

YOUR
SMILE.
OUR
VISION.

Call 1800 337 003

info@sdi.com.au

www.sdi.com.au

www.polawhite.com.au

New South Wales

Elli Mouhtouris 0400 849 017

Aurelie Lequimener 0418 210 341

Western Australia

Toni Doherty 0448 086 925

South Australia/Northern Territory

Jade Knowles 0417 080 771

Victoria

Jo Coyle 0419 313 765

Snezana Josifoski 0418 567 758

Shelby MacDonald 0408 174 626

Tasmania

Snezana Josifoski 0418 567 758

Queensland

Natalie Johnson 0409 807 290

Dora Zanin 0428 910 805

*Offers are valid 1st September - 31st October 2025. Offers cannot be combined. All prices include GST. Free goods will be of equal or lesser value. Purchases must be on one order to qualify for free goods. The above-mentioned offers are not applicable on corporate, government accounts or those customers who have received discounted pricing. Offers are valid in Australia only.