



HIGHLIGHTS

March | April 2025



TOOTH WHITENING

TAKE-HOME WHITENING

4+1*
BUY 4, GET 1 FREE

10+3*
BUY 10, GET 3 FREE



POLA DAY

HYDROGEN PEROXIDE BASED ADVANCED TRAY TOOTH WHITENING SYSTEM

- Available in **6% 7.5%** and **9.5% hydrogen peroxide** applications
- From 30 minutes once a day

4 X 1.3g MINI KIT	MRP \$41.47
10 X 1.3g SYRINGE KIT	MRP \$101.57
50 X 3g DISPENSER PACK	MRP \$564.30



POLA NIGHT

CARBAMIDE PEROXIDE BASED ADVANCED TRAY TOOTH WHITENING SYSTEM

- Available in **10% 16% 18%** and **22% carbamide peroxide**
- From 45 minutes once a day

4 X 1.3g MINI KIT	MRP \$41.47
10 X 1.3g SYRINGE KIT	MRP \$101.57
50 X 3g DISPENSER PACK	MRP \$564.30



POLA DAY CP

SUPER FAST TRAY TOOTH WHITENING SYSTEM!

- Available in **35% carbamide peroxide**
- From just 15 minutes per day

4 X 1.3g MINI KIT	MRP \$41.47
10 X 1.3g SYRINGE KIT	MRP \$101.57
50 X 3g DISPENSER PACK	MRP \$564.30



POLA LUMINATE

BRUSH ON & GO!

- Available in **6% hydrogen peroxide**, low viscosity fast drying gel
- No trays required - quick and easy professional brush on tooth whitener
- Apply twice a day for just 30 minutes

DENTIST PRICE \$27.30
RECOMMENDED PATIENT PRICE \$75.00



POLA LIGHT

COMPLETE YOUR WHITENING OFFERING WITH POLA LIGHT

- Available in **6% and 9.5% hydrogen peroxide**
- No custom trays required
- Reduced patient consult time
- Entry level whitening option
- Great front desk conversation starter

DENTIST PRICE \$109.15
RECOMMENDED PATIENT PRICE \$189.00



POLA FOR ALIGNERS

STRAIGHTEN AND BRIGHTEN TOGETHER

- Available in **10% carbamide peroxide**
- The greatest patient comfort in a whitening kit
- Contains potassium nitrate and fluoride to reduce sensitivity for long term use
- High viscosity gel ensures it can be easily and securely placed into the aligner tray, staying in place

POLA FOR ALIGNERS MRP \$87.10



Dr. David Carralero Plaza

Pola Day 6%



TOOTH WHITENING

POLA RAPID SUPER. FAST. WHITENING.



Pola Rapid is the new in office whitening treatment that safely whitens teeth in a 24 minute application time.

SUPERIOR HANDLING AND FAST REMOVAL

Pola Rapid's higher viscosity gel enables faster application and removal due to its specially formulated non-stick attributes. Pola Rapid can be extruded quickly and easily and remains on teeth when applied.

ENHANCED BLUE GEL FOR FASTER APPLICATION

Pola Rapid uses an advanced blue gel for ease of visibility and application. The easy to apply non-stick gel enables faster application and removal.

MINIMAL CHAIR TIME WITH IMMEDIATE RESULTS

A unique **38% hydrogen peroxide** formulation that releases active ingredients immediately. In only one visit, Pola Rapid initiates a powerful whitening treatment designed to whiten teeth in 24 minutes. No pre-desensitiser step is required.

POLA RAPID 3 PATIENT KIT	MRP \$464.48
POLA RAPID 3 PATIENT KIT WITH OPTRAGATE	MRP \$494.73
POLA RAPID 1 PATIENT KIT	MRP \$164.63
POLA RAPID 1 PATIENT KIT WITH OPTRAGATE	MRP \$176.14



Dr Bob Gibbins - Australia



BEFORE



AFTER

3+1^x
BUY 3, GET 1 FREE

EQUIPMENT

IN-CHAIR WHITENING BUNDLE!

Purchase the Radii Xpert for the special price of

\$2100^x

AND RECEIVE THE FOLLOWING FREE*



x1

BLEACH ARCH ATTACHMENT



x1

POLA RAPID 3 PATIENT KIT WITH OPTRAGATE



x1

POLA STAND 2



SAVE
\$2190.35*



RADI XPERT MRP \$2248.90

CEMENTS

RIVA CEM AUTOMIX


RESIN MODIFIED GLASS IONOMER LUTING CEMENT

SAVE
\$21.06*

SUPERIOR BOND, HIGH STRENGTH & LESS WASTE

Riva Cem Automix is a self-curing, radiopaque, fluoride releasing paste/paste resin modified glass-ionomer luting cement. It is indicated for the permanent cementation of metal and ceramic restorations. The automix tip attachment offers convenience over traditional hand mixed cement systems.

SPECIAL PRICE OF
\$99.00*



» The dispensing tips are great and allow me to use less cement, which reduces waste. «

🇺🇸 Dr Byron Davis, DDS, MS Okatie, SC



RIVA CEM AUTOMIX
SYRINGE REFILL 2x4g MRP \$120.06
MIXING TIPS 50PK MRP \$75.23

*Limit of 2 per clinic

ADHESIVES

ZIPBOND

THE ALL-IN-ONE ADHESIVE SOLUTION FOR EVERY SITUATION

- Consistent bond strengths to dentine and enamel
- Wet and dry dentine bonding
- Proprietary combination of adhesive monomers including MDP
- Predictable bond strengths to indirect substrates
- Effective hybridisation
- Exceptional control
- Highly aesthetic
- Pleasant aroma
- Fluoride releasing
- Easy to follow application
- Convenient delivery

ZIPBOND REFILL 5mL BOTTLE MRP \$155.19

DISCOUNTED PRICE
\$99.00



29 CLINICAL EVALUATORS
862 TOTAL USES
93% CLINICAL RATING

- Dental Advisor rates Zipbond ++++ 1/2 or **93% clinical rating**
- 862 uses and **zero post op sensitivity**
- 97% of Dental Advisor evaluators rated it either equal or better than their existing adhesive
- **100% would recommend Zipbond** to a colleague

DEFENCE
AGAINST
CLINICAL
VARIABILITY


SAVE
\$56.19*

ACCESSORIES

POINTS

DISPOSABLE BRUSH APPLICATORS

- Points disposable brush applicators are used to quickly transfer small amounts of material into areas that are difficult to reach. Points are resistant to common solvents in dental materials.
- Solvent resistant, non-absorbent, non-linting fibers
- Boxes of 400 in medium, fine or superfine



BOX OF 400 POINTS MRP \$47.19

BUY 4, GET 3 FREE

ETCHANT

SUPER ETCH

37% PHOSPHORIC ACID ETCHANT GEL FOR ENAMEL AND DENTINE

- Stays where placed.
- Does not run onto soft tissue
- Washes away easily with water
- Easily controlled gel which flows under direction, not gravity
- Blue tint for easy placement visibility
- Thixotropic gel to improve handling



SUPER ETCH BULK KIT MRP \$102.64
SUPER ETCH 12GM SYRINGE MRP \$66.21
SUPER ETCH GEL 15mL BOTTLE MRP \$44.32
SUPER ETCH JUMBO MRP \$220.57

BUY 3, GET 1 FREE

SDI | STELA

THE FUTURE OF COMPOSITES

- ✓ Unlimited depth of cure (self cure)
- ✓ Gap-free interface
- ✓ Only 2 steps: 15 sec prep
- ✓ High strength
- ✓ Choose: capsule or automix



BEFORE
Cavity prepared and ready to be isolated with rubber dam



AFTER
Final aspect after finishing and polishing

Images courtesy of Prof Dr Rocio Lazo - Peru

2+1^x
BUY 2, GET 1 FREE



STELA AUTOMIX 1x 8g SYR REFILL	MRP \$191.26
STELA AUTOMIX INTRO KIT	MRP \$231.66
STELA CAPSULE REFILL 50 CAPS	MRP \$278.08
STELA CAPSULES INTRO KIT	MRP \$138.61

SDI | ULTRAMAT X

THE EASIER, QUIETER CAPSULE MIXER



THE QUIETEST MIXER. USE WITH ONE HAND

The Ultramat X is 30% quieter than the nearest leading capsule mixer*. The Ultramat X's Easy arm is designed with a state-of-the-art clamp design, engineered with safety and practicality. Insert and remove capsules with one hand eliminating glove and finger pinch.



USE WITH ONE HAND



Single hand use

30% quieter than other mixers**

Full capsule compatibility

5 Year warranty

ULTRAMAT X MRP \$1724.80

*Source: Waveform Acoustics assessment (2024). Measured 300mm from length / height of each unit.



BUY ULTRAMAT X RECEIVE ANY BOX OF GIC 50 CAPS **FREE**

ULTRAMAT 2 \$999*

SAVE
\$536.91*

HIGH SPEED MULTI-USE TRITURATOR

- Low vibration and extremely quiet
- Accurate and consistent trituration
- For mixing most encapsulated dental materials
- Safe and simple to use



5 YEAR LIMITED WARRANTY

ULTRAMAT 2 MRP \$1535.91

PLASTVAC \$999*

SAVE
\$538.71*

TRAYFORM VACUUM MACHINE

- Powerful motor
- Quick and practical
- Suitable for round and square thermoplastic materials of different dimensions



PLASTVAC MRP \$1537.71



YOUR SMILE.
OUR VISION.

Call 1800 337 003
info@sdi.com.au
www.sdi.com.au
www.polawhite.com.au



New South Wales
Elli Mouhtouris 0400 849 017
Allison Ohlback 0400 305 512
Western Australia
Toni Doherty 0448 086 925
South Australia/Northern Territory
Andriana Spithas 0417 080 771

Victoria
Jo Coyle 0419 313 765
Snezana Josifoski 0418 567 758
Shelby MacDonald 0408 174 626
Tasmania
Snezana Josifoski 0418 567 758
Queensland
Natalie Johnson 0409 807 290
Dora Zanin 0428 910 805

*Offers are valid 1st March - 30th April 2025. Offers cannot be combined. All prices include GST. Free goods will be of equal or lesser value. The above-mentioned offers are not applicable on corporate, government accounts or those customers who have received discounted pricing. Offers are valid in Australia only.