



# FEATURES

Jan | Feb | Mar 2026



SDI (North America) Inc.  
[www.sdi.com.au/en-ca](http://www.sdi.com.au/en-ca)

# STELA

## THE FUTURE OF COMPOSITES

Stela is a self cure flowable composite offering a gap-free interface, high strength, unlimited depth of cure and a simplified application process for easy to place restorations.

### COMPETITOR COMPARISON

	STELA	FILTEK ONE BULK FILL	ACTIVA BIOACTIVE	BULK EZ
Compressive Strength (MPa)	331*	300	239	271
Flexural Strength (MPa)	143*	157	88	121
Material Type	Flowable	Packable	Flowable	Flowable
Depth of Cure	Unlimited	5mm	4mm	Unlimited
Preparation	Simplified	Standard	Standard	Standard

Source: Manufacturer's brochures and SDI Research & Development Department  
\*Stela Automix



### UNLIMITED DEPTH OF CURE

Self cures fully at any depth, eliminating the need for incremental layering and reducing polymerization uncertainty.



### GAP-FREE INTERFACE

Stela cures towards the tooth, reducing voids and enhancing marginal adaptation.



### HIGH STRENGTH

Engineered for long-term durability with a combination of high compressive and flexural strengths.



### SIMPLE 2 STEP PROCESS

Reduce in-chair time by eliminating etching and light curing from placement protocols.



### BIOACTIVE PROPERTIES

Unlike many light cured composites, Stela releases fluoride and supports mineral retention for long-term tooth health.

### CAPSULE OR SYRINGE YOUR CHOICE



# POLA

## TAKE-HOME WHITENING

**BUY 5, GET 2 FREE\***

### POLA DAY

9.5% HYDROGEN PEROXIDE

- Fast acting: 30 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



### POLA NIGHT

10%, 16% and 22% CARBAMIDE PEROXIDE

- From 45 minutes to overnight wear
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



### POLA DAY CP

35% CARBAMIDE PEROXIDE

- Super fast: Just 15 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



### POLA LIGHT

9.5% HYDROGEN PEROXIDE and  
22% CARBAMIDE PEROXIDE

- From just 20 minutes a day
- Portable, rechargeable and easy to wear



**BUY 2, GET 1 FREE\***

## IN-OFFICE WHITENING

### POLA OFFICE+

Pola Office+, with its 37.5% Hydrogen Peroxide formulation, is designed to whiten teeth in 24 minutes

- Minimal chair time - Improved productivity and profitability
- Contains Potassium Nitrate to maximize patient comfort



SDI

## RIVA LIGHT CURE & RIVA SELF CURE



**BUY 2, GET 1 FREE\***



### RIVA LIGHT CURE/HV

LIGHT CURED RESIN REINFORCED GLASS IONOMER RESTORATIVE MATERIAL

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralizes and fights against secondary caries formation
- Exceptional handling; easy to sculpt and non-sticky
- **Choice of viscosity - High Viscosity formulation is flowable enough to inject and highly sculptable**
- **Unique, patented, orange light safe tip, prevents curing in the tip and significantly extends working time**
- Easy Activation Capsule - 70% easier to activate
- BPA free



### RIVA SELF CURE/HV

BULK FILL GLASS HYBRID RESTORATIVE SYSTEM

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralizes and prevents secondary caries formation
- Exceptional handling; non-sticky, non-slumping and sculptable formulation
- **Choice of viscosity - High Viscosity formulation is flowable enough to inject and highly sculptable**
- Easy Activation Capsule - 70% easier to activate
- BPA & HEMA free

## ALL COMPOSITES ~ 3+1<sup>x</sup> OR 14+7<sup>x</sup>

\*Excludes Wave and Stela range

BUY 3, GET 1 FREE

BUY 14, GET 7 FREE

### UNIVERSAL COMPOSITES



#### LUNA 2

- Non-stick handling that maintains its sculpted shape
- Chameleon effect for accurate shade matching
- High strength & esthetics
- BPA free for health-conscious patients



#### AURA EASY

- 4 shades cover 16 Vita® shades
- Non-stick and easy to sculpt
- Wear resistant with high strength
- Easy handling for effortless placement

### FLOWABLE COMPOSITES



#### LUNA FLOW

- Effortlessly flows to margins with thixotropic properties
- High strength, radiopacity and wear resistance
- Accurate shade matching



#### AURA EASYFLOW

- 4 shades cover 16 Vita® shades
- Thixotropic for maximum control during application
- High flexural strength and radiopacity

### BULK FILL COMPOSITES



#### AURA BULK FILL

- 5mm depth of cure
- Non-stick, waxy consistency
- Low volumetric shrinkage

PLEASE SCAN FOR MORE DETAILED INFORMATION ON THE FULL RANGE OF SDI COMPOSITES.



ALL PROMOTIONS AVAILABLE THROUGH THE FOLLOWING SDI DISTRIBUTION PARTNERS



**SDI (North America) Inc.**  
1279 Hamilton Parkway  
Itasca, IL 60143

**Toll Free 1 800 228 5166**  
usa.canada@sdi.com.au  
www.sdi.com.au

**Redeem your free goods by sending your invoice to:**  
usa.canada@sdi.com.au

\*Please allow 2-4 weeks for delivery. Please note offers cannot be combined. Offers cannot be used in conjunction with contract pricing. Valid on invoices dated January 1st - March 30th, 2026. All claims for free goods must be made within 30 days from the end of promotion date. **Purchases must be on one order to qualify for free goods.** Free products must be of equal or lesser value. Offers are valid in Canada only.