



HIGHLIGHTS

Oct | Nov | Dec 2025

NEW ZEALAND



WHITENING

Call Courtney Gerlach: 021 400 323

POLA

ADVANCED TOOTH WHITENING SYSTEM

Take-home whitening



POLA DAY

HYDROGEN PEROXIDE BASED
ADVANCED TRAY TOOTH
WHITENING SYSTEM

- Available in **7.5%** and **9.5% hydrogen peroxide**
- From 30 minutes once a day

4 X 1.3g MINI KIT
10 X 1.3g SYRINGE KIT
50 X 3g DISPENSER PACK



POLA NIGHT

CARBAMIDE PEROXIDE BASED
ADVANCED TRAY TOOTH
WHITENING SYSTEM

- Available in **10%, 16% and 22% carbamide peroxide**
- From 45 minutes once a day

4 X 1.3g MINI KIT
10 X 1.3g SYRINGE KIT
50 X 3g DISPENSER PACK



POLA DAY CP

SUPER FAST TRAY TOOTH
WHITENING SYSTEM!

- Available in **35% carbamide peroxide**
- From just 15 minutes per day tray wear time

4 X 1.3g MINI KIT
10 X 1.3g SYRINGE KIT
50 X 3g DISPENSER PACK



POLA LUMINATE

BRUSH ON & GO!

- Available in **6% hydrogen peroxide** low viscosity fast drying gel
- No trays required - quick and easy professional brush on tooth whitener

POLA LUMINATE

POLA FOR ALIGNERS

STRAIGHTEN AND BRIGHTEN TOGETHER

- Available in **10% carbamide peroxide**
- High viscosity gel ensures it can be easily and securely placed into the aligner tray, staying in place



POLA FOR ALIGNERS 4 X 3g

POLA LIGHT

6% and 9.5% HYDROGEN PEROXIDE

- No custom trays required
- Reduced patient consult time
- Entry level whitening option
- Great front desk conversation starter

POLA LIGHT 3 X 3g SYRINGES



All Whitening **3+1^{*}** BUY 3 GET 1 FREE | **10+5^{*}** BUY 10 GET 5 FREE

In-office whitening

POLA RAPID

38% HYDROGEN PEROXIDE

- Fast and simple to use:
- From 3 x 8 minute applications
- Non-stick gel for faster application and removal
- Enhanced blue gel for superior application visibility



POLA RAPID 3 PATIENT KIT
POLA RAPID 3 PATIENT KIT W/O
POLA RAPID 1 PATIENT KIT
POLA RAPID 1 PATIENT KIT W/O



STELA | LUNA 2 | LUNA FLOW - NOW 30% OFF



STELA
THE FUTURE
OF COMPOSITES

50 years of dental material experience has engineered a two step restorative that is stronger, gap free, easier to place and perfectly cures every time.

- ✓ Unlimited depth of cure (self cure)
- ✓ Gap-free interface
- ✓ Only 2 steps: 15 sec prep
- ✓ High strength
- ✓ Choose: capsule or automix

BPA &
HEMA
FREE

STELA AUTOMIX 1 SYR REFILL
STELA AUTOMIX INTRO KIT
STELA CAPSULE REFILL 50 CAP
STELA CAPSULES INTRO KIT
STELA PRIMER REFILL 5ML BT



BPA
FREE
HEMA
FREE

SUPERIOR STRENGTH
& AESTHETICS

LUNA 2

- UNIVERSAL COMPOSITE
- Improved non-stick, non-slump handling and reliable mechanical properties
 - Logical Shade Matching Technology for easier and accurate shade selection
 - BPA free
 - Very high radiopacity of 250% AI
 - Flexural strength similar to dentine for Class III, IV and V restorations and proximal boxes
 - Opalescent and fluorescent like the tooth for lifelike results

LUNA 2 4g SYRINGE
LUNA 2 20 X 0.25g COMPLETS
LUNA 2 INTRO KIT



Photos courtesy:
Dr Neus Sancho Lecina - SPAIN




HEMA
FREE
Dental Product Shopper
Recommended Product
catapult
Trusted by Dental Experts

PERFECT FLOW
PERFECT STAY


LUNA FLOW

- UNIVERSAL FLOWABLE COMPOSITE
- Exceptionally high compressive strength of 421 MPa and 128 MPa of flexural strength
 - Flows on demand, holds its shape and does not string for maximum control over application
 - Outstanding radiopacity of 265% AI for faster and more precise diagnosis
 - Logical Shade Matching Technology for easier and accurate shade selection
 - Two viscosity options for all clinical indications: Luna Flow offers medium/high viscosity and Luna Flow LV offers low viscosity

LUNA FLOW/LV 2G SYRINGE
LUNA FLOW/LV BULK KIT 5 X 2G SYRINGES
LUNA FLOW/LV INTRO KIT 4 X 2G SYRINGES



Photos courtesy:
Dr Carla Cohn - CANADA



*Flow at 90° angle for 40 seconds



BUT 10 LUNA 2 SYRINGES
OR COMPLETS
GET A RADII XPERT

FREE*

**Offer cannot be combined

Professional LED Curing Light

Smarter design. Optimal cure every time.

Target assist precision curing

Radii Xpert's unique target assist technology enables you to correctly position the light and be confident of the angulation before curing. A light source must be held close and perpendicular to a restoration to avoid incomplete polymerization¹. By activating the positioning light and releasing once over the tooth, you can be sure that curing only begins when the light is in position and restorations will receive a complete cure. Ensure you are on target every time with Radii Xpert.



TARGET ASSIST
PRECISION CURING



COLLIMATED
LIGHT



SMART TECHNOLOGY
FOR EASE OF USE



ERGONOMICALLY
WELL-BALANCED GRIP



VALUE
\$3023.00*



ULTRAMAT 2

20% DISCOUNT

RRP \$1867.00

HIGH SPEED MULTI-USE TRITURATOR

- Low vibration and extremely quiet
- Accurate and consistent trituration
- For mixing most encapsulated dental materials
- Safe and simple to use

HIGH SPEED,
MULTI-USE
TRITURATOR

SAVE
\$374.40*

**5 YEAR
LIMITED
WARRANTY**



HENRY SCHEIN®
DENTAL

Rely on Us™

free phone 0800 808 855
orders@henryschein.co.nz
www.henryschein.co.nz

*All prices exclude GST. Offers cannot be combined. Valid on invoices dated October 1st - December 31st 2025. All claims for free goods must be made within 30 days from the end of promotion date. The above-mentioned offers are not applicable on corporate, government accounts or those customers who have received discounted pricing. Bonus goods must be of equal or lesser value. Offers are valid in New Zealand only.

Courtney Gerlach
Sales Product Specialist -
New Zealand
SDI Limited
m: +64 21 400 323
e: Courtney.Gerlach@sdi.com.au

SDI

YOUR
SMILE.
OUR
VISION.