



FEATURES

Jan | Feb | Mar 2026



SDI (North America) Inc.
www.sdi.com.au/en-us

STELA

THE FUTURE OF COMPOSITES

Stela is a self cure flowable composite offering a gap-free interface, high strength, unlimited depth of cure and a simplified application process for easy to place restorations.

COMPETITOR COMPARISON

	STELA	FILTEK ONE BULK FILL	ACTIVA BIOACTIVE	BULK EZ
Compressive Strength (MPa)	331*	300	239	271
Flexural Strength (MPa)	143*	157	88	121
Material Type	Flowable	Packable	Flowable	Flowable
Depth of Cure	Unlimited	5mm	4mm	Unlimited
Preparation	Simplified	Standard	Standard	Standard

Source: Manufacturer's brochures and SDI Research & Development Department
*Stela Automix



UNLIMITED DEPTH OF CURE

Self cures fully at any depth, eliminating the need for incremental layering and reducing polymerization uncertainty.



GAP-FREE INTERFACE

Stela cures towards the tooth, reducing voids and enhancing marginal adaptation.



HIGH STRENGTH

Engineered for long-term durability with a combination of high compressive and flexural strengths.



SIMPLE 2 STEP PROCESS

Reduce in-chair time by eliminating etching and light curing from placement protocols.



BIOACTIVE PROPERTIES

Unlike many light cured composites, Stela releases fluoride and supports mineral retention for long-term tooth health.

CAPSULE OR SYRINGE YOUR CHOICE



POLA

SDI

TAKE-HOME WHITENING

BUY 5, GET 2 FREE*

POLA DAY

7.5% and 9.5% HYDROGEN PEROXIDE

- Fast acting: 30 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



POLA NIGHT

10%, 16% and 22% CARBAMIDE PEROXIDE

- From 45 minutes to overnight wear
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



POLA DAY CP

35% CARBAMIDE PEROXIDE

- Super fast: Just 15 minutes a day
- High viscosity, neutral pH gel contains Potassium Nitrate to maximize patient comfort



POLA LIGHT

9.5% HYDROGEN PEROXIDE and 22% CARBAMIDE PEROXIDE

- From just 20 minutes a day
- Portable, rechargeable and easy to wear



POLA FOR ALIGNERS

10% and 35% CARBAMIDE PEROXIDE

- Straighten and brighten together
- The greatest patient comfort in a whitening kit



IN-OFFICE WHITENING

BUY 2, GET 1 FREE*

POLA RAPID

38% HYDROGEN PEROXIDE

- Fast and simple to use
- Non-stick gel for faster application and removal
- Enhanced blue gel for superior application visibility



POLA RAPID 1 & 3 PATIENT KITS

ALL COMPOSITES

*Excludes Stela and Wave range

3+1^x OR 14+7^x

BUY 3, GET 1 FREE

BUY 14, GET 7 FREE

UNIVERSAL COMPOSITES



LUNA 2

Luna 2 is SDI's standout universal composite, with a unique blend of fillers to deliver outstanding optical and mechanical properties.

- Non-stick handling that maintains its sculpted shape
- Chameleon effect for accurate shade matching
- High Strength & Esthetics
- BPA free for health-conscious patients



AURA EASY

Aura Easy is a universal composite available in 4 simplified shades to ensure easy restorations with reduced inventory.

- 4 shades cover 16 Vita® shades
- Non-stick and easy to sculpt
- Wear resistant with high strength
- Easy handling for effortless placement

FLOWABLE COMPOSITES



LUNA FLOW

Luna Flow is an extra strong flowable composite available in two viscosities that offer excellent mechanical and optical properties.

- Effortlessly flows to margins with thixotropic properties
- High strength, radiopacity and wear resistance
- Accurate shade matching



AURA EASYFLOW

Aura Easyflow is a universal flowable composite available in 4 simplified shades with optimized features for every layer of use.

- 4 shades cover 16 Vita® shades
- Thixotropic for maximum control during application
- High flexural strength and radiopacity

BULK FILL COMPOSITES



AURA BULK FILL

Aura Bulk Fill is a posterior composite with a natural looking finish available in just one universal shade for easy and quick restorations. The waxy consistency is packable and non-sticky to optimize handling.

- 5mm depth of cure
- Non-stick, waxy consistency
- Low volumetric shrinkage

PLEASE SCAN FOR MORE DETAILED INFORMATION ON THE FULL RANGE OF SDI COMPOSITES.



SDI GLASS IONOMERS

BUY 2, GET 1 FREE*



RIVA SELF CURE/HV

BULK FILL GLASS HYBRID RESTORATIVE SYSTEM

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralizes and prevents secondary caries formation
- Exceptional handling; non-sticky, non-slumping and sculptable formulation
- **Choice of viscosity - High Viscosity formulation is flowable enough to inject and highly sculptable**
- Easy Activation Capsule - 70% easier to activate
- BPA & HEMA free



RIVA LIGHT CURE/HV

LIGHT CURED RESIN REINFORCED GLASS IONOMER RESTORATIVE MATERIAL

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralizes and fights against secondary caries formation
- Exceptional handling; easy to sculpt and non-sticky
- **Choice of viscosity - High Viscosity formulation is flowable enough to inject and highly sculptable**
- **Unique, patented, orange light safe tip, prevents curing in the tip and significantly extends working time**
- Easy Activation Capsule - 70% easier to activate
- BPA free

SDI RIVA STAR AQUA

BUY 1, GET 1 FREE*

NEXT GENERATION SDF: RIVA STAR AQUA

SDI continues to innovate SDF solutions by releasing Riva Star Aqua, a next generation water based silver fluoride solution that is also as effective as an SDF solution, with additional patient benefits.

Similar to Riva Star (SDF), Riva Star Aqua (AgF) is also a non-invasive 2 step patented system used to desensitize tooth pain in a few drops. Riva Star Aqua is an aqueous silver fluoride solution without the Ammonia base.

Riva Star Aqua Bottle Kit only



REDUCE STAINING

Riva Star Aqua Bottle Kit
1 x 1.5mL Riva Star Aqua Step 1 Bottle
1 x 3.0mL Riva Star Step 2 Bottle

8800532



YOUR SMILE. OUR VISION.

SDI (North America) Inc.
1279 Hamilton Parkway
Itasca, IL 60143

Toll Free 1 800 228 5166
usa.canada@sdicom.au
www.sdi.com.au

Redeem your free goods by sending your invoice to:
usa.canada@sdicom.au

facebook.com/sdinainc instagram.com/sdinainc

*Please allow 2-4 weeks for delivery. Please note offers cannot be combined. Offers cannot be used in conjunction with contract pricing. Valid on invoices dated January 1st - March 30th, 2026. All claims for free goods must be made within 30 days from the end of promotion date. Purchases must be on one order to qualify for free goods. Free products must be of equal or lesser value. Offers are valid in USA only.