



**29** CLINICAL EVALUATORS

**862** TOTAL USES

**93%** CLINICAL RATING

**Key features:** Universal bonding agent • Can be used for indirect or direct restorations in any etching technique • No dual-cured activator required


## Description

**ZIPBOND** is a light-cured, single-component universal dental adhesive that is:

- Indicated for direct and indirect restorations.
- Fluoride releasing.
- Compatible with self-etch, total-etch, or selective-etch techniques without the use of a separate dual-cured activator.

## Indications

- Bonding of light, dual, and self-cured composite or compomer restorations.
- Sealing of cavity preparations prior to direct and indirect restoration placement.
- Core build ups.
- Bonding of indirect restorations when combined with a resin-based cement.
- Root surface desensitization.
- Bonding of pit and fissure sealants.
- Protective varnish for glass ionomer restorations.
- Repair of composite and compomer restorations.



## Clinical Tips

- Try it as a desensitizer, as well as for bonding.
- The light color makes it great for anterior esthetics.

**"ZERO POST-OP SENSITIVITY WITH ANY OF THE PATIENTS I USED ZIPBOND ON."**

## Evaluators' Comments

*"Light-colored, thin material. It flowed great and you could tell it got in all the crevices."*

*"Made even coats possible when applying multiple layers."*


*"This was a very easy and versatile bonding agent. I like to keep things simple and efficient and this material allows me to do just that."*

*"ZIPBOND is an excellent desensitizer. It is better than our current products."*


*"A few times, it was very difficult to remove a sectional or circumferential matrix band due to the adhesive bond."*

*"It worked well in all situations in 92 restorations."*

### Evaluation Summary:



### Compared to Competitive Products:



### Consultants who would:

**100%** Recommend to a colleague

### Consultants who would want to stock in their office:

**24%** Yes, instead of current product

**62%** Yes, in addition to current product

**7%** I might want to order this product for certain cases