



YOUR  
SMILE.  
OUR  
VISION.

# LUNA COMPOSITE FAMILY

## BLEACH SHADES EQUIVALENCE TO VITA SYSTEM 3D-MASTER®

### LUNA UNIVERSAL COMPOSITE

Available in two bleach shades:

- B – Bleach shade: equivalent to 0M3 Vita®
- XB - Extra Bleach shade: equivalent to 0M2 Vita®

### LUNA 2 UNIVERSAL COMPOSITE (BPA free)

### LUNA FLOW (Regular Viscosity flowable composite)

### LUNA FLOW LV (Low Viscosity flowable composite)

Available in two bleach shades, being 2XB the lightest of all:

- XB - Extra Bleach shade: equivalent to 0M2 Vita®
- 2XB – 2x Extra Bleach shade: equivalent to 0M1 Vita®

SDI | LUNA  
 SDI | LUNA 2  
 SDI | LUNA FLOW  
 SDI | LUNA FLOW LV

**XB - Extra Bleach**

Vita® 0M2

SDI | LUNA 2  
 SDI | LUNA FLOW  
 SDI | LUNA FLOW LV

**2XB – 2x Extra Bleach**

Vita® 0M1



SDI | LUNA

**B – Bleach**

Vita® 0M3



YOUR  
SMILE.  
OUR  
VISION.

SDI LIMITED ACN 008 075 581 | ABN 27 008 075 581  
 HEAD OFFICE 3-15 BRUNSDON STREET, BAYSWATER. PO BOX 314, VICTORIA, 3153 AUSTRALIA.  
 TOLL FREE 1800 337 003 | TELEPHONE +61 3 8727 7111  
 FAX +61 3 8727 7222 | info@sdi.com.au | www.sdi.com.au