



LUNA COMPOSITE FAMILY

BLEACH SHADES EQUIVALENCE TO VITA SYSTEM 3D-MASTER®

LUNA UNIVERSAL COMPOSITE

Available in two bleach shades:

- B - Bleach shade: equivalent to OM3 Vita®
- XB - Extra Bleach shade: equivalent to OM2 Vita®

LUNA 2 UNIVERSAL COMPOSITE [BPA free]

LUNA FLOW [Regular viscosity flowable composite]

LUNA FLOW LV [Low viscosity flowable composite]

Available in two bleach shades, being 2XB the lightest of all:

- XB - Extra Bleach shade: equivalent to OM2 Vita®
- 2XB - 2x Extra Bleach shade: equivalent to OM1 Vita®

SDI | LUNA

SDI | LUNA 2

SDI | LUNA FLOW

SDI | LUNA FLOW LV

XB - Extra Bleach

Vita® OM2

SDI | LUNA

B - Bleach

Vita® OM3



SDI | LUNA 2

SDI | LUNA FLOW

SDI | LUNA FLOW LV

2XB - 2 x Extra Bleach

Vita® OM1

