

# SDI PRODUCTS



AMALGAM  
TOOTH WHITENING  
ACCESSORIES  
GLASS IONOMERS  
LED CURING LIGHTS  
ADHESIVES  
SEALANTS



YOUR SMILE. OUR VISION.

# COMPOSITES & SEALANTS

HIGH STRENGTH. LOW SHRINKAGE | PREVENTING DECAY BEFORE IT STARTS

## aura bulk fill >



- One shade, one fill
- High compressive & wear strength - no need for capping layer
- 5mm depth of cure
- Excellent aesthetics due to novel curing feature & chameleon effect
- Low shrinkage
- Quick & easy to place without sticking or slumping



## aura >



- 8 dentine shades and 3 enamel shades to reproduce the natural colours of teeth
- A simple layering technique with zero guesswork
- Superior strength and wear resistance
- High polishability and gloss
- Non-stick and easy to sculpt
- Ultra High Density (UHD) filler with low-shrinkage resin
- High flexural strength



## ice >



ANTERIOR/POSTERIOR  
LIGHT CURED  
RADIOPAQUE NANO  
HYBRID COMPOSITE

- Tooth-like appearance - excellent translucency and natural fluorescence
- Non-stick, non-slump handling



## rok >

POSTERIOR PACKABLE  
COMPOSITE

- Condensable non-stick handling
- Does not set in ambient light or under operative light
- High wear resistance



## aura eASY >



4 SHADES = SIMPLE, BRILLIANT  
RESTORATIONS, VERY TIME!

- 4 shades cover 8 VITA classical A1-D4® shades
- Increased wear resistance and strength
- Smooth high gloss finish every time
- Non-sticky consistency makes it simple to sculpt
- Low volumetric shrinkage to offset marginal stress
- 4 universal shades to minimise inventory
- No slumping and holds its shape
- Available in syringes and complets



## wave, wave mv wave hv >



FLUORIDE RELEASING, FLOWABLE  
RADIOPAQUE COMPOSITES

- High compressive strength
- Nano technology
- Three viscosities



## conseal f >



LIGHT CURED FLUORIDE  
RELEASING RESIN PIT  
AND FISSURE SEALANT

- No Bis-Phenol A or it's derivatives
- Accurate and controlled delivery with 27 gauge tips
- Extremely low viscosity - quickly and deeply penetrate pits and fissures



## conseal clear >

LIGHT CURED PIT AND  
FISSURE RESIN SEALANT

- No Bis-Phenol A or it's derivatives
- Accurate, easy and controlled delivery with 27 gauge tips
- Extremely low viscosity to quickly and deeply penetrate pits and fissures





# GLASS IONOMERS

## THE REVOLUTIONARY GIC SYSTEM

### riva self cure HV >

HIGH VISCOSITY, SELF CURING, TOOTH REMINERALIZING GLASS IONOMER RESTORATIVE MATERIAL

- Bulk fill to save time
- Non-stick handling for easy placement
- Strong dentin replacement



### riva light cure HV >

HIGH VISCOSITY, LIGHT CURED, RESIN REINFORCED GLASS IONOMER RESTORATIVE MATERIAL

- Non-stick handling for easy placement
- No bonding agent needed
- Easy to manipulate and sculpt



### riva self cure >



SELF CURING, TOOTH REMINERALIZING GLASS IONOMER RESTORATIVE MATERIAL

- Bulk fill to save time
- no bonding agent needed!
- BPA & HEMA free



### riva light cure >



LIGHT CURED, RESIN REINFORCED GLASS IONOMER RESTORATIVE MATERIAL

- No bonding agent needed
- Composite like esthetics
- BPA free



### riva protect >



SELF CURING CONVENTIONAL GLASS IONOMER USED AS A SEALANT, LINER AND SURFACE PROTECTOR

- Works in a moist environment – indicated for partially erupted teeth and challenging patients



### riva silver >

MERCURY FREE SELF CURING CONVENTIONAL GLASS IONOMER RESTORATIVE MATERIAL

- Chemically adheres to tooth structure
- High fluoride releasing and recharging ability
- Easy handling



### riva applicator 2 >

- Allows more material to extrude from the capsules
- Autoclavable - easy to clean
- Fits comfortably in smaller hands
- Robust & long lasting



RIVA PRODUCTS ARE AVAILABLE IN BOXES OF 50 CAPSULES OR POWDER / LIQUID KITS

# AMALGAM & EQUIPMENT

SUPERIOR STRENGTH AND HANDLING | HIGH POWERED LED CURING LIGHTS

## permite >

HIGH SILVER NON-GAMMA 2  
ADMIX ALLOY

- Zero microleakage
- Very high early compressive strength
- Not affected by moisture



## radii plus > 2008 REALITY

NOT JUST AN LED CURING LIGHT

- 5 year limited warranty
- Pulse technology = High intensity 1500mW/cm<sup>2</sup> plus low heat emission
- Stays charged for 3h 20 min = 1200 ten second cures



## lojic+ >

PLATINUM MODIFIED,  
HIGH SILVER, HIGH STRENGTH,  
SPHERICAL AMALGAM

- High early compressive strength
- Very low microleakage
- Superior handling characteristics



## radii-cal >

LIGHTWEIGHT, RADICALLY  
POWERFUL LED CURING LIGHT

- Pulse technology = High intensity 1200mW/cm<sup>2</sup> plus low heat emission
- Stays charged for 3h 20 min = 1200 ten second cures
- 2 year limited warranty\*



## gs-80 >

ECONOMICAL NON-GAMMA 2  
ADMIX ALLOY

- High early compressive strength
- Very low microleakage
- Optimal handling characteristics



## ultramatt 2 >

HIGH SPEED MULTI-USE  
TRITURATOR

- 5 year limited warranty
- Quiet, low vibration
- Accurate and consistent trituration



AMALGAM PRODUCTS ARE  
AVAILABLE IN BOXES OF 50  
OR JARS OF 500 CAPSULES

## LED radiometer >

MONITOR FOR CURING  
LIGHT INTEGRITY

- Measures the intensity of curing lights between 400 and 525 nanometers
- Immediate digital reading



# ADHESIVES & ACCESSORIES

AN ADHESIVE FOR EVERY INDICATION | APPLICATORS, TIPS AND MORE

## riva bond LC >



THE FIRST ADHESIVE TO COUNTER COMPOSITE SHRINKAGE STRESS

- Unique Stress Reduction Technology™
- Compensates and absorbs composite shrinkage stress



## go! >



7TH GENERATION ONE STEP, ONE BOTTLE, SELF ETCHING LIGHT CURED ADHESIVE SYSTEM

- No Bis-Phenol A or its derivatives
- Etch, prime and bond in only 35 seconds



## super etch >

37% PHOSPHORIC ACID ETCH GEL

- No interference with adhesive
- Low or high viscosity gel
- Thixotropic gel which flows under direction, not gravity



## points & dispenser >

DISPOSABLE BRUSH APPLICATORS



- Solvent resistant, non-absorbent, non-linting fibers
- Soft, flexible material – bendable brush applicators that retain their bent angle

**Dispenser:**

- Simple turn dispensing
- Economically priced



## stae >



5TH GENERATION TOTAL ETCH - FLUORIDE RELEASING, DENTIN/ ENAMEL ADHESIVE

- High bond strength
- No Bis-Phenol A or its derivatives
- Easy to use - bonds to moist and dry tooth surfaces



## polishing paste >

COMPOSITE POLISHING PASTE

- Compatible with composites, compomers, amalgams, porcelain, ceramic and gold restorations
- No anatomical reduction
- Enhances esthetic finish



## riva conditioner >

POLYACRYLIC ACID CONDITIONER

- Partially removes the smear layer
- Optimizes adhesion of the GIC
- 2.6% polyacrylic acid
- 10mL bottle



## riva coat >

LIGHT CURED COATING MATERIAL FOR GLASS IONOMER AND RESIN PRODUCTS

- Prevents saliva and water contamination to exposed surfaces of glass ionomer cements during the first 24 hours immediately after placement
- 5mL bottle



# TOOTH WHITENING

## ADVANCED TOOTH WHITENING SYSTEMS

### pola day >

HYDROGEN PEROXIDE BASED TAKE HOME ADVANCED TOOTH WHITENING SYSTEM

- From 30 minutes once a day
- Available in 7.5% and 9.5% hydrogen peroxide



### pola office+ >

THE WORLD'S FASTEST BLEACH

- Can be used with or without a light
- Fast and simple to use: 3 x 8 minute applications
- 37.5% hydrogen peroxide



### pola night >

CARBAMIDE PEROXIDE BASED TAKE HOME ADVANCED TOOTH WHITENING SYSTEM

- From 45 minutes once a day
- Available in 10%, 16% and 22% carbamide peroxide



### pola office >

IN-OFFICE TOOTH WHITENING SYSTEM

- Can be used with or without a light
- 30 minutes from preparation to final application
- 35% hydrogen peroxide



### pola day CP >

GREAT RESULTS IN ONLY 15 MINUTES A DAY!

- 35% carbamide peroxide
- Contains fluoride & desensitizing additives
- 2 packaging choices: 4 or 50 syringe kit



### pola luminate >

NEW

BRUSH ON & GO!

- No trays required - quick and easy professional brush on tooth whitener
- 6% hydrogen peroxide, low viscosity fast drying gel
- Apply twice a day for just 30 minutes
- Over 60 applications in a 5.5mL tube
- Contains fluoride to assist with remineralization and reduce sensitivity
- Long lasting spearmint flavor
- Fits easily into a pocket or handbag
- Ideal touch up procedure to maintain tooth whitening between in office or take-home whitening sessions



### pola zing >

IN OFFICE "JUMP START" TOOTH WHITENING SYSTEM

- Shorten the take home tooth whitening process
- 35% carbamide peroxide
- 4 syringe kit



### soothe >

DESENSITIZING GEL

- Instant relief – use before or after whitening process
- 6% potassium nitrate and 0.1% fluoride
- 4 and 50 syringe kits





# MARKETING MATERIALS

FREE



POSTER POLADAY  
Reorder M100274



POSTER POLANIGHT  
Reorder M100275



POSTER POLAOFFICE+  
Reorder M100276



POSTER POLAOFFICE+ AND POLADAY  
Reorder M100279



POSTER POLADAY AND POLANIGHT  
Reorder M100280



POSTER POLAOFFICE+ AND  
POLANIGHT  
Reorder M100281



POLA PATIENT HANDOUTS  
Reorder M100283



POLA PATIENT HANDOUT STAND  
Reorder M100290



POLA BEFORE & AFTER STAND  
Reorder M100288

# CEMENTS

WE HAVE ALL YOUR CEMENT REQUIREMENTS COVERED

## riva luting >

SELF CURING CONVENTIONAL  
GLASS IONOMER LUTING CEMENT

- Great for PFM crowns, metal posts, orthodontic bands, implants
- BPA & HEMA free
- Very low film thickness – flows easily between restoration and tooth



## riva luting plus >

SELF CURING RESIN MODIFIED  
GLASS IONOMER LUTING CEMENT

- Chemically adheres to tooth and substrates
- BPA free
- High bond strengths
- High fluoride releasing and recharging ability



## seT PP >



DUAL CURE, SELF ADHERING  
RESIN CEMENT

- No Bis-Phenol A or HEMA
- No need to etch – minimizes post operative sensitivity
- Easy to clean up
- 7g syringe refills and kits



## seT >



DUAL CURE, SELF ADHERING  
RESIN CEMENT

- No Bis-Phenol A or HEMA
- High bond strength
- Easy clean up
- Box of 50 capsules



# VISIT [WWW.SDI.COM.AU](http://WWW.SDI.COM.AU)

Pola Office+ the world's fastest bleach: Based on the total treatment time compared to that of all other competitor's printed literature. Radium Plus: 5 year warranty on all parts except battery and full arch bleaching LED attachment which is 2 years. Radium-Cal: 2 year manufacturing warranty on parts and battery. Pola Office+: THE DENTAL ADVISOR, Vol. 25, No. 9, November 2008. Pola Office: THE DENTAL ADVISOR, Vol. 23, No. 3, April 2006. Soothe: THE DENTAL ADVISOR, Vol. 23, No. 6, July/Aug 2006. Go!: THE DENTAL ADVISOR, Vol. 24, No.10, Dec 2007. Stae: THE DENTAL ADVISOR, Vol. 15, 1998. Points: THE DENTAL ADVISOR, Vol. 20, 2003. Permite: THE DENTAL ADVISOR, Vol. 5, No. 2, March/April 1995. Radium Plus: THE DENTAL ADVISOR, Vol. 23, No. 3, April 2006. Lojic+: THE DENTAL ADVISOR, Vol. 4, No. 5, Sep/Oct 1994. Ultramat 2: THE DENTAL ADVISOR editors choice Vol 28, No. 9, Nov 2011. Riva Light Cure: THE DENTAL ADVISOR, Vol.27, No. 6, July-Aug 2010. Conseal F: THE DENTAL ADVISOR, Vol. 20, 2003. Ice: THE DENTAL ADVISOR, Vol. 21, No. 6, 2004. Wave (wave only): THE DENTAL ADVISOR, Vol. 15, 1998. seT: THE DENTAL ADVISOR, Vol. 25, No. 8, October 2008. seT PP: THE DENTAL ADVISOR, Vol.27, No. 7, Sept 2010.



YOUR SMILE. OUR VISION.

[www.sdi.com.au](http://www.sdi.com.au)

Australia 1 800 337 003  
Austria 00800 022 55 734  
Brazil 0800 770 1735  
France 00800 022 55 734  
Germany 0800 100 5759  
Ireland 01 886 9577

Italy 800 780 625  
New Zealand 0800 734 034  
Spain 00800 022 55 734  
United Kingdom 00800 022 55 734  
USA & Canada 1 800 228 5166

MX00012 G



+D036MX00012H